

MABSOL[®] Biotinylated Proteins

Bring pre-labeled biotinylated proteins directly to your bench



ELISA



FAC Sorting



Biopanning



SPR



Enrichment



More

PD-1

TNF-alpha

HER2

CD3E

FcRn

BCMA

CD16a

Content


- I. Introduction**
- II. Application**
- III. Product Series**
- IV. Case Studies**
- V. PrecisionAvi Series Product List**
- VI. UltraLys Series Product List**
- VII. Testimonials**



OX40 ROR1 PD-L1 TIGIT EGFR CD19 CTLA-4

Introduction

Biotin is commonly used as a protein tag to facilitate the detection, purification, and immobilization of the proteins.



The bond between biotin and its binding partner Avidin (or Streptavidin) is unique in the following ways:

- Strong (K_a of 10^{15} M^{-1})
- Specific
- Multi-moiety
- Stable
- Minimal interference

Fig. 1 Biotin-avidin binding structure

With characteristics mentioned above, biotin-avidin (streptavidin) system is now considered a versatile independent technology in following applications.

Application

■ ELISA

Biotinylated proteins can be used in ELISA as two-way antibodies for both capture and detection with high specificity and detection sensitivity.

■ FACS

Biotinylated proteins can be used along fluorophore-tagged SA to detect/isolate cells expressing particular surface markers.

■ Biopanning

It is a technique often used for the selection of phage displays during antibody drug development. Biotinylated proteins can be used with SA-coated magnetic beads/surface in biopanning with higher coating density and uniformed antigen presentation.

■ SPR

It is a standard method used by pharmaceutical researchers to study protein binding kinetics. Biotinylated proteins can be used along with Biacore Sensor Chip SA for SPR analysis with low baseline drift and noise, and low activity loss during surface regeneration cycle.

■ Immuno-capture and Enrichment

Biotinylated proteins can be used to isolate antibodies from plasma or other biological fluid for subsequent analyses with high sensitivity, and the processed sample can be easily analyzed in quantitative mass spectrometry.

Product Series


MABSOL® biotinylated protein collection includes two unique and complimentary product series, the PrecisionAvi series built upon the Avitag™ technology, and the UltraLys series produced using the in-house developed chemical labeling method. These products are made with every attention to details.

■ PrecisionAvi Series

An exclusive collection of ready-to-use Avitag™ biotinylated proteins

The products in this series are exclusively produced using the Avitag™ technology. Briefly, a unique 15 amino acid peptide, the Avi tag, is introduced into the recombinant protein during expression vector construction. The single lysine residue in the Avi tag is enzymatically biotinylated by the *E. Coli* biotin ligase BirA.


This single-point enzymatic labeling technique brings many advantages for commonly used binding assays. The biotinylation happens on the lysine residue of Avi tag, and therefore does NOT interfere with the target protein's natural binding activities. In addition, when immobilized on an avidin-coated surface, the protein orientation is uniform because the position of the Avi tag in the protein is precisely controlled since biotin is labeled at single lysine in the Avi tag.



■ UltraLys Series

A unique series of chemically labeled biotinylated proteins with ultra-sensitivity

The products in this series are produced using our in-house developed chemical labeling approach. The primary amines in the side chains of lysine residues and the N-terminus of protein are conjugated with biotins.



Chemical labeling usually results in multiple biotin attachment on a single protein molecule, which could potentially leads to higher detection sensitivity.

Case Studies

■ SPR: Antibody Optimization with Biotinylated FcRn

The half-life and efficacy of a therapeutic antibody largely depends on its Fc fragment. To obtain candidate antibody with desired pharmaceutical properties, the interaction between Fc fragment and Fc receptors needs to be evaluated. SPR is a common assay used for determining binding affinity between Fc and Fc receptors.

Efficient coating and chip regeneration present a major challenge in this approach. As shown here, the combo of ACROBiosystems' Avitag™ biotin-labeled FcRn and GE's sensor Chip-SA delivers satisfactory result in a binding assay against Herceptin.




Fig. 2 Biotinylated Human FcRn / FCGRT & B2M, Avi Tag (Avitag™) (Cat. # FCM-H82W4) coupled to SA coated sensor chip can bind Herceptin with an affinity constant of 2.24 μ M as determined in a SPR assay (Biacore 2000). The data is provided by Biaffin GmbH & Co KG, Germany.

■ Immuno-capture with Biotinylated TNF-alpha

TNF-alpha antibody is widely used in treating autoimmune diseases. In order to optimize the administration, it's very important to monitor the serum concentration of TNF-alpha antibody. Our collaborator AbSciex has shown that, using ACROBiosystems's biotinylated TNF-alpha (Cat. # TNA-H821R), TNF-alpha antibody in the serum could be easily detected by immune-capture coupled with quantitative MS/MS.




Fig. 3 Chromatograms spiked Adalimumab in human plasma: blank plasma, 5, 50, and 500 ng/ml

Biotinylated TNF-alpha is first coupled to streptavidin coated magnetic beads, and then capture TNF-alpha antibody in the serum. After that, MS/MS is applied to quantify TNF-alpha antibody. Immuno-capture with biotinylated TNF-alpha significantly increases the sensitivity by ten folds.

■ Inhibitor Screening with Biotinylated PD-1

ACROBiosystems has developed a PD1-PD-L1 inhibitor screening assay kit (Cat. # EP-101) for rapid and high throughput screening of candidate inhibitory antibodies or small molecules of the PD1 pathway.

This inhibitor screening ELISA pair is designed to facilitate the identification and characterization of new PD-1 pathway inhibitors. The assay takes advantage of our in-house developed binding of biotinylated human PD-1 to immobilized human PD-L1 in a functional ELISA assay, and employs a simple colorimetric sandwich ELISA platform. Briefly, we provide you with a human PD-1-Biotin protein, a human PD-L1 protein, an anti-PD-1 neutralizing antibody (as method verified Std.), and streptavidin-HRP reagent. Both biotinylated PD-1 and PD-L1 proteins are expressed in the HEK293 cells.




Fig. 4 Immobilized human PD-L1 protein at 2 µg/mL (100 µL/well) can bind biotinylated human PD-1 with a linear range of 0.038 - 0.6 µg/mL when detected by Streptavidin-HRP.




Fig. 5 Inhibition of PD-1-PD-L1 binding by Anti-PD-1 neutralizing antibody is measured using the PD-1 [Biotinylated] : PD-L1 Inhibitor Screening ELISA Assay Pair (Cat. # EP-101).

PrecisionAvi Series (Avitag™) Product List

Molecule	Cat. No.	Species	Structure	Size
4-1BB	41B-H82F7	Human	4-1BB (24-186) Fc Avi	25ug, 200ug,
4-1BB	41B-H82E3	Human	4-1BB (24-186) Avi His	25ug, 200ug,
B7-1	B71-H82F2	Human	B7-1 (35-242) Fc Avi	25ug, 200ug,
B7-2	CD6-H82E2	Human	B7-2(26-247) Avi His	25ug, 200ug,
B7-2	CD6-H82F5	Human	B7-2 (26-247) Fc Avi	25ug, 200ug,
B7-H4	B74-H82E2	Human	B7-H4(29-258) Avi His	25ug, 200ug,
B7-H5	B75-H82E1	Human	B7-H5 (33-194) Avi His	25ug, 200ug,
BAFF	BAF-H82F3	Human	Avi Fc BAFF (134-285)	25ug, 200ug,
BAFF	BAF-H82Q2	Human	His Avi BAFF (134-285)	25ug, 200ug,
BAFFR	BAR-M82F0	Mouse	BAFFR (10-71) Fc Avi	25ug, 200ug,

PrecisionAvi Series (Avitag™) Product List

Molecule	Cat. No.	Species	Structure	Size
BCMA	BC7-H82F0	 Human	BCMA(1-54) Fc Avi	25ug, 200ug,
BTLA	BTA-H82F8	 Human	BTLA (31-134) Fc Avi	25ug, 200ug,
CD155	CD5-H82F6	 Human	CD155 (21-343) Fc Avi	25ug, 200ug,
CD27	TN7-H82F6	 Human	CD27 (21-192) Fc Avi	25ug, 200ug,
CD28	CD8-H82F2	 Human	CD28 (19-152) Fc Avi	25ug, 200ug,
CD40	CD0-H82E8	 Human	CD40 (21-193) Avi His	25ug, 200ug,
CD40	CD0-H82F2	 Human	CD40 (21-193) Fc Avi His	25ug, 200ug,
CD47	CD7-H82F6	 Human	CD47 (19-139) Fc Avi	25ug, 200ug,
CD47	CD7-H82E9	 Human	CD47(19-139) His Avi	25ug, 200ug,
CTLA-4	CT4-H82F3	 Human	CTLA-4 (37-162) Fc Avi	25ug, 200ug,
CTLA-4	CT4-C82E5	 Cynomolgus	CTLA-4 (37-160) Avi His	25ug, 200ug,
DNAM-1	DN1-H82F9	 Human	DNAM-1 (19-247) Fc Avi	25ug, 200ug,
EGFR	EGR-H82E7	 Human	EGFR (25-645) Avi His	25ug, 200ug,
EGFR (vIII)	EGR-H82E0	 Human	EGFR(25-29) EGFR (298-645) Avi His	25ug, 200ug,
Fc gamma RI / CD64	CD4-M52E8	 Mouse	CD64 (25-297) Avi His	25ug, 200ug,











PrecisionAvi Series (Avitag™) Product List

Molecule	Cat. No.	Species	Structure	Size
Fc gamma RIIA / CD32a	CDA-H82E6	 Human	^{H131} CD32a (36-218) Avi His	25ug, 200ug,
Fc gamma RIIA / CD32a	CDA-H82E7	 Human	^{R131} CD32a (35-217) Avi His	25ug, 200ug,
Fc gamma RIIIB / CD32b	CDB-H82E0	 Human	CD32b (46-217) Avi His	25ug, 200ug,
Fc gamma RIIIA / CD16a	CDA-H82E8	 Human	^{F158} CD16a (17-208) Avi His	25ug, 200ug,
Fc gamma RIIIA / CD16a	CDA-H82E9	 Human	^{V158} CD16a (17-208) Avi His	25ug, 200ug,
Fc gamma RIIIB / CD16b	CDB-H82E1	 Human	CD16b (17-200) Avi His	25ug, 200ug,
FcRn	FCM-H82W4	 Human	FCGRT (24-297) Avi His B2M (21-119) Strep II	25ug, 200ug,
FcRn	FCM-M82W6	 Mouse	FCGRT (22-297) Avi His B2M (21-119) Twin-Strep	25ug, 200ug,
FcRn	FCM-R82W7	 Rat	FCGRT (23-298) Avi His B2M (21-119) Strep II	25ug, 200ug,
FcRn	FCM-C82W5	 Cynomolgus / Rhesus	FCGRT (24-297) Avi His B2M (21-119) Strep II	25ug, 200ug,
G-CSF R	GCR-H82E4	 Human	G-CSF R (25-621) Avi His	25ug, 200ug,
Hemagglutinin (HA)	HA1-V82E1	 Influenza	Hemagglutinin (HA) (19-523) Avi His	25ug, 200ug,
Hemagglutinin (HA)	HA1-V82E2	 Influenza	Hemagglutinin (HA) (19-338) Avi His	25ug, 200ug,
Hemagglutinin (HA)	HA1-V82E4	 Influenza	Hemagglutinin (HA) (17-339) Avi His	25ug, 200ug,





PrecisionAvi Series (Avitag™) Product List

Molecule	Cat. No.	Species	Structure	Size
HGF R	MET-H82E1	 Human	HGF R (25-932) Avi His	25ug, 200ug,
HVEM	HV4-H82F1	 Human	HVEM (39-202) Fc Avi	25ug, 200ug,
IGF-II	IG2-H82F9	 Human	Avi Fc IGF-II (25-91)	25ug, 200ug,
IgG Fc	IG1-H82E2	 Human	Human IgG1 Fc (100-330) Avi His	25ug, 500ug,
IgG Fc	IGA-M8210	 Mouse	Mouse IgG2a Fc (98-330) Avi	25ug, 500ug,
IL-17 RA	ILA-H82F1	 Human	IL-17 RA (33-320) Fc Avi	25ug, 200ug,
IL-17 RA	ILR-H82E0	 Human	IL-17 RA (33-320) Avi His	25ug, 200ug,
IL23A & IL12B	ILB-H82W6	 Human	His Avi IL23A (20-189) IL12B (23-328)	25ug, 200ug,
IL-4 R alpha	ILR-H82E9	 Human	IL-4 R alpha (26-232) Avi His	25ug, 200ug,
IL-7 R alpha	IL7-H82F8	 Human	IL-7 R alpha (21-236) Fc Avi	25ug, 200ug,
OX40	OX0-H82F7	 Human	OX40 (29-216) Fc Avi	25ug, 200ug,
PCSK9	PC9-H82E7	 Human	PCSK9 (31-692) Avi His	25ug, 200ug,
PCSK9	PCY-H82E7	 Human	^{D374Y} PCSK9 (31-692) Avi His	25ug, 200ug,
PCSK9	PC9-M82E1	 Mouse	PCSK9 (35-694) Avi His	25ug, 200ug,
PD-1	PD1-H82E4	 Human	PD-1 (25-167) Avi His	25ug, 200ug,





PrecisionAvi Series (Avitag™) Product List

Molecule	Cat. No.	Species	Structure	Size
PD-1	PD1-H82F2	 Human	PD-1 (25-167) Fc Avi His	25ug, 200ug,
PD-1	PD1-M82F4	 Mouse	PD-1 (25-167) Fc Avi	25ug, 200ug,
PD-L1	PD1-H82F3	 Human	PD-L1 (19-238) Fc Avi His	25ug, 200ug,
PD-L1	PD1-H82E5	 Human	PD-L1 (19-238) Avi His	25ug, 200ug,
PD-L1	PD1-M82F5	 Mouse	PD-L1 (19-238) Fc Avi	25ug, 200ug,
PD-L2	PD2-H82E8	 Human	PD-L2 (20-219) Avi His	25ug, 200ug,
PD-L2	PD2-H82F6	 Human	PD-L2 (20-219) Fc Avi	25ug, 200ug,
ROR1	RO1-H82F4	 Human	ROR1 (30-403) Fc Avi	25ug, 200ug,
SCF	SCF-H82E1	 Human	SCF (26-190) Avi His	25ug, 200ug,
Siglec-3	CD3-H82E7	 Human	CD33 (18-259) Avi His	25ug, 200ug,
TIGIT	TIT-H82F1	 Human	TIGIT (22-141) Fc Avi	25ug, 200ug,
TIGIT	TIT-H82E5	 Human	TIGIT (22-141) Avi His	25ug, 200ug,
TNF-alpha (HPLC-verified)	TNA-H82E3	 Human	TNF-alpha (77-233) His Avi	25ug, 200ug,
TNF-alpha (HPLC-verified)	TNA-M82E9	 Mouse	TNF-alpha (80-235) His Avi	25ug, 200ug,
TNFR1	TN1-H82E3	 Human	TNFR1 (22-211) Avi His	25ug, 200ug,

PrecisionAvi Series (Avitag™) Product List

Molecule	Cat. No.	Species	Structure	Size
VEGF R2	KDR-H82E5	 Human	VEGF R2 (20-764) Avi His	25ug, 200ug,
VEGF121	VE1-H82E7	 Human	Avi His VEGF121 (27-147)	25ug, 200ug,
VEGF164	VE4-M82Q3	 Mouse	His Avi VEGF164 (27-190)	25ug, 200ug,
VEGF165	VE5-H82Q0	 Human	His Avi VEGF165 (27-191)	25ug, 200ug,

UltraLys Series (Chemical labeling) Product List

Molecule	Cat. No.	Species	Structure	Size
B7-H4	B74-H8222	 Human	B7- H4 (29-258) His	25ug, 200ug
CD19	CD9-H8259	 Human	CD19 (20-291) Fc	25ug, 200ug
CD3E & CD3D	CDD-H82W0	 Human	CD3E (23-126) Fc His	25ug, 200ug
			CD3D (22-105) Fc Flag	
CD3 epsilon	CDE-H8223	 Human	CD3 epsilon (23-126) His	25ug, 200ug
CX3CL1	CX1-H8221	 Human	CX3CL1 (25-341) His	25ug, 200ug
EphB4	EP4-H8229	 Human	EphB4 (16-539) His	25ug, 200ug
EpCAM	EPM-H8223	 Human	EpCAM (24-265) His	25ug, 200ug
EpCAM	EPM-H8254	 Human	EpCAM (24-265) Fc	25ug, 200ug
ErbB3	ER3-H8223	 Human	ErbB3 (20-643) His	25ug, 200ug
FcRn	FCM-H8286	 Human	FCGRT (24-297) His	25ug, 200ug
			B2M (21-119) Strep II	
FGF basic	BFF-H8117	 Human	FGF basic (143-288)	50ug, 500ug
Growth Hormone R	GHR-H8222	 Human	Growth Hormone R (27-264) His	50ug, 200ug
GM-CSF	GMF-H8214	 Human	GM-CSF (18-144)	25ug, 200ug
GPA33	GP3-H8224	 Human	GPA33 (22-235) His	25ug, 200ug

UltraLys Series (Chemical labeling) Product List

Molecule	Cat. No.	Species	Structure	Size
Her2	HE2-H822R	 Human	Her2 (23-652) His	25ug, 200ug
IL-6	IL6-H8218	 Human	IL-6 (30-212)	25ug, 200ug
IL-1 alpha	ILA-H8213	 Human	IL-1 alpha (113-271)	25ug, 200ug
Mesothelin	MSN-H8223	 Human	Mesothelin (296-580) His	25ug, 200ug
Mesothelin	MSN-H826x	 Human	Fc Mesothelin (296-580)	25ug, 200ug
PD-1	PD1-M8259	 Mouse	PD-1 (25-167) Fc	25ug, 200ug
Protein L	RPL-P814R	N/A	N/A	500ug, 2mg
SCF	SCF-M8228	 Mouse	SCF (26-189) His	25ug, 200ug
SOST	SOT-H8245	 Human	His SOST (24-213)	25ug, 200ug
TFPI	TFI-H8226	 Human	TFPI (29-282) His	25ug, 200ug
Transferrin R	TFR-H8243	 Human	His Transferrin R (89-760)	25ug, 200ug
TNF-alpha	TNA-H8211	 Human	TNF-alpha (77-233)	25ug, 200ug
TNF-alpha	TNA-H821R	 Human	TNF-alpha (77-233)	25ug, 200ug
VEGF165	VE5-H8210	 Human	VEGF165 (27-191)	25ug, 200ug

Testimonials

"I have been very satisfied with ACROBiosystems' products and customer service. Their technical support team has been extraordinarily proactive in making sure that our needs are being fully addressed. They always work closely with us to help us find the products that suited our research the best. This includes offering an array of samples for free for our preliminary testing. The demonstration of such awareness of subtle issues surrounding our research was quite unexpected and highly impressive.

Vidal F.de la Cruz, Founder & Principal of Fomento Pharma, LLC

"Rubicon Biotechnology considers ACROBiosystems a preferred vendor to supply us with high quality proteins for our drug development programs, including using their proteins in our in vivo studies. We have contracted with ACROBiosystems to manufacture high purity custom proteins and they proposed, and met, aggressive delivery times that allowed Rubicon to start some critical proof-of-concept studies. Acro's customer service was extremely knowledgeable and responsive and kept us up to date on their progress. We never felt out of the loop. We are extremely satisfied with ACROBiosystems and would recommend them to anyone in the protein discovery and development fields.

Richard Richieri, Founding Partner & Chief Operating Officer, Rubicon Biotechnology

"TNF-alpha immunocapture is easy to automate and offers a very specific sample preparation by isolating only the compound of interest. Furthermore the high specificity generates very low and robust MS/MS signal to noise ratios and has allowed us to gain a factor of ten in sensitivity. Thus our mass spectrometer which hitherto was not sensitive enough, can now detect levels of concentration of our therapeutic monoclonal antibody which are compatible with the clinical diagnosis.

Jean-François Jourdil, University Hospital of Grenoble

"In my ten years of research experience, I have found that the protein products produced by ACROBiosystems are equal to none. They are exceptional in their purity and activity and, as such, have greatly enhanced the quality of my research. It has been a pleasure to find a great multitude of proteins, most with a variety of tags, all offered by a single company. In addition, the efficiency and cordiality of their customer service matches the excellence of their products, resulting in one very satisfied customer!

Adam D. Friedman, University of North Carolina at Chapel Hill

CD3E

HER2

CD16a

BCMA

FcRn

PD-1

TNF-alpha

TIM-3 CTLA-4 4-1BB
Immune Checkpoint
Biotinylated HER2
B7-H4 TIGIT 4-1BB Ligand
CD40 GTR
PD-L2 HER2
DNAM-1 B7-1
TIM-3 FcRn
TNF- alpha LAG-3
CD19 OX40 Ligand B7-H2
PD-1 CD47 PCSK9
Immune Checkpoint
Biotin-labeled
VEGF165
FcRn HER2 ICOS
B7-H4 BTLA CD27 CD48
TIM-3 CTLA-4
Biotin-labeled
HER2 PD-L2 B7-H4 DNAM-1
VEGF165 FcRn PD-L1
BTLA CTLA-4
DNAM-1 PD-L2 FcRn
GTR Ligand Biotinylated PD-1
4-1BB
TNF-alpha PD-L2