



FuGENE SI Quick Protocol

The following quick protocol will generate enough mix to transfect one 35mm vessel, five wells of a 24-well plate, or twenty-five wells of a 96-well plate (at the recommended starting ratio of 0.3ul FuGENE® SI + 1pmol total siRNA per well of 96 well plate, please consult Table 1 in users guide for other vessel sizes and scalings)

Preparing the FuGENE SI® Transfection Reagent

- 1. Before use, allow the vial of FuGENE® SI Transfection Reagent, siRNA stock, and serumfree medium to reach room temperature
- 2. Mix by inverting or vortexing briefly. If a precipitate is visible, briefly warm at 37 degrees C then cool to room temperature

General Transfection Protocol

- 1. Seed cells to be 25-50% confluent.
- 2. Dilute 7.5ul of FuGENE SI reagent in 117.5ul of pre-warmed serum free medium, vortex 1 sec.
- 3. Dilute 2.5ul of 10uM siRNA stock in 122.5ul of pre-warmed serum free medium, vortex 1 sec.
- 4. Combine 125ul of FuGENE SI dilution from step #2 with 125ul of siRNA dilution from step #3.
- 5. Incubate FuGENE SI + siRNA complex for 5 minutes at room temperature. (up to 15 min)
- 6. Add FuGENE SI + siRNA complex to cells in drop-wise manner, swirl or shake to mix
 - i. 35mm vessel: Add all 250ul to vessel
 - ii. 24-well plate: Add 50ul per well
 - iii. 96-well plate: Add 10ul per well
- 7. Incubate transfected cells for 24-72 hours
- 8. Analyze transfected cells.
- 9. See additional protocol information in full users guide available at www.fugene.com
- 10. For additional support please contact us at contact@fugene.com or visit us at www.fugene.com



Learn more at: www.fugene.com