

MOLECULAR SIGNATURE®

## DNA Damage (8-OHdG) ELISA kit

Colorimetric detection of 8-hydroxy-2-deoxy Guanosine  
Catalog No. SKT-120



Discovery through partnership | Excellence through  
quality

### Product Name

DNA Damage (8-OHdG) ELISA kit

### Description

Colorimetric detection of 8-hydroxy-2-deoxy Guanosine

### Species Reactivity

Species Independent

### Platform

Microplate

### Sample Types

Cell lysates, Plasma, Sample matrices, Urine

### Detection Method

Colorimetric Assay

### Assay Type

Competitive ELISA (Enzyme-linked Immunosorbent Assay)

### Utility

ELISA Kit for 8-OHdG detection in samples.

### Sensitivity

0.59 ng/mL

### Assay Range

0.94 - 60 ng/mL

### Precision

Intra-Assay Precision: Three samples of known concentration were assayed thirty times on one plate; the intra-assay coefficient of variation of the DNA Damage ELISA has been determined to be <5%. Inter-Assay Precision: Three samples of known concentration were assayed thirty times in three individual assays; the inter-assay coefficient of variation of the DNA Damage ELISA has been determined to be <5%.

---

**Incubation Time**

---

1 hour

---

**Number Of Samples**

---

39 samples in duplicate

---

**Other Resources**

---

Kit Booklet , MSDS , Calculations Worksheet

---

**Field Of Use**

---

Not for use in humans. Not for use in diagnostics or therapeutics. For in vitro research use only.

---

## Properties

---

---

**Storage Temperature**

4°C and -20°C

---

**Shipping Temperature**

---

Blue Ice

---

**Product Type**

---

ELISA Kits

---

**Specificity**

---

Cross-reactivity: 8-Hydroxy-2-deoxy Guanosine (8-OHdG): 100%. 8-Hydroxy Guanosine (8-OHG): 23%. 8-Hydroxy Guanine (8-oxoG): 23%. Guanosine: <0.01%.

---

**Assay Overview**

---

1. Prepare standard and samples in the Sample and Standard Diluent. 2. Add 50 µL of prepared standards and samples in triplicate to appropriate wells. 3. Add 50 µL of the diluted antibody preparation to the appropriate wells. 4. Cover plate with Plate Cover and incubate at room temperature (20-25°C) for 1 hour. 5. Wash plate 4 times with 1X Wash Buffer. 6. Add 100 µL of TMB Substrate to each well. 7. Cover plate and develop the plate in the dark at room temperature for 30 minutes. 8. Add 100 µL of Stop Solution to each well. 9.

Measure absorbance on a plate reader at 450 nm. 10. Plot the standard curve and calculate sample concentrations.

---

## Kit Overview

---

### Component No.

#### Item

#### Quantity / Size

---

#### SKC-120A

---

8-hydroxy-2-deoxy Guanosine : BSA Coated Plate

---

1 Plate

---

#### SKC-120C

---

8-hydroxy-2-deoxy Guanosine Standard

---

1 vial/ 100uL

---

#### SKC-120F

---

8-hydroxy-2-deoxy Guanosine HRP Conjugated Monoclonal Antibody

---

1 vial/75uL

---

#### SKC-0001

---

Sample and Standard Diluent

---

1 vial/50mL

---

#### SKC-0002

---

8-hydroxy-2-deoxy Guanosine Antibody Diluent

---

1 vial/13mL

---

#### SKC-0003

---

Wash Buffer Concentrate

---

1 vial/50mL

---

#### SKC-0004

---

TMB Substrate

---

1 vial/13mL

---

---

**SKC-0005**

---

Stop Solution

---

1 vial/13mL

---

**SKC-0009**

---

Plate Cover

---

2 covers

---

---

**Cite This Product**

---

DNA Damage (8-OHdG) ELISA Kit (StressMarq Biosciences Inc., Victoria BC CANADA, Catalog # SKT-120)

---

## Biological Description

---

**Alternative Names**

8-OH-dG ELISA Kit, 8OHG ELISA Kit, 8OG ELISA Kit, 8 hydroxyguanine ELISA Kit, 8-OHdG ELISA Kit, DNA Damage ELISA Kit

---

**Research Areas**

---

Cancer, Cell Signaling, Oxidation, Oxidative Stress, Post-translational Modifications

---

**Scientific Background**

---

8-hydroxy-2-deoxy Guanosine (8-OH-dG) is produced by the oxidative damage of DNA by reactive oxygen and nitrogen species and serves as an established marker of oxidative stress (1-4). Hydroxylation of guanosine occurs in response to both normal metabolic processes and a variety of environmental factors (i.e., anything that increases reactive oxygen and nitrogen species). Increased levels of 8-OH-dG are associated with the aging process as well as with a number of pathological conditions including cancer, diabetes, and hypertension(5-9). In complex samples such as plasma, cell lysates, and tissues, 8-OH-dG can exist as either the free nucleoside or incorporated in DNA. Once the blood enters the kidney, free 8-OH-dG is readily filtered into the urine, while larger DNA fragments remain in the bloodstream. Because of the complexity of plasma samples, urine is a more suitable matrix for the measurement of free 8-OH-dG than plasma. Urinary levels of 8-OH-dG range between 2.7-13 ng/mg creatine, while plasma levels of free 8-OH-dG have been reported to be between 4-21 pg/ml as determined by LC-MS (10-11).



---


**References**

---


1. Maxey K.M., Maddipati K.R., Birkmeier J. (1992) J Clin Immunoassay 15: 116-120.
  2. Pradelles P., Grassi J., Maclouf J. (1990) Methods Enzymol. 187: 24-34.
  3. Maclouf J., Grassi J., Pradelles P. (1987) Dev Immunoassay Tech Meas eicosanoids.
  4. Lin H., et al. (2004) Biochem J. 380: 541-548.
  5. Bogdanov M.B., et al. (1999) Free Radic Biol Med. 27(5/6): 647-666.
  6. Lee J., et al. (2005) Hypertension 45: 986-990.
  7. Leinonen, J., et al. (1997) FEBSLett. 417: 150-152.
  8. Endo K., et al. (2006) J. Atheroscler. Thromb. 13:68-75.
  9. Kuo H., et al. (2007) Mutat Res. 631:62-68.
  10. Shen J., et al. (2007) Cancer 109: 574-580.
  11. Beckman K.B., Ames B.N. (1997) J Biol Chem 272: 19633-19636.
  12. Epe B., et al. (1996) Nucleic Acids Res 24: 4105-4110.
  13. Spencer J.P.E., et al. (1995) FEBS Lett 374: 233-236.
  14. Floyd R.A. (1990) FASEB J 4: 2587-2597.
- 

## Product Images


---




Typical Standard Curve for the DNA Damage (8-OHdG) ELISA kit (Enzyme-Linked Immunosorbent Assay) StressXpress® – SKT-120. Assay Type: Competitive ELISA. Detection Method: Colorimetric Assay. Assay Range: 0.94 – 60 ng/ml.




## Chemical Equation of the Oxidation of Guanosine



## Diagram of the 8-OHdG Competitive ELISA



## Urine Spike Assay




## Diagram of the Preparation of the 8-OHdG Standards

	1	2	3	4	5	6	7	8	9	10	11	12
A	S1	S1	S1	Blk	Blk	Blk	8	8	8	16	16	16
B	S2	S2	S2	1	1	1	9	9	9	17	17	17
C	S3	S3	S3	2	2	2	10	10	10	18	18	18
D	S4	S4	S4	3	3	3	11	11	11	19	19	19
E	S5	S5	S5	4	4	4	12	12	12	20	20	20
F	S6	S6	S6	5	5	5	13	13	13	21	21	21
G	S7	S7	S7	6	6	6	14	14	14	22	22	22
H	S8	S8	S8	7	7	7	15	15	15	23	23	23

S1 – S7: 60 to 0.94 ng/ml Standards  
S8: Zero Standard

Blk: Blank  
1 – 23: Samples

Diagram of the Triplicate Sample Plate Format



StressMarq  
Biosciences Inc.

Discovery through partnership | Excellence through quality

**DNA/RNA (8-OHdG) Damage ELISA Kit Worksheet**

Catalog No. SKT-120-965/SKT-120-4805

118-1537 Millside Ave | Victoria BC V8T 3C1 CANADA

Tel: 250-294-9880 | Fax: 250-294-9025 | Email: info@stressmarq.com

www.stressmarq.com

**STEP 1:** Fill in your Blank Optical Density (OD) values. Mean OD will automatically calculate the mean of the four values. The Blank Mean OD will be automatically subtracted from your subsequent OD values provided in Steps 2 & 3.


Step 1				
For Duplicate				
	Well 1 OD <sub>mean</sub>	Well 2 OD <sub>mean</sub>	Well 3 OD <sub>mean</sub>	Mean OD <sub>mean</sub>
Blank				

**STEP 2:** Fill in your Standard OD values. Calculated [8-OHdG] will calculate the [8-OHdG] of each standard used based on the equation that your standards have generated. These should be very similar to the [8-OHdG] from your standard values.

Standard Curve						
For Duplicate						
Standard No.	[8-OHdG] (ng/ml)	Well 1 OD <sub>mean</sub>	Well 2 OD <sub>mean</sub>	Well 3 OD <sub>mean</sub>	Mean OD <sub>mean</sub>	Calculated [8-OHdG] (ng/ml)
1	60.00					
2	30.00					
3	15.00					
4	7.50					
5	3.75					
6	1.88					
7	0.94					

Constants for  $y = ax^3 + bx^2 + cx + d$

a	b	c	d	r <sup>2</sup>



**STEP 3:** Fill in your sample OD values, and dilution factor. If you diluted your samples 1:20, then input 20 into the Dilution Factor column. The Corrected [8-OHdG] will calculate the concentration of your samples before dilution by accounting for your

Unknown Samples					
	Well 1 OD <sub>mean</sub>	Well 2 OD <sub>mean</sub>	Well 3 OD <sub>mean</sub>	Mean OD <sub>mean</sub>	Corrected [8-OHdG] (ng/ml)

Calculations Worksheet | 1 |

Preview of the Calculations Worksheet

## Product Citations (53)

### Other Citations

#### Association of serum total bilirubin and plasma 8-OHdG in HIV/AIDS patients

Kolgiri, V., Nagar, V. and Patil, V. (2018) Interventional Medicine and Applied Science [Epub ahead of print].



**PubMed ID:** N/A    **Reactivity:** Human    **Applications:** Blood

**The biological age linked to oxidative stress modifies breast cancer aggressiveness.**

Saez Freire, M.D.M. et al. (2018) Free Radic Biol Med. [Epub ahead of print].

**PubMed ID:** 29550329    **Reactivity:** Mouse    **Applications:** Liver

**SOCS1-targeted therapy ameliorates renal and vascular oxidative stress in diabetes via STAT1 and PI3K inhibition.**

Lopez-Sanz, L. et al. (2018) Lab Invest. [Epub ahead of print].

**PubMed ID:** 29540859    **Reactivity:** Mouse    **Applications:** Blood and urine

**Effects of ozone repeated short exposures on the airway/lung inflammation, airway hyperresponsiveness and mucus production in a mouse model of ovalbumin-induced asthma.**

Bao, W. et al. (2018) Biomed Pharmacother. [Epub ahead of print].

**PubMed ID:** 29499403    **Reactivity:** Mouse    **Applications:** Lung

**Ginsenoside Rb1 prevents high glucose-induced Schwann cell injury through the mitochondrial apoptosis pathway.**

Li, X. et al. (2018) J Tradit Chin Med. 37(6): 746-755.

**PubMed ID:** N/A    **Reactivity:** Rat    **Applications:** Schwann cells

**Urinary 8-hydroxydeoxyguanosine is a better biomarker of aging in non-smokers.**

Zhao, Z.X. et al. (2018) Tradit. Med. Mod. Med.

**PubMed ID:** N/A    **Reactivity:** Human    **Applications:** Urine

**IL-6 receptor blockade ameliorates diabetic nephropathy via inhibiting inflammasome in mice.**

Wu, R. et al. (2018) Metabolism. [Epub ahead of print].

**PubMed ID:** 29336982    **Reactivity:** Mouse    **Applications:** Urine

**Analysis of Oxidative DNA Damage in HIV-Positive Pregnant Women.**

Rajekar, R.B. (2017) IJIRMS. 02(10):2455-8737.

**PubMed ID:** N/A    **Reactivity:** Human    **Applications:** Plasma

**Gene expression, oxidative stress and senescence of primary coronary endothelial cells exposed to postprandial serum of healthy adult and elderly volunteers after oven cooked meat meals.**

Laura, C. et al. (2017) Mediators Inflamm. [Epub ahead of print].

**PubMed ID:** N/A    **Reactivity:** Human    **Applications:** HCAECs

**A multi-biomarker analysis of the antioxidant efficacy of Parkinson's disease therapy.**

Colamartino, M. et al. (2017) Toxicol In Vitro. 47:1-7.

**PubMed ID:** 29080800    **Reactivity:** Human    **Applications:** buffy coat

**Association of Metabolic Syndrome and Oxidative DNA Damage in HIV/AIDS Patients.**

Kolgiri, V., Nagar, V. and Patil, V. (2017) Ind J Clin Biochem. 1669:18-26.

**PubMed ID:** 28535982    **Reactivity:** Human    **Applications:** Plasma

**Antiretroviral Therapy-induced Insulin Resistance and Oxidative Deoxy Nucleic Acid Damage in Human Immunodeficiency Virus-1 Patients.**

VK Honnapurmath, V.W.P. (2017) Indian J Endocr Metab. 21(2):316-321.

**PubMed ID:** 28459032    **Reactivity:** Human    **Applications:** Blood

**Urinary excretion of uric acid, allantoin, and 8-OH-Deoxyguanosine in uricase-knockout mice.**

Inazawa, K. et al. (2016) Nucleosides Nucleotides Nucleic Acids. 35(10-12):559-565.

**PubMed ID:** 27906613    **Reactivity:** Mouse    **Applications:** Urine

**Protein carbonyl content: a novel biomarker for aging in HIV/AIDS patients.**

Kolgiri, V. and Patil, V.W. (2016) Braz J Infect Dis. pii: S1413-8670(16)30432-9.

**PubMed ID:** 27821249    **Reactivity:** Human    **Applications:** Plasma

**Coagulin-I ameliorates TLR4 induced oxidative damage and immune response by regulating mitochondria and NOX-derived ROS.**

Reddy, SS (2016) Toxicol Appl Pharmacol. 309:87-100.

**PubMed ID:** 27568862    **Reactivity:** Mouse    **Applications:** Liver homogenates

**Preventive Effect of Resveratrol against Brain Mitochondria DNA Damage, Lipid Peroxidation, Inflammation and Seizures Induced by Kainic Acid in Mice.**

Hussein, S. A., Abdel-mageid, A.D., Abd-Elhamed, O. M., Amin, A. and Al harthy, H.S. (2016) Int J Phar Sci. 6(4): 1634-1646.

**PubMed ID:** N/A    **Reactivity:** Mouse    **Applications:** Brain tissue

**Neuroprotective Effect of Curcumin on Kainic Acid Model of Epilepsy in Male Swiss**

#### **Albino Mice.**

Hussein, S. A., Abdel-mageid, A.D., Abd-Elhamed, O. M., Amin, A. and Al harthy, H.S. (2016) Int J Chem Nat Sci. 4(4): 447-460.

**PubMed ID:** N/A    **Reactivity:** Mouse    **Applications:** Brain tissue

#### **DNA damage and the activation of the p53 pathway mediate alterations in metabolic and secretory functions of adipocytes.**

Vergoni, B. et al. (2016) Diabetes.

**PubMed ID:** 27388216    **Reactivity:** Human    **Applications:** Adipocytes

#### **Discovering biomarkers for antidepressant response: protocol from the Canadian biomarker integration network in depression (CAN-BIND) and clinical characteristics of the first patient cohort.**

Lam, R.W. et al. (2016) BMC Psychiatry. 16(1):105.

**PubMed ID:** 27084692    **Reactivity:** Human    **Applications:** Blood and Urine

#### **Podocyte-specific Nox4 deletion affords renoprotection in a mouse model of diabetic nephropathy.**

Jha, J.C. et al. (2015) Diabetologia. 59(2):379-89.

**PubMed ID:** 26508318    **Reactivity:** Mouse    **Applications:** Urine

#### **Does chronic raise of metal ion levels induce oxidative DNA damage and hypoxia-like response in patients with metal-on-metal hip resurfacing?**

Savarino, L. et al. (2015) J Biomed Mater Res B Appl Biomater. 105(2):460-466.

**PubMed ID:** 26477446    **Reactivity:** Human    **Applications:** Serum

#### **Higher urinary Levels of 8-hydroxy-2'-deoxyguanosine are associated with a worse RANKL/OPG ratio in postmenopausal women with osteopenia.**

Cervellati, C. et al. (2015) Oxid Med Cell Long. 2016:6038798.

**PubMed ID:** 26635910    **Reactivity:** Human    **Applications:** Urine

#### **Turmerone enriched standardized Curcuma longa extract alleviates LPS induced inflammation and cytokine production by regulating TLR4-IRAK1-ROS-MAPK-NFκB axis.**

Rana, M. et al. (2015) J Functional Foods. 16 (2015): 152-163.

**PubMed ID:** N/A    **Reactivity:** Mouse    **Applications:** Liver homogenates

#### **Synthesis, characterization and biological activity of some unsymmetrical Schiff base transition metal complexes.**

Esmadi, F.T. et al. (2015) Drug Chem Toxicol. :1-7.

**PubMed ID:** 25791998    **Reactivity:** Human    **Applications:** Lymphocytes

**Late-intervention study with ebselen in an experimental model of type 1 diabetic nephropathy.**

Tan, S.M., Sharma, A., Stefanovic, N. and de Haan, J.B. (2015) Free Radic Res. 49(3):219-27.

**PubMed ID:** 25465090    **Reactivity:** Mouse    **Applications:** Urine

**Immediate and delayed effects of growth conditions on ageing parameters in nestling zebra finches.**

Reichert, S. et al. (2015) J Exp Biol. 218(Pt 3):491-9.

**PubMed ID:** 25524985    **Reactivity:** Zebra finch    **Applications:** Blood cells

**A derivative of Bardoxolone methyl, dh404, in an inverse dose-dependent manner, lessens diabetes-associated atherosclerosis and improves diabetic kidney disease.**

Tan, S. M., et al. (2014) Diabetes. 63(9):3091-103.

**PubMed ID:** 24740568    **Reactivity:** Mice    **Applications:** Urine

**Elevation impacts the balance between growth and oxidative stress in coal tits.**

Stier, A., et al. (2014) Oecologia. 175(3):791-800.

**PubMed ID:** 24805201    **Reactivity:** Bird (coal tits)    **Applications:** Red blood cells

**Evaluation of vitamin B12 effects on DNA damage induced by paclitaxel.**

Alzoubi, K., Khabour O., Khader M., Mhaidat N. and Al-Azzam, S. (2014) Drug Chem Toxicol. 37(3):276-80.

**PubMed ID:** 24215581    **Reactivity:** Human    **Applications:** Blood

**Effect of Therapeutic Hypothermia on DNA Damage and Neurodevelopmental Outcome Among Term Neonates with Perinatal Asphyxia: A Randomized Controlled Trial.**

Gane, B.D., et al. (2014) J Trop Pediatr. 60(2):134-40.

**PubMed ID:** 24343823    **Reactivity:** Human    **Applications:** Blood

**Disruption of pro-oxidant and antioxidant systems with elevated expression of the ubiquitin proteasome system in the cachectic heart muscle of nude mice.**

Hinch, E. C. A., Sullivan-Gunn, M. J., Vaughan V. C., McGlynn M. A., Lewandowski, Paul A. (2013) J Cachexia Sarcopenia Muscle. 4(4):287-93.

**PubMed ID:** 24030522    **Reactivity:** Mouse    **Applications:** Tumor tissue

**Oxidative Damage of DNA Confers Resistance to Cytosolic Nuclease TREX1 Degradation and Potentiates STING-Dependent Immune Sensing.**

Gehrke, N. et al. (2013) Immunity. 39(3):482-95.

**PubMed ID:** 23993650    **Reactivity:** Mouse    **Applications:** RMA cells (T-cell tumor cells)

**Differential effects of docosahexanoic acid (DHA) on preterm and term placental pro-oxidant/anti-oxidant balance.**

Stark, M. et al. (2013) Reproduction. 146(3):243-51.

**PubMed ID:** 23813449    **Reactivity:** Human    **Applications:** Placenta

**The Modified Selenenyl Amide, M-hydroxy Ebselen, Attenuates Diabetic Nephropathy and Diabetes-Associated Atherosclerosis in ApoE/GPx1 Double Knockout Mice.**

Tan, S.M., et al. (2013) PLoS ONE. 8(7): e69193.

**PubMed ID:** 23874911    **Reactivity:** Mouse    **Applications:** Urine (diabetic and non-diabetic)

**TiO<sub>2</sub> Nanoparticle Exposure and Illumination during Zebrafish Development: Mortality at Parts per Billion Concentrations.**

Bar-Ilan, O. et al. (2013) Environ Sci Technol. 47 (9): 4726-4733.

**PubMed ID:** 23510150    **Reactivity:** Zebrafish    **Applications:** DNA extract

**Poor sleep in PCOS; is melatonin the culprit?**

Shreeve, N. et al. (2013) Hum Reprod. 28 (5): 1348-1353.

**PubMed ID:** 23438443    **Reactivity:** Human    **Applications:** Urine

**Grasshoppers' adaptation to elevated radioactivity in the Chernobyl exclusion zone.**

Mortensen, L.H. (2013) Roskilde University Dissertation.

**PubMed ID:** N/A    **Reactivity:** Grasshopper    **Applications:** Intestines

**Consumption of a low glycaemic index diet in late life extends lifespan of Balb/c mice with differential effects on DNA damage.**

Nankervis, S.A., Mitchell, J.M., Charchar, F.J., McGlynn, M.A. and Lewandowski, P.A. (2013) Longev Healthspan. 2(1):4.

**PubMed ID:** 24472560    **Reactivity:** Mouse    **Applications:** WBC's

**Light-induced cell detachment for cell sheet technology.**

Hong, Y. et al. (2013) Biomaterials. 34 (1): 11-18.

**PubMed ID:** 23069710    **Reactivity:** Mouse    **Applications:** MC3T3-E1 cells

**Puerarin prevents high glucose-induced apoptosis of Schwann cells by inhibiting oxidative stress.**

Wu, Y., Xue, B., Li, X. and Liu, H. (2012) Neural Regen Res. 7 (33): 2583-2591.

**PubMed ID:** 25368634    **Reactivity:** Mouse    **Applications:** Schwann cells

**The protective effect of Alpha lipoic acid on Schwann cells exposed to constant or intermittent high glucose.**

Sun, L. et al. (2012) Biochem Pharmacol. 84 (7): 961-973.

**PubMed ID:** 22796564    **Reactivity:** Mouse    **Applications:** Schwann cells

**Ginsenoside Rb1 relieves glucose fluctuation induced oxidative stress and apoptosis in Schwann cells.**

Xue, B. et al. (2012) Neural Regen Res. 2012 (30): 2340-2346.

**PubMed ID:** 25538758    **Reactivity:** Mouse    **Applications:** Schwann cells

**Eicosapentaenoic Acid and Oxypurinol in the Treatment of Muscle Wasting in a Mouse Model of Cancer Cachexia.**

Vaughan, V.C., Sullivan-Gunn, M., Hinch, E., Martin, P., Lewandowski, P.A. (2012) PLoS ONE. 7(9): e45900.

**PubMed ID:** 23029301    **Reactivity:** Mouse    **Applications:** Gastrocnemius muscle

**Titanium dioxide nanoparticles produce phototoxicity in the developing zebrafish.**

Bar-Ilan, O. et al. (2012) Nanotoxicology. 6 (6): 670-679.

**PubMed ID:** 21830861    **Reactivity:** Zebrafish    **Applications:** Embryos

**Evaluation of vitamin B12 effects on DNA damage induced by pioglitazone.**

Alzoubi, K., Khabour, O., Hussain, N., Al-azzam, S. and Mhaidat, N. (2012) Mutat Res. 748 (1-2): 48-51.

**PubMed ID:** 22790087    **Reactivity:** Human    **Applications:** Blood

**Protective Effects of Salvianolic Acid B on Schwann Cells Apoptosis Induced by High Glucose.**

Sun, L. et al. (2012) Neurochemical Res. 37 (5): 996-1010.

**PubMed ID:** 22252725    **Reactivity:** Mouse    **Applications:** Schwann cells

**Oxidative Stress in HPV-Driven Viral Carcinogenesis: Redox Proteomics Analysis of**

### **HPV-16 Dysplastic and Neoplastic Tissues.**

De Marco, F. et al. (2012) PLoS ONE. 7(3): e34366.

**PubMed ID:** 22470562    **Reactivity:** Human    **Applications:** Ecto-cervical tissue

### **Inhibitory effects of Salvianolic acid B on apoptosis of Schwann cells and its mechanism induced by intermittent high glucose.**

Sun, L. et al. (2012) Life Sci. 90 (3-4): 99-108.

**PubMed ID:** 22036624    **Reactivity:** Mouse    **Applications:** Schwann cells

### **A qPCR-based assay to quantify oxidized guanine and other FPG-sensitive base lesions within telomeric DNA.**

O'Callaghan, N., Baack, N., Sharif, R. and French, M. (2011) Biotechniques. 51 (6): 403-411.

**PubMed ID:** 22150331    **Reactivity:** Human    **Applications:** Spleen B lymphoblastoid cells

### **The effect of cocoa supplementation on hepatic steatosis, reactive oxygen species and LFABP in a rat model of NASH.**

Janevski, M., Antonas, K.N., Sullivan-Gunn, M.J., McGlynn, M.A. and Lewandowski, P.A. (2011) Comp Hepatol. 10(1):10.

**PubMed ID:** 22081873    **Reactivity:** Rat    **Applications:** Liver tissue

### **Reduction of Oxidative Damage Reflects a Better Kidney Transplantation Outcome.**

La Manna, G. et al. (2011) Am J Nephrol. 34, 496-504.

**PubMed ID:** 22041478    **Reactivity:** Human    **Applications:** PBMC

### **Integrated exposure assessment of sewage workers to genotoxicants: an urinary biomarker approach and oxidative stress evaluation.**

Zabadi, H.A. et al. (2011) Environ Health. 10 (23).

**PubMed ID:** 21435260    **Reactivity:** Human    **Applications:** Urine

### **Transgenic Mice Expressing Cyclooxygenase-2 in Hepatocytes Reveal a Minor Contribution of This Enzyme to Chemical Hepatocarcinogenesis.**

Izquierdo, C.L. et al. (2011) Am J Pathol. 178 (3): 1361-1373.

**PubMed ID:** 21356386    **Reactivity:** Mouse    **Applications:** COX-2 Tg cells

## **Reviews**

---

Based on validation through cited publications.

SSSSS

**StressMarq Biosciences**

January 17, 2017: