



ELISA and Assay Kits




Western Blot



Sandwich ELISA




Immunohistochemistry




Immunofluorescence



Immunoprecipitation



Flow Cytometry



ELISA Kits

1

- ELISA Kits
- Disease ELISA Kits
- Immunoglobulin ELISA Kits
- Environmental Risk ELISA Kits
- Food Safety ELISA Kits

Assay Kits


21

- Activity Assay Kits
- Apoptosis Assay Kits
- Cell Damage and Oxidative Stress Assay Kits
- Cell Proliferation and Cytotoxicity Assay Kits
- Epigenetic Assay Kits
- Fluorescent Protein Assay Kits
- Metabolism Assay Kits
- Signal Transduction Assay Kits
- Others Assay Kits


Other Kits

28


- Ab Concentration Kits
- Ab Labeling Kits
- Ab Purification Kits
- IHC Kits
- IP & ChIP Kits
- Protein Expression Kits




KA0432




KA0202




KA0457



KA1071



KA0502




KA1252

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA0486	A2M (Human) ELISA Kit	0.002 to 0.5 ug/mL	2 ng/mL
KA0487	A2M (Human) ELISA Kit	1.25 to 40 ug/mL	1 ug/mL
KA1043	A2m (Mouse) ELISA Kit	2 to 500 ng/mL	0.1 ng/mL
KA1044	A2m (Mouse) ELISA Kit	1.56 to 100 ug/mL	1.5 ug/mL
KA2051	A2m (Rabbit) ELISA Kit	6.25 to 400 ng/mL	N/A
KA1931	A2m (Rat) ELISA Kit	6.25 to 400 ng/mL	N/A
KA1126	ABCA1 (Human) ELISA Kit	78 to 5000 pg/mL	39 pg/mL
KA1158	Abeta1-40 (Human) ELISA Kit	0.78 to 50 nmol/mL	0.25 nmol/mL
KA1174	Abeta1-42 (Human) ELISA Kit	15.6 to 1,000 pg/mL	7.8 pg/mL
KA0432	ACE (Human) ELISA Kit	156 to 10000 pg/mL	5 pg/mL
KA1133	ACHE Ab (Human) ELISA Kit	156 to 10000 pg/mL	N/A
KA0041	ACP5 (Human) ELISA Kit	0.2 to 4 U/L	0.4 U/L
KA1130	ACTA2 (Human) ELISA Kit	1.56 to 100 ng/mL	0.7 ng/mL
KA0916	ACTH (Human) ELISA Kit	5 to 165 pg/mL	0.46 pg/mL
KA0917	ACTH (Mouse/Rat) ELISA Kit	5 to 700 pg/mL	1 pg/mL
KA1235	Activated PROC (Human) ELISA Kit	156 to 10000 pg/mL	39 pg/mL
KA1239	ADAM10 (Human) ELISA Kit	0.156 to 10 ng/mL	0.039 ng/mL
KA1128	ADMA (Human) ELISA Kit	0.312 to 20 nmol/mL	0.078 nmol/mL
KA0202	AFP (Human) ELISA Kit	0 to 300 ng/mL	2 ng/mL
KA1024	AFP (Human) ELISA Kit	0.156 to 10 ng/mL	0.05 ng/mL
KA1187	AFP-L3 (Human) ELISA Kit	0.312 to 20 ng/mL	0.078 ng/mL
KA1764	AGER (Human) ELISA Kit	N/A	3 pg/mL
KA1835	Ager (Rat) ELISA Kit	N/A	300 pg/mL
KA0457	AGT (Human) ELISA Kit	0.063 to 1 ng/mL	0.03 ng/mL
KA1679	AGT (Human/Mouse/Rat) ELISA Kit	1 to 100 ng/mL	49.29 pM
KA0478	AHSG (Human) ELISA Kit	6.25 to 400 ng/mL	5 ng/mL
KA1071	AHSG (Human) ELISA Kit	0.312 to 20 ng/mL	< 10 pg/mL
KA0454	Albumin (Human) ELISA Kit	1.563 to 100 ug/mL	1.5 ug/mL
KA0455	Albumin (Human) ELISA Kit	3.125 to 200 ng/mL	3 ng/mL
KA0488	Albumin (Mouse) ELISA Kit	0.098 to 25 ug/mL	0.3 ug/mL
KA0489	Albumin (Mouse) ELISA Kit	1.56 to 400 ng/mL	1.5 ng/mL
KA0497	Albumin (Pig) ELISA Kit	1.56 to 100 ng/mL	1 ng/mL
KA1027	Albumin (Pig) ELISA Kit	1.56 to 100 ug/mL	0.3 ug/mL
KA0501	Albumin (Rat) ELISA Kit	0.391 to 100 ug/mL	0.3 ug/mL
KA0502	Albumin (Rat) ELISA Kit	0.781 to 200 ng/mL	0.7 ng/mL
KA1692	ALCAM (Human) ELISA Kit	N/A	15 pg/mL
KA1224	ALDH (Human) ELISA Kit	1.56 to 100 U/L	0.78 U/L
KA1883	Aldosterone ELISA Kit	N/A	10 pg/mL
KA1965	Alpha 1-Anti-Chymotrypsin (Human) ELISA Kit	6.25 to 400 ng/mL	N/A
KA1406	Alpha1-Microglobulin (Human) ELISA Kit	1 to 10 ng/mL	0.6 ng/mL
KA1209	alpha-GST (Human) ELISA Kit	0.78 to 50 ng/mL	N/A
KA1190	ALPP (Human) ELISA Kit	15.6 to 1,000 pg/mL	3.9 pg/mL
KA1132	Amiodarone (Human) ELISA Kit	1.56 to 100 ng/mL	0.39 ng/mL
KA1136	AMS/AMY (Human) ELISA Kit	15.6 to 1000 U/L	6.2 U/L
KA1898	Androstenedione ELISA Kit	N/A	0.04 ng/mL
KA1252	ANG (Human) ELISA Kit	78 to 5000 pg/mL	12 pg/mL
KA1693	Angiostatin (Human) ELISA Kit	N/A	20 ng/mL
KA1866	ANGPT1 (Human) ELISA Kit	N/A	30 pg/mL
KA1867	ANGPT2 (Human) ELISA Kit	N/A	10 pg/mL
KA1694	ANGPTL4 (Human) ELISA Kit	N/A	20 pg/mL
KA0150	Annexin V Ab (Human) ELISA Kit	6.3 to 400 ng/mL	1.2 ng/mL
KA0158	AOC3 (Human) ELISA Kit	31.3 to 2000 pg/mL	19 pg/mL
KA1841	APCS (Human) ELISA Kit	N/A	0.15 ng/mL
KA1976	Apcs (Mouse) ELISA Kit	12.5 to 200 ng/mL	N/A
KA1934	Apcs (Rat) ELISA Kit	2.5 to 160 ng/mL	N/A
KA1681	APLN (Human/Mouse/Rat) ELISA Kit	0.1 to 1000 ng/mL	29.1 pg/mL
KA0460	APOA1 (Human) ELISA Kit	0.625 to 20 ug/mL	0.6 ug/mL
KA0462	APOA1 (Human) ELISA Kit	3.13 to 200 ug/mL	3 ng/mL
KA0461	APOA2 (Human) ELISA Kit	0.004 to 4 ug/mL	0.004 ug/mL
KA1028	APOB (Human) ELISA Kit	0.008 to 0.250 ug/mL	0.008 ug/mL
KA0463	APOC1 (Human) ELISA Kit	0.063 to 4 ug/mL	0.06 ug/mL
KA1997	APOC2 (Human) ELISA Kit	N/A	0.15 ug/mL
KA0465	APOC3 (Human) ELISA Kit	0.002 to 2 ug/mL	0.005 ug/mL
KA1398	APOC3 (Monkey) ELISA Kit	0.078 to 5 ug/mL	0.1 ug/mL
KA1030	Apoc3 (Rat/Mouse) ELISA Kit	0.5 to 8 ug/mL	0.5 ug/mL

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA1031	APOE (Human) ELISA Kit	31 to 2000 ng/mL	30 ng/mL
KA0982	APOH (Human) ELISA Kit	0.625 to 40 ng/mL	0.5 ng/mL
KA1201	AQP1 (Human) ELISA Kit	0.156 to 10 ng/mL	0.04 ng/mL
KA1202	AQP2 (Human) ELISA Kit	0.156 to 10 ng/mL	0.039 ng/mL
KA1127	AR (Human) ELISA Kit	0.312 to 20 ng/mL	0.11 ng/mL
KA1695	AREG (Human) ELISA Kit	N/A	10 pg/mL
KA1775	Areg (Mouse) ELISA Kit	N/A	2 pg/mL
KA1895	ARG1 (Human) ELISA Kit	N/A	0.5 ng/mL
KA0301	Arg8-Vasopressin ELISA Kit	4.10 to 1000 pg/mL	3.39 pg/mL
KA1141	ARRB2 (Human) ELISA Kit	0.156 to 10 ng/mL	N/A
KA1696	AXL (Human) ELISA Kit	N/A	2 pg/mL
KA1776	Axl (Mouse) ELISA Kit	N/A	1 pg/mL
KA1689	AZGP1 (Human/Mouse/Rat) ELISA Kit	0.1 to 1000 ng/mL	21 pg/mL
KA0222	B2M (Human) ELISA Kit	0 to 10 ug/mL	0.1 ug/mL
KA1045	B2M (Human) ELISA Kit	0.78 to 200 ng/mL	0.7 ng/mL
KA2021	B2m (Rat) ELISA Kit	25 to 800 ng/mL	N/A
KA1231	BACE1 (Human) ELISA Kit	1.56 to 100 U/L	0.39 U/L
KA1184	BALP (Human) ELISA Kit	6.25 to 400 U/L	3.12 U/L
KA0147	BCL2 (Human) ELISA Kit	0.50 to 32 ng/mL	0.5 ng/mL
KA0329	BDNF (Human) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA0331	Bdnf (Mouse) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA0330	Bdnf (Rat) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA1150	Beta-endorphin (Human) ELISA Kit	15.6 to 1,000 pg/mL	3.9 pg/mL
KA0441	BIRC5 (Human) ELISA Kit	62.5 to 4000 pg/mL	2 pg/mL
KA0540	BMP2 (Human) ELISA Kit	62.5 to 4000 pg/mL	2 pg/mL
KA0542	Bmp2 (Mouse) ELISA Kit	62.5 to 4000 pg/mL	2 pg/mL
KA0541	Bmp2 (Rat) ELISA Kit	62.5 to 4000 pg/mL	2 pg/mL
KA0543	BMP4 (Human) ELISA Kit	62.5 to 4000 pg/mL	4 pg/mL
KA1501	BMP5 (Human) ELISA Kit	62.5 to 4000 pg/mL	N/A
KA1698	BMP6 (Human) ELISA Kit	N/A	150 pg/mL
KA1699	BMP7 (Human) ELISA Kit	N/A	10 pg/mL
KA0988	BSG (Human) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA1006	BSG (Human) ELISA Kit	0.78 to 50 ng/mL	0.30 ng/mL
KA1700	BTC (Human) ELISA Kit	N/A	30 pg/mL
KA1186	C1QTNF3 (Human) ELISA Kit	4.68 to 300 ng/mL	1.2 ng/mL
KA1923	C3 (Dog) ELISA Kit	6.25 to 400 ng/mL	N/A
KA1925	C3 (Guinea Pig) ELISA Kit	3.125 to 200 ng/mL	N/A
KA0466	C3 (Human) ELISA Kit	0.47 to 30 ug/mL	0.5 ug/mL
KA1020	C3 (Human) ELISA Kit	0.625 to 40 ng/mL	0.6 ng/mL
KA1926	C3 (Mouse) ELISA Kit	3.125 to 200 ng/mL	N/A
KA1932	C3 (Rat) ELISA Kit	12.5 to 800 ng/mL	N/A
KA1021	C4 (Human) ELISA Kit	N/A	0.08 ng/mL
KA1159	CA (Human) ELISA Kit	1.56 to 100 IU/mL	0.8 IU/mL
KA0205	CA125 (Human) ELISA Kit	0 to 400 U/mL	5 U/mL
KA0206	CA15-3 (Human) ELISA Kit	0 to 240 U/mL	5 U/mL
KA0207	CA19-9 (Human) ELISA Kit	0 to 600 U/mL	10 U/mL
KA1226	CALB2 (Human) ELISA Kit	0.78 to 50 ng/mL	0.39 ng/mL
KA1682	CARTPT (Human/Mouse/Rat) ELISA Kit	1 to 100 ng/mL	497.87 pM
KA1955	Casein (Bovine) ELISA Kit	3.125 to 200 ng/mL	N/A
KA0184	CASP8 (Human) ELISA Kit	0.16 to 10 ng/mL	0.1 ng/mL
KA0185	CASP9 (Human) ELISA Kit	1.6 to 100 ng/mL	0.4 ng/mL
KA1862	CCK (Human/Mouse/Rat) ELISA Kit	0.1 to 1000 pg/mL	3.86 pg/mL
KA1725	CCL1 (Human) ELISA Kit	N/A	4 pg/mL
KA0334	CCL11 (Human) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA0240	Ccl11 (Mouse) ELISA Kit	31.3 to 2000 pg/mL	3.6 pg/mL
KA1801	Ccl12 (Mouse) ELISA Kit	N/A	0.5 pg/mL
KA1745	CCL13 (Human) ELISA Kit	N/A	3 pg/mL
KA1748	CCL15 (Human) ELISA Kit	N/A	3 pg/mL
KA1724	CCL16 (Human) ELISA Kit	N/A	8 pg/mL
KA1506	CCL17 (Human) ELISA Kit	31.2 to 2000 pg/mL	N/A
KA1757	CCL18 (Human) ELISA Kit	N/A	2 pg/mL
KA1749	CCL19 (Human) ELISA Kit	N/A	2 pg/mL
KA1806	Ccl19 (Mouse) ELISA Kit	N/A	5 pg/mL
KA0169	CCL2 (Human) ELISA Kit	15.6 to 1000 pg/mL	2.3 pg/mL
KA0255	Ccl2 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	2.2 pg/mL




KA0441




KA0540




KA0543




KA1501




KA0206




KA1506




KA0389




KA0386




KA0388



KA1073



KA1069




KA0668

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA1831	Ccl2 (Rat) ELISA Kit	N/A	15 pg/mL
KA1832	Ccl2 (Rat) ELISA Kit	N/A	15 pg/mL
KA0389	CCL20 (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1805	Ccl20 (Mouse) ELISA Kit	N/A	3 pg/mL
KA2000	Ccl20 (Rat) ELISA Kit	15.6 to 1000 pg/mL	< 1 pg/mL
KA1774	Ccl21a (Mouse) ELISA Kit	N/A	20 pg/mL
KA0386	CCL22 (Human) ELISA Kit	31.2 to 2000 pg/mL	1 pg/mL
KA0387	Ccl22 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1751	CCL23 (Human) ELISA Kit	N/A	7 pg/mL
KA1709	CCL24 (Human) ELISA Kit	N/A	2 pg/mL
KA1782	Ccl24 (Mouse) ELISA Kit	N/A	4 pg/mL
KA1766	CCL25 (Human) ELISA Kit	N/A	200 pg/mL
KA1813	Ccl25 (Mouse) ELISA Kit	N/A	15 pg/mL
KA1781	Ccl27a (Mouse) ELISA Kit	N/A	15 pg/mL
KA1701	CCL28 (Human) ELISA Kit	N/A	25 pg/mL
KA0388	CCL3 (Human) ELISA Kit	31.2 to 2000 pg/mL	10 pg/mL
KA1804	Ccl3 (Mouse) ELISA Kit	N/A	8 pg/mL
KA1747	CCL4 (Human) ELISA Kit	N/A	2.5 pg/mL
KA0411	CCL5 (Human) ELISA Kit	31.2 to 2000 pg/mL	3 pg/mL
KA2001	Ccl5 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	< 1 pg/mL
KA0519	CCL7 (Human) ELISA Kit	30 to 1600 pg/mL	30 pg/mL
KA1744	CCL8 (Human) ELISA Kit	N/A	1.5 pg/mL
KA0446	CD14 (Human) ELISA Kit	5 to 50 ng/mL	5 ng/mL
KA0447	Cd14 (Mouse) ELISA Kit	5 to 50 ng/mL	5 ng/mL
KA0172	CD28 (Human) ELISA Kit	0.8 to 50 ng/mL	0.18 ng/mL
KA1780	Cd36 (Mouse) ELISA Kit	N/A	30 pg/mL
KA0161	CD40 (Human) ELISA Kit	7.8 to 500 pg/mL	6.92 pg/mL
KA1061	CD40 (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1062	Cd40 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0144	CD40LG (Human) ELISA Kit	0.16 to 10 ng/mL	0.062 ng/mL
KA0551	CD40LG (Human) ELISA Kit	62.5 to 4000 pg/mL	6 pg/mL
KA0239	CD40lg (Mouse) ELISA Kit	0.31 to 20 ng/mL	0.14 ng/mL
KA0119	CD44std (Human) ELISA Kit	0.12 to 4 ng/mL	0.015 ng/mL
KA0129	CD44var (v5) (Human) ELISA Kit	0.32 to 10 ng/mL	0.14 ng/mL
KA0120	CD44var (v6) (Human) ELISA Kit	0.32 to 10 ng/mL	0.05 ng/mL
KA1073	CD5L (Human) ELISA Kit	1.57 to 100 ng/mL	0.745 ng/mL
KA0133	CD8 (Human) ELISA Kit	31.3 to 2000 U/mL	6.9 U/mL
KA0991	Cd80 (Mouse) ELISA Kit	62.5 to 4000 pg/mL	10 pg/mL
KA0433	CDH1 (Human) ELISA Kit	156 to 10000 pg/mL	15 pg/mL
KA1069	Cdh1 (Mouse) ELISA Kit	156 to 10000 pg/mL	15 pg/mL
KA1135	CDH2 (Human) ELISA Kit	312 to 20000 pg/mL	78 pg/mL
KA0668	CDH3 (Human) ELISA Kit	62.5 to 4000 pg/mL	2 pg/mL
KA0152	CDH5 (Human) ELISA Kit	0.16 to 10 ng/mL	0.15 ng/mL
KA0271	CDH5 (Monkey) ELISA Kit	0.19 to 12.5 U/mL	0.1 U/mL
KA0204	CEA (Human) ELISA Kit	0 to 120 ng/mL	1 ng/mL
KA1702	CEACAM1 (Human) ELISA Kit	N/A	20 pg/mL
KA1152	CETP (Human) ELISA Kit	3.12 to 200 ng/mL	1.56 ng/mL
KA1691	CFD (Human) ELISA Kit	N/A	4 pg/mL
KA0212	CG (Human) ELISA Kit	0 to 300 mIU/mL	2 mIU/mL
KA1160	CGRP (Human) ELISA Kit	15.6 to 1,000 pg/mL	3.9 pg/mL
KA1884	CHGA (Human) ELISA Kit	N/A	10 ng/mL
KA2073	CKM (Pig) ELISA Kit	3.125 to 50 ng/mL	N/A
KA2072	Ckm (Rat) ELISA Kit	3.125 to 50 ng/mL	N/A
KA1188	CLU (Human) ELISA Kit	1.56 to 100 ng/mL	0.39 ng/mL
KA1985	Clu (Mouse) ELISA Kit	3.9 to 250 ng/mL	N/A
KA0042	Clu (Rat) ELISA Kit	4 to 128 ng/mL	0.7 ng/mL
KA0096	CML/N ^ε (carboxymethyl) lysine ELISA Kit	80 to 5000 ng/mL	0.063 ug/mL
KA1703	CNTF (Human) ELISA Kit	N/A	0.2 ng/mL
KA1824	Cntf (Rat) ELISA Kit	N/A	10 pg/mL
KA1708	COL18A1 (Human) ELISA Kit	N/A	10 pg/mL
KA2100	COLL2-1 ELISA Kit	80 to 2000 pM	40 pM
KA2098	COLL2-1 NO2 ELISA Kit	100 to 1000 pM	15 pM
KA0021	COMP (Human) ELISA Kit	4 to 128 ng/mL	0.5 ng/mL
KA0467	Complement C4 (Human) ELISA Kit	0.313 to 5 ug/mL	0.05 ug/mL
KA0308	Complement C4a des Arg (Human) ELISA Kit	0.78 to 200 ng/mL	0.76 ng/mL

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA0468	Corticosterone ELISA Kit	0.391 to 100 ng/mL	0.3 ng/mL
KA0918	Cortisol (Human) ELISA Kit	20 to 800 ng/mL	0.44 ng/mL
KA0919	Cortisol (Sheep) ELISA Kit	1 to 60 ng/mL	0.1 ng/mL
KA1885	Cortisol ELISA Kit	0.66 to 20 ng/mL	0.3 ng/mL
KA1901	Cortisol ELISA Kit	N/A	2.0 ng/mL
KA0469	CP (Human) ELISA Kit	0.63 to 40 ug/mL	0.1 ug/mL
KA0470	CP (Human) ELISA Kit	1.25 to 320 ng/mL	1 ng/mL
KA1033	Cp (Rat) ELISA Kit	0.78 to 50 ug/mL	0.5 ug/mL
KA1034	Cp (Rat) ELISA Kit	0.31 to 20 ng/mL	0.3 ng/mL
KA1004	Creatinine ELISA Kit	1.25 to 20 mg/dl	0.625 mg/dl
KA2064	Crp (Cat) ELISA Kit	3.125 to 200 ng/mL	N/A
KA1973	CRP (Dog) ELISA Kit	3.125 to 200 ng/mL	N/A
KA2065	CRP (Horse) ELISA Kit	10 to 320 ng/mL	N/A
KA0238	CRP (Human) ELISA Kit	0 to 0.10 mg/L	0.1 mg/L
KA1052	CRP (Human) ELISA Kit	59.4 to 3800 pg/mL	28.6 pg/mL
KA1983	CRP (Monkey) ELISA Kit	1.17 to 75 ng/mL	N/A
KA1974	Crp (Mouse) ELISA Kit	0.39 to 25 ng/mL	N/A
KA1920	CRP (Pig) ELISA Kit	6.25 to 100 ng/mL	N/A
KA1975	CRP (Rabbit) ELISA Kit	6.25 to 200 ng/mL	N/A
KA1035	Crp (Rat) ELISA Kit	0.78 to 50 ng/mL	0.5 ng/mL
KA1036	Crp (Rat) ELISA Kit	0.078 to 20 ug/mL	0.07 ug/mL
KA1233	CS (Human) ELISA Kit	0.78 to 50 ng/mL	0.195 ng/mL
KA0384	CSF1 (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0385	Csf1 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0170	CSF2 (Human) ELISA Kit	7.8 to 500 pg/mL	2.9 pg/mL
KA0241	Csf2 (Mouse) ELISA Kit	3.9 to 250 pg/mL	1.98 pg/mL
KA0272	Csf2 (Rat) ELISA Kit	7.8 to 500 pg/mL	3.5 pg/mL
KA0985	Csf2 (Rat) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1718	CSF3 (Human) ELISA Kit	N/A	1 pg/mL
KA0983	Csf3 (Mouse) ELISA Kit	31.2 to 2000 pg/mL	4 pg/mL
KA0022	CST3 (Human) ELISA Kit	200 to 10000 ng/mL	0.2 ng/mL
KA2026	CST3 (Monkey) ELISA Kit	0.031 to 2 ng/mL	N/A
KA1060	Cst3 (Mouse) ELISA Kit	312 to 20000 pg/mL	10 pg/mL
KA2022	Cst3 (Rat) ELISA Kit	2.63 to 20 ng/mL	N/A
KA0434	CTF1 (Human) ELISA Kit	31.2 to 2000 pg/mL	10 pg/mL
KA0165	CTLA4 (Human) ELISA Kit	0.16 to 10 ng/mL	0.13 ng/mL
KA0552	CTSB (Human) ELISA Kit	156 to 10000 pg/mL	5 pg/mL
KA0992	Ctsb (Mouse) ELISA Kit	156 to 10000 pg/mL	5 pg/mL
KA2002	CTSD (Human) ELISA Kit	156 to 10000 pg/mL	< 10 pg/mL
KA0156	CTSL1 (Human) ELISA Kit	3.1 to 50 ng/mL	1.71 ng/mL
KA1223	CTX-I (Human) ELISA Kit	0.156 to 10 ng/mL	0.039 ng/mL
KA1225	CTX-II (Human) ELISA Kit	78 to 5000 pg/mL	19.5 pg/mL
KA0131	Cu/ZnSOD (Human) ELISA Kit	0.08 to 5 ng/mL	0.04 ng/mL
KA0545	CX3CL1 (Human) ELISA Kit	156 to 10000 pg/mL	50 pg/mL
KA0663	Cx3cl1 (Mouse) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA1826	Cx3cl1 (Rat) ELISA Kit	N/A	5 pg/mL
KA1827	Cx3cl1 (Rat) ELISA Kit	N/A	5 pg/mL
KA1722	CXCL1 (Human) ELISA Kit	N/A	50 pg/mL
KA0553	Cxcl1 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0554	Cxcl1 (Rat) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA2004	CXCL10 (Human) ELISA Kit	31.2 to 2000 pg/mL	< 1 pg/mL
KA0259	Cxcl10 (Mouse) ELISA Kit	7.8 to 500 pg/mL	6.5 pg/mL
KA1739	CXCL11 (Human) ELISA Kit	N/A	2 pg/mL
KA1798	Cxcl11 (Mouse) ELISA Kit	N/A	6 pg/mL
KA0664	Cxcl12 (Mouse) ELISA Kit	62.5 to 4000 pg/mL	10 pg/mL
KA1697	CXCL13 (Human) ELISA Kit	N/A	1.5 pg/mL
KA1778	Cxcl13 (Mouse) ELISA Kit	N/A	2 pg/mL
KA0556	Cxcl16 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1505	Cxcl2 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	N/A
KA1870	CXCL5 (Human) ELISA Kit	N/A	10 pg/mL
KA1799	Cxcl5 (Mouse) ELISA Kit	N/A	35 pg/mL
KA1828	Cxcl5 (Rat) ELISA Kit	N/A	20 pg/mL
KA1829	Cxcl5 (Rat) ELISA Kit	N/A	15 pg/mL
KA0555	CXCL9 (Human) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA1803	Cxcl9 (Mouse) ELISA Kit	N/A	3 pg/mL




KA0384




KA1060




KA2002



KA0545




KA0553




KA0556


Note:
 1. Assay Buffer: 0.30 pmol/mL (non-acetylated); 0.039 pmol/mL (acetylated). HCl: 0.39 pmol/mL (non-acetylated); 0.037 pmol/mL (acetylated)




KA0514




KA0544




KA0442



KA1070



KA1987



KA1391

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA0320	cyclic AMP Complete ELISA Kit	0.078 to 20 pmol/mL	0.039 pmol/mL ¹
KA0321	cyclic GMP Complete ELISA Kit	0.08 to 50 pmol/mL	0.043 pmol/mL ²
KA1207	Cyclosporin A (Human) ELISA Kit	0.312 to 20 nmol/mL	0.078 nmol/mL
KA0311	Cysteinyl Leukotriene ELISA Kit	78.1 to 2500 pg/mL	26.6 pg/mL
KA0160	Cytochrome c (Human) ELISA Kit	0.08 to 5 ng/mL	0.05 ng/mL
KA0532	DBI (Human) ELISA Kit	25 to 1600 pg/mL	25 pg/mL
KA1705	DCN (Human) ELISA Kit	N/A	1.5 pg/mL
KA0315	DHEA ELISA Kit	12.21 to 50000 pg/mL	2.9 pg/mL
KA0920	DHEA-S (Human) ELISA Kit	0.1 to 10 ug/mL	0.02 ug/mL
KA1902	DHEA-S ELISA Kit	N/A	25 pg/mL
KA1003	Diacetylspermidine ELISA Kit	18.75 to 600 nM	9.375 nM
KA1002	Diacetylspermine ELISA Kit	12.5 to 200 nM	6.25 nM
KA1886	Dihydrotestosterone ELISA Kit	N/A	6 pg/mL
KA1869	DKK1 (Human) ELISA Kit	N/A	100 pg/mL
KA1706	DKK3 (Human) ELISA Kit	N/A	3 pg/mL
KA1707	DKK4 (Human) ELISA Kit	N/A	10 pg/mL
KA1887	Dopamine ELISA Kit	114 to 1917 ng/mL ³	4.5 ng/mL ⁴
KA0141	DPP4 (Human) ELISA Kit	15.6 to 500 ng/mL	7.26 ng/mL
KA0514	EGF (Human) ELISA Kit	20 to 1000 pg/mL	20 pg/mL
KA0544	Egf (Mouse) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1842	ELA2 (Human) ELISA Kit	N/A	0.15 ng/mL
KA0442	ENG (Human) ELISA Kit	156 to 10000 pg/mL	5 pg/mL
KA1910	ENO2 (Human) ELISA Kit	N/A	4 ng/mL
KA1156	Epinephrine (Human) ELISA Kit	3.12 to 200 pg/mL	0.78 pg/mL
KA1882	Epinephrine ELISA Kit	4.6 to 81.4 ng/mL ⁵	0.3 ng/mL ⁶
KA1877	Epinephrine/Norepinephrine ELISA Kit	1.3 to 81.4 ng/mL ⁷	0.3 ng/mL urine ⁸
KA1880	Epinephrine/Norepinephrine/Dopamine ELISA Kit	1.3 to 81.4 ng/mL ⁹	0.3 ng/mL urine ¹⁰
KA0175	EPO (Human) ELISA Kit	1.6 to 100 mIU/mL	0.4 mIU/mL
KA1998	Epo (Mouse) ELISA Kit	46.9 to 3000 pg/mL	< 15 pg/mL
KA0118	ERBB2 (Human) ELISA Kit	0.16 to 10 ng/mL	0.06 ng/mL
KA1070	ERBB2 (Human) ELISA Kit	62.5 to 4000 pg/mL	10 pg/mL
KA1710	ERBB3 (Human) ELISA Kit	N/A	4 pg/mL
KA1783	Ereg (Mouse) ELISA Kit	N/A	10 pg/mL
KA0116	E-selectin (Human) ELISA Kit	0.8 to 50 ng/mL	0.33 ng/mL
KA0260	E-selectin (Monkey) ELISA Kit	1.9 to 60 U/mL	0.4 U/mL
KA0234	Estradiol (Human) ELISA Kit	0 to 1000 pg/mL	10 pg/mL
KA0297	Estradiol ELISA Kit	29.3 to 30000 pg/mL	28.5 pg/mL
KA1907	Estradiol ELISA Kit	0.8 to 93 pg/mL	0.4 pg/mL
KA0316	Estriol ELISA Kit	122 to 500000 pg/mL	59.6 pg/mL
KA1908	Estrone ELISA Kit	N/A	10.0 pg/mL
KA1153	ET (Human) ELISA Kit	0.0156 to 1 EU/mL	0.01 EU/mL
KA0473	F10 (Human) ELISA Kit	1.56 to 100 ng/mL	1 ng/mL
KA1400	F11 (Human) ELISA Kit	2 to 20 ng/mL	0.7 ng/mL
KA1037	F12 (Human) ELISA Kit	3.13 to 100 ng/mL	0.05 ng/mL
KA0474	F13 (Human) ELISA Kit	2.5 to 160 ng/mL	0.05 ng/mL
KA1407	F2 (Pig) ELISA Kit	1 to 100 ng/mL	0.03 ng/mL
KA2099	F3F ELISA Kit	12.5 to 400 pM	8 pM
KA0471	F7 (Human) ELISA Kit	6.25 to 200 ng/mL	6 ng/mL
KA1837	F7 (Human) ELISA Kit	N/A	1.4 ng/mL
KA0472	F9 (Human) ELISA Kit	1.56 to 100 ng/mL	1.5 ng/mL
KA2020	FABP1 (Human) ELISA Kit	1 to 32 ng/mL	N/A
KA1389	Fabp1 (Mouse) ELISA Kit	N/A	7.49 pg/mL
KA1987	Fabp1 (Rat) ELISA Kit	15.63 to 1000 ng/mL	N/A
KA2049	Fabp1 (Rat) ELISA Kit	12.5 to 400 ng/mL	N/A
KA1990	FABP3 (Dog) ELISA Kit	0.3125 to 20 ng/mL	N/A
KA2071	FABP3 (Monkey) ELISA Kit	3.9 to 250 ng/mL	N/A
KA1390	Fabp3 (Mouse) ELISA Kit	N/A	19.51 pg/mL
KA1993	FABP3 (Pig) ELISA Kit	0.78 to 25 ng/mL	N/A
KA2060	FABP3 (Rabbit) ELISA Kit	1.56 to 100 ng/mL	N/A
KA1995	Fabp3 (Rat) ELISA Kit	0.078 to 5 ng/mL	N/A
KA1056	Fabp4 (Mouse) ELISA Kit	39 to 2500 pg/mL	34.4 pg/mL
KA1055	Fabp4 (Rat) ELISA Kit	0.156 to 10 ng/mL	163.4 pg/mL
KA1391	Fabp5 (Mouse) ELISA Kit	N/A	6.91 pg/mL
KA0148	FAS (Human) ELISA Kit	15.6 to 1000 pg/mL	13.2 pg/mL
KA0159	FASLG (Human) ELISA Kit	0.16 to 10 ng/mL	0.07 ng/mL

Note:


- Assay Buffer: 0.42 pmol/mL (non-acetylated); 0.043 pmol/mL (acetylated). HCl: 0.604 pmol/mL (non-acetylated); 0.059 pmol/mL (acetylated)
- 114 to 1917 ng/mL urine, 1000 to 8000 pg/mL plasma
- 4.5 ng/mL urine, 20 pg/mL plasma
- 4.6 to 81.4 ng/mL urine, 92 to 545 pg/mL plasma
- 0.3 ng/mL urine, 10 pg/mL plasma

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA1825	Faslg (Rat) ELISA Kit	N/A	90 pg/mL
KA0134	FCER2 (Human) ELISA Kit	12.5 to 400 U/mL	6.8 U/mL
KA1941	Ferritin (Mouse) ELISA Kit	12.5 to 400 ng/mL	N/A
KA1949	Ferritin (Rat) ELISA Kit	12.5 to 400 ng/mL	N/A
KA1129	FGB (Human) ELISA Kit	0.156 to 10 nmol/L	0.078 nmol/L
KA0527	FGF2 (Human) ELISA Kit	78 to 8000 pg/mL	78 pg/mL
KA1777	Fgf2 (Mouse) ELISA Kit	N/A	8 pg/mL
KA1849	FGF21 (Human) ELISA Kit	N/A	0.03 ng/mL
KA1711	FGF4 (Human) ELISA Kit	N/A	50 pg/mL
KA1712	FGF4 (Human) ELISA Kit	N/A	50 pg/mL
KA1713	FGF6 (Human) ELISA Kit	N/A	6 pg/mL
KA1714	FGF7 (Human) ELISA Kit	N/A	0.5 pg/mL
KA0337	FGF9 (Human) ELISA Kit	62.5 to 4000 pg/mL	15 pg/mL
KA1401	Fibrinogen (Dog) ELISA Kit	20 to 200 ng/mL	12 ng/mL
KA1402	Fibrinogen (Dog) ELISA Kit	1 to 10 ug/mL	0.6 ug/mL
KA0475	Fibrinogen (Human) ELISA Kit	0.16 to 40 ug/mL	0.1 ug/mL
KA0476	Fibrinogen (Human) ELISA Kit	1.25 to 80 ng/mL	1 ng/mL
KA0492	Fibrinogen (Mouse) ELISA Kit	3.13 to 200 ng/mL	2 ng/mL
KA1399	Fibrinogen (Mouse) ELISA Kit	0.63 to 40 ug/mL	0.63 ug/mL
KA1964	Fibrinogen (Pig) ELISA Kit	6.125 to 400 ng/mL	N/A
KA0505	Fibrinogen (Rat) ELISA Kit	3.125 to 200 ng/mL	2 ng/mL
KA1405	Fibrinogen (Rat) ELISA Kit	500 to 10000 ng/mL	300 ng/mL
KA0162	FLT1 (Human) ELISA Kit	0.16 to 10 ng/mL	0.0295 ng/mL
KA0665	FLT1 (Human) ELISA Kit	31.2 to 2000 pg/mL	4 pg/mL
KA1819	Flt1 (Mouse) ELISA Kit	N/A	8 pg/mL
KA1784	Flt3l (Mouse) ELISA Kit	N/A	15 pg/mL
KA1716	FLT3LG (Human) ELISA Kit	N/A	1.5 pg/mL
KA1773	FLT4 (Human) ELISA Kit	N/A	15 pg/mL
KA0188	FN1 (Human) ELISA Kit	0.31 to 20 ng/mL	0.1 ng/mL
KA0491	Fn1 (Mouse) ELISA Kit	3.13 to 200 ng/mL	2 ng/mL
KA1253	Fn1 (Mouse) ELISA Kit	156 to 10000 pg/mL	15 pg/mL
KA0339	Fn1 (Rat) ELISA Kit	156 to 10000 pg/mL	15 pg/mL
KA0210	Free beta-CG (Human) ELISA Kit	0 to 50 ng/mL	0.25 ng/mL
KA1211	Free Fatty Acid (Human) ELISA Kit	0.47 to 30 nmol/mL	0.12 nmol/mL
KA0209	Free PSA (Human) ELISA Kit	0 to 25 ng/mL	0.1 ng/mL
KA0213	FSH (Human) ELISA Kit	0 to 200 mIU/mL	2.5 mIU/mL
KA1715	FSTL3 (Human) ELISA Kit	N/A	250 pg/mL
KA1859	FURIN (Human) ELISA Kit	N/A	15 pg/mL
KA0305	GAST (Human) ELISA Kit	39.1 to 10000 pg/mL	7.27 pg/mL
KA0319	Gast (Rat) ELISA Kit	78.1 to 5000 pg/mL	78.1 pg/mL
KA0479	GC (Human) ELISA Kit	0.025 to 100 ug/mL	0.02 ug/mL
KA1933	Gc (Rat) ELISA Kit	3.9 to 250 ng/mL	N/A
KA2018	Gc (Rat) ELISA Kit	3.125 to 100 ng/mL	N/A
KA1871	GDF15 (Human) ELISA Kit	N/A	2 pg/mL
KA0984	GDNF (Human) ELISA Kit	31.2 to 2000 pg/mL	4 pg/mL
KA0341	Gdnf (Rat) ELISA Kit	31.2 to 2000 pg/mL	4 pg/mL
KA0024	GFAP (Human) ELISA Kit	0.25 to 25 ng/mL	0.33 ng/mL
KA0215	GH (Human) ELISA Kit	0 to 50 ng/mL	5 ng/mL
KA0521	GH (Human) ELISA Kit	0.1 to 12 ng/mL	0.1 ng/mL
KA1189	GHR (Human) ELISA Kit	7.8 to 500 pmol/ml	N/A
KA1909	Glutamate ELISA Kit	N/A	0.3 ug/mL
KA1137	GNLY (Human) ELISA Kit	31.2 to 2000 pg/mL	7.8 pg/mL
KA1175	GPC3 (Human) ELISA Kit	15.6 to 1,000 pg/mL	3.9 pg/mL
KA0530	GPX1 (Human) ELISA Kit	1.56 to 100 ng/mL	1.56 ng/mL
KA0318	GSK-3-beta (total) ELISA Kit	78.1 to 5000 pg/mL	74.4 pg/mL
KA2076	Gsta (Rat) ELISA Kit	3.125 to 100 ng/mL	N/A
KA0186	GZMA (Human) ELISA Kit	0.8 to 600 pg/mL	0.4 pg/mL
KA0187	GZMB (Human) ELISA Kit	0.66 to 480 pg/mL	0.2 pg/mL
KA2011	Haptoglobin (Horse) ELISA Kit	18.75 to 600 ng/mL	N/A
KA2052	Haptoglobin (Horse) ELISA Kit	3.13 to 100 ng/mL	N/A
KA1243	HAVCR1 (Human) ELISA Kit	0.78 to 50 ng/mL	0.39 ng/mL
KA1064	Havcr1 (Mouse) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA1410	Havcr1 (Rat) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA1723	HBEGF (Human) ELISA Kit	N/A	20 pg/mL
KA2019	Hemoglobin (Horse) ELISA Kit	12.5 to 400 ng/mL	N/A




KA0527




KA0665




KA0188



KA1933




KA0341




KA1064


Note:
 7. Epinephrine: 1.3 to 81.4 ng/mL, Norepinephrine: 4.9 to 339 ng/mL
 8. Epinephrine: 0.3 ng/mL urine, 10 pg/mL plasma, Norepinephrine: 1.5 ng/mL urine, 50 pg/mL plasma
 9. Epinephrine: 1.3 to 81.4 ng/mL, Norepinephrine: 4.9 to 339 ng/mL, Dopamine: 14 to 917 ng/mL
 10. Epinephrine: 0.3 ng/mL urine, 10 pg/mL plasma, Norepinephrine: 1.5 ng/mL urine, 50 pg/mL plasma, Dopamine: 4.5 ng/mL urine, 20 pg/mL plasma




KA0211




KA0343




KA1981



KA0345



KA0350




KA0353

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA1942	Hemoglobin (Mouse) ELISA Kit	12.5 to 400 ng/mL	N/A
KA1989	Hemoglobin (Rat) ELISA Kit	6.25 to 200 ng/mL	N/A
KA1192	Hepc (Human) ELISA Kit	0.187 to 12 ng/mL	0.094 ng/mL
KA0211	H-ferritin (Human) ELISA Kit	0 to 1000 ng/mL	5 ng/mL
KA0343	HGF (Human) ELISA Kit	31.2 to 2000 pg/mL	5 pg/mL
KA1785	Hgf (Mouse) ELISA Kit	N/A	100 pg/mL
KA1786	Hgf (Mouse) ELISA Kit	N/A	100 pg/mL
KA1247	HIF1A (Human) ELISA Kit	6.25 to 400 ng/mL	1.56 pg/mL
KA1171	Histamine (Human) ELISA Kit	1.56 to 100 ng/mL	0.78 ng/mL
KA1888	Histamine ELISA Kit	N/A	75 ng/g
KA1019	Hmox1 (Rat) ELISA Kit	0.195 to 12.5 ng/ml	0.036 ng/mL
KA1242	Homocysteine (Human) ELISA Kit	15.6 to 1000 pmol/mL	3.9 pmol/mL
KA1850	HP (Bovine) ELISA Kit	N/A	3.7 ng/mL
KA1853	HP (Canine) ELISA Kit	N/A	0.3 ng/mL
KA1981	HP (Cat) ELISA Kit	0.39 to 25 ng/mL	N/A
KA2015	HP (Cat) ELISA Kit	6.25 to 400 ng/mL	N/A
KA0483	HP (Human) ELISA Kit	1.56 to 400 ng/mL	0.3 ng/mL
KA0533	HP (Human) ELISA Kit	0.78 to 50 ng/mL	0.78 ng/mL
KA2012	HP (Monkey) ELISA Kit	1.56 to 25 ng/mL	N/A
KA1921	Hp (Mouse) ELISA Kit	7.8 to 125 ng/mL	N/A
KA1930	HP (Rabbit) ELISA Kit	12.5 to 200 ng/mL	N/A
KA1403	Hp (Rat) ELISA Kit	1000 to 15000 ng/mL	20 ng/mL
KA1922	Hp (Rat) ELISA Kit	7.8 to 125 ng/mL	N/A
KA1852	HP (Swine) ELISA Kit	N/A	1.8 ng/mL
KA1230	HPSE (Human) ELISA Kit	1.56 to 100 U/L	0.78 U/L
KA0481	HPX (Human) ELISA Kit	0.078 to 5 ug/mL	0.07 ug/mL
KA0482	HPX (Human) ELISA Kit	6.25 to 400 ng/mL	6 ng/mL
KA1927	Hpx (Mouse) ELISA Kit	12.5 to 400 ng/mL	N/A
KA0981	HSP27 (Human) ELISA Kit	0.078 to 80 ng/mL	0.1 ng/mL
KA1404	HSP47 (Human) ELISA Kit	5 to 50 ng/mL	3 ng/mL
KA1018	HSPB1 (Human) ELISA Kit	0.39 to 25 ng/mL	0.39 ng/mL
KA1065	HSPB1 (Human) ELISA Kit	62.5 to 4000 pg/mL	5 pg/mL
KA1843	HSPD1 (Human) ELISA Kit	N/A	2 ng/mL
KA1212	Hyp (Human) ELISA Kit	62.5 to 4000 pmol/mL	15.6 pmol/mL
KA0113	ICAM1 (Human) ELISA Kit	6.25 to 100 ng/mL	2.17 ng/mL
KA0261	ICAM1 (Monkey) ELISA Kit	6.25 to 100 U/mL	1.05 U/mL
KA0345	Icam1 (Mouse) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA0346	Icam1 (Rat) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA0142	ICAM2 (Human) ELISA Kit	0.75 to 24 U/mL	0.3 U/mL
KA0128	ICAM3 (Human) ELISA Kit	0.8 to 50 ng/mL	0.38 ng/mL
KA0126	IFNA1 (Human) ELISA Kit	7.8 to 500 pg/mL	3.16 pg/mL
KA0127	IFN-alpha Ab (Human) ELISA Kit	3.2 to 200 ng/mL	1.38 ng/mL
KA0135	IFNG (Human) ELISA Kit	1.6 to 100 pg/mL	0.99 pg/mL
KA0242	Ifng (Mouse) ELISA Kit	15.6 to 1000 pg/mL	5.3 pg/mL
KA0282	IFNG (Pig) ELISA Kit	31.3 to 2000 pg/mL	3.3 pg/mL
KA0273	Ifng (Rat) ELISA Kit	31.3 to 2000 pg/mL	9.9 pg/mL
KA0139	IFNW1 (Human) ELISA Kit	4.7 to 300 pg/mL	1.51 pg/mL
KA0349	IGF1 (Human) ELISA Kit	62.5 to 4000 pg/mL	10 pg/mL
KA0493	Igf1 (Mouse) ELISA Kit	0.16 to 10 ng/mL	0.1 ng/mL
KA0546	Igf1 (Rat) ELISA Kit	62.5 to 4000 pg/mL	5 pg/mL
KA0350	IGF2 (Human) ELISA Kit	62.5 to 4000 pg/mL	2 pg/mL
KA1791	Igf2 (Mouse) ELISA Kit	N/A	1.5 ng/mL
KA1872	IGFBP1 (Human) ELISA Kit	N/A	5 pg/mL
KA1726	IGFBP2 (Human) ELISA Kit	N/A	20 pg/mL
KA1787	Igfbp2 (Mouse) ELISA Kit	N/A	0.6 ng/mL
KA0353	IGFBP3 (Human) ELISA Kit	156.2 to 10000 pg/mL	10 pg/mL
KA1788	Igfbp3 (Mouse) ELISA Kit	N/A	70 pg/mL
KA1873	IGFBP4 (Human) ELISA Kit	N/A	0.25 ng/mL
KA1727	IGFBP5 (Human) ELISA Kit	N/A	2 ng/mL
KA1789	Igfbp5 (Mouse) ELISA Kit	N/A	5 pg/mL
KA1728	IGFBP6 (Human) ELISA Kit	N/A	150 pg/mL
KA1790	Igfbp6 (Mouse) ELISA Kit	N/A	6 pg/mL
KA1200	IL1 (Human) ELISA Kit	15.6 to 1000 pg/mL	3.9 pg/mL
KA0125	IL10 (Human) ELISA Kit	3.2 to 200 pg/mL	0.99 pg/mL
KA0262	IL10 (Monkey) ELISA Kit	3.9 to 250 U/mL	1.32 U/mL

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA0243	IL10 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	5 pg/mL
KA1038	IL10 (Mouse) ELISA Kit	31 to 2000 pg/mL	100 pg/mL
KA0274	IL10 (Rat) ELISA Kit	15.6 to 1000 pg/mL	1.5 pg/mL
KA1731	IL10RB (Human) ELISA Kit	N/A	5 pg/mL
KA0526	IL11 (Human) ELISA Kit	4 to 1600 pg/mL	4 pg/mL
KA1792	IL11 (Mouse) ELISA Kit	N/A	1 ng/mL
KA0647	IL12 (Human) ELISA Kit	7.8 to 500 pg/mL	7 pg/mL
KA0143	IL12 (p70) (Human) ELISA Kit	3.2 to 200 pg/mL	2.1 pg/mL
KA0263	IL12 (p70) (Monkey) ELISA Kit	3.9 to 250 U/mL	1.6 U/mL
KA0244	IL12 (p70) (Mouse) ELISA Kit	15.6 to 1000 pg/mL	10 pg/mL
KA0245	IL12 (p70/ p40 alone) (Mouse) ELISA Kit	31.3 to 2000 pg/mL	1.9 pg/mL
KA0374	IL12B (Human) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA0137	IL13 (Human) ELISA Kit	1.6 to 100 pg/mL	0.73 pg/mL
KA0246	IL13 (Mouse) ELISA Kit	7.8 to 500 pg/mL	2.8 pg/mL
KA1408	IL13 (Rat) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0524	IL13RA2 (Human) ELISA Kit	100 to 6400 pg/mL	10 pg/mL
KA0376	IL15 (Human) ELISA Kit	15.6 to 1000 pg/mL	3 pg/mL
KA1793	IL15 (Mouse) ELISA Kit	N/A	0.5 ng/mL
KA0520	IL16 (Human) ELISA Kit	20 to 3200 pg/mL	20 pg/mL
KA0180	IL17A (Human) ELISA Kit	1.6 to 100 pg/mL	0.5 pg/mL
KA1732	IL17A (Human) ELISA Kit	N/A	10 pg/mL
KA0247	IL17a (Mouse) ELISA Kit	7.8 to 500 pg/mL	1.6 pg/mL
KA1001	IL17a (Rat) ELISA Kit	1.6 to 100 pg/mL	1 pg/mL
KA1794	IL17b (Mouse) ELISA Kit	N/A	12 pg/mL
KA1733	IL17F (Human) ELISA Kit	N/A	20 pg/mL
KA1000	IL17f (Mouse) ELISA Kit	62.5 to 4000 pg/mL	10 pg/mL
KA1734	IL17RA (Human) ELISA Kit	N/A	10 pg/mL
KA0561	IL18 (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0283	IL18 (Pig) ELISA Kit	39 to 2500 pg/mL	23.1 pg/mL
KA0439	IL18 (Rat) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1735	IL18BP (Human) ELISA Kit	N/A	20 pg/mL
KA0146	IL1A (Human) ELISA Kit	3.2 to 200 pg/mL	0.66 pg/mL
KA0249	IL1a (Mouse) ELISA Kit	7.8 to 500 pg/mL	3.9 pg/mL
KA0275	IL1a (Rat) ELISA Kit	15.6 to 1000 pg/mL	4.1 pg/mL
KA0356	IL1B (Human) ELISA Kit	1.56 to 100 pg/mL	0.15 pg/mL
KA0357	IL1b (Mouse) ELISA Kit	7.8 to 500 pg/mL	1 pg/mL
KA1820	IL1B (Pig) ELISA Kit	N/A	6 pg/mL
KA1502	IL1b (Rat) ELISA Kit	15.6 to 1000 pg/mL	N/A
KA1066	IL1RN (Human) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA0130	IL2 (Human) ELISA Kit	19 to 1200 pg/mL	9.9 pg/mL
KA0264	IL2 (Monkey) ELISA Kit	19.5 to 1250 U/mL	9.9 U/mL
KA0250	IL2 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	5.3 pg/mL
KA0276	IL2 (Rat) ELISA Kit	31.3 to 2000 pg/mL	9.6 pg/mL
KA0986	IL2 (Rat) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0194	IL21 (Human) ELISA Kit	78 to 5000 pg/mL	20 pg/mL
KA1795	IL21 (Mouse) ELISA Kit	N/A	6 pg/mL
KA0195	IL22 (Human) ELISA Kit	31.3 to 2000 pg/mL	5 pg/mL
KA0251	IL22 (Mouse) ELISA Kit	7.8 to 500 pg/mL	1.6 pg/mL
KA0183	IL-23 (Human) ELISA Kit	78 to 5000 pg/mL	10 pg/mL
KA1409	IL23a (Mouse) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA1736	IL28A (Human) ELISA Kit	N/A	15 pg/mL
KA1796	IL28b (Mouse) ELISA Kit	N/A	2 pg/mL
KA0997	IL29 (Human) ELISA Kit	15.6 to 1000 pg/mL	2 pg/mL
KA0122	IL-2R (Human) ELISA Kit	0.31 to 20 ng/mL	0.27 ng/mL
KA0522	IL2RA (Human) ELISA Kit	10 to 2000 pg/mL	10 pg/mL
KA0360	IL3 (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0361	IL3 (Mouse) ELISA Kit	7.8 to 500 pg/mL	1 pg/mL
KA0998	IL31 (Human) ELISA Kit	7.8 to 500 pg/mL	1 pg/mL
KA0999	IL33 (Human) ELISA Kit	4 to 250 pg/mL	1 pg/mL
KA1143	IL35 (Human) ELISA Kit	25 to 1600 pg/mL	6.2 pg/mL
KA0132	IL4 (Human) ELISA Kit	7.8 to 500 pg/mL	1.3 pg/mL
KA0252	IL4 (Mouse) ELISA Kit	3.9 to 250 pg/mL	1.98 pg/mL
KA0277	IL4 (Rat) ELISA Kit	1.6 to 100 pg/mL	0.2 pg/mL
KA0167	IL5 (Human) ELISA Kit	7.8 to 500 pg/mL	1.45 pg/mL
KA0987	IL5 (Human) ELISA Kit	15.6 to 1000 pg/mL	2 pg/mL




KA1408




KA0561




KA0356




KA1066




KA0360




KA0987




KA1067




KA2009




KA1063



KA0380



KA0559




KA0438

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA0253	IL5 (Mouse) ELISA Kit	7.8 to 500 pg/mL	3.3 pg/mL
KA0123	IL6 (Human) ELISA Kit	1.6 to 100 pg/mL	0.92 pg/mL
KA0265	IL6 (Monkey) ELISA Kit	1.9 to 120 U/mL	1.25 U/mL
KA0254	IL6 (Mouse) ELISA Kit	31.3 to 2000 pg/mL	7.92 pg/mL
KA1821	IL6 (Pig) ELISA Kit	N/A	45 pg/mL
KA0278	IL6 (Rat) ELISA Kit	31.3 to 2000 pg/mL	12 pg/mL
KA1503	IL6 (Rat) ELISA Kit	62.5 to 4000 pg/mL	N/A
KA0124	IL6R (Human) ELISA Kit	0.08 to 5 ng/mL	0.01 ng/mL
KA0523	IL6R (Human) ELISA Kit	31.2 to 2000 pg/mL	31.2 pg/mL
KA0525	IL7 (Human) ELISA Kit	1.5 to 50 pg/mL	1 pg/mL
KA1504	IL7 (Human) ELISA Kit	15.6 to 1000 pg/mL	N/A
KA1067	IL7 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0115	IL8 (Human) ELISA Kit	15.6 to 1000 pg/mL	2 pg/mL
KA0266	IL8 (Monkey) ELISA Kit	19.5 to 1250 U/mL	2.5 U/mL
KA1822	IL8 (Pig) ELISA Kit	N/A	10 pg/mL
KA1737	IL9 (Human) ELISA Kit	N/A	0.1 ng/mL
KA1797	IL9 (Mouse) ELISA Kit	N/A	0.1 ng/mL
KA1155	IMA (Human) ELISA Kit	15.6 to 1000 nmol/L	N/A
KA0328	INHBA (Human) ELISA Kit	62.5 to 4000 pg/mL	12 pg/mL
KA1690	INHBA (Human) ELISA Kit	N/A	15 pg/mL
KA0538	Inhba (Mouse) ELISA Kit	62.5 to 4000 pg/mL	10 pg/mL
KA0539	Inhba (Rat) ELISA Kit	62.5 to 4000 pg/mL	12 pg/mL
KA1683	Inhibin B (Human/Mouse/Rat) ELISA Kit	1 to 10000 ng/mL	3.46 ng/mL
KA0477	ITGA2B/ITGB3 (Human) ELISA Kit	1.25 to 40 ng/mL	1 ng/mL
KA1968	Kappa (Human) ELISA Kit	6.25 to 200 ng/mL	N/A
KA1812	Kitl (Mouse) ELISA Kit	N/A	5 pg/mL
KA1844	Kitl (Mouse) ELISA Kit	N/A	0.06 ng/mL
KA1022	KITLG (Human) ELISA Kit	0.156 to 40 ng/mL	0.02 ng/mL
KA2009	KLK1 (Human) ELISA Kit	156 to 10000 pg/mL	< 10 pg/mL
KA1063	KLK3 (Human) ELISA Kit	312 to 20000 pg/ml	10 pg/mL
KA1039	KNG1 (Human) ELISA Kit	1.56 to 100 ng/mL	1.5 ng/mL
KA1040	KNG1 (Human) ELISA Kit	125 to 4000 ng/mL	50 ng/mL
KA1935	Kng111 (Rat) ELISA Kit	6.25 to 400 ng/mL	N/A
KA1205	KS (Human) ELISA Kit	0.312 to 20 ng/mL	0.078 ng/mL
KA0178	L-12/IL-23 p40 (Human) ELISA Kit	31.3 to 2000 pg/mL	10 pg/mL
KA1179	LALBA (Human) ELISA Kit	0.625 to 40 ng/mL	0.156 ng/mL
KA1969	Lambda (Human) ELISA Kit	6.25 to 200 ng/mL	N/A
KA1999	Laminin (Human) ELISA Kit	156 to 10000 pg/mL	< 10 pg/mL
KA0380	Laminin (Mouse) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA0379	Laminin (Rat) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA0448	LBP (Human) ELISA Kit	5 to 50 ng/mL	5 ng/mL
KA0449	Lbp (Mouse) ELISA Kit	5 to 50 ng/mL	5 ng/mL
KA2023	LCN2 (Dog) ELISA Kit	3.125 to 100 ng/mL	N/A
KA0043	LCN2 (Human) ELISA Kit	2 to 100 ng/mL	0.02 ng/mL
KA1051	LCN2 (Human) ELISA Kit	46.9 to 3000 pg/mL	26.7 pg/mL
KA0559	Lcn2 (Mouse) ELISA Kit	156 to 10000 pg/mL	5 pg/mL
KA0560	Lcn2 (Rat) ELISA Kit	156 to 10000 pg/mL	15 pg/mL
KA1388	Lcn2 (Rat) ELISA Kit	N/A	6.02 pg/mL
KA1250	Leukotriene B4 (Human) ELISA Kit	6.25 to 400 pg/mL	3.12 pg/mL
KA0993	Lgals1 (Mouse) ELISA Kit	156 to 10000 pg/mL	< 5 pg/mL
KA0168	LGALS3 (Human) ELISA Kit	0.39 to 25 ng/mL	0.12 ng/mL
KA0140	LGALS3BP (Human) ELISA Kit	12.5 to 200 ng/mL	0.92 ng/mL
KA1717	LGALS7 (Human) ELISA Kit	N/A	20 pg/mL
KA0214	LH (Human) ELISA Kit	0 to 200 mIU/mL	1 mIU/mL
KA1740	LIF (Human) ELISA Kit	N/A	15 pg/mL
KA0149	LIFR (Human) ELISA Kit	0.08 to 5 ng/mL	0.052 ng/mL
KA1241	LIPC (Human) ELISA Kit	156 to 10000 pg/mL	39 pg/mL
KA1041	LPA (Human) ELISA Kit	0.78 to 50 ng/mL	0.8 ng/mL
KA0117	L-selectin (Human) ELISA Kit	0.39 to 25 ng/mL	0.19 ng/mL
KA0438	LTA (Human) ELISA Kit	62.5 to 4000 pg/mL	5 pg/mL
KA0484	LTF (Human) ELISA Kit	0.625 to 40 ng/mL	0.1 ng/mL
KA1284	Lysozyme ELISA Kit	0.78 to 50 ng/mL	N/A
KA1743	LYVE1 (Human) ELISA Kit	N/A	8 pg/mL
KA0485	LYZ (Human) ELISA Kit	0.313 to 10 ng/mL	0.1 ng/mL
KA1042	LYZ (Human) ELISA Kit	90 to 6000 ng/mL	50 ng/mL

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA1982	MB (Dog) ELISA Kit	3.91 to 250 ng/mL	N/A
KA2068	Mb (Guinea Pig) ELISA Kit	1.563 to 100 ng/mL	N/A
KA0237	MB (Human) ELISA Kit	0 to 1000 ng/mL	5 ng/mL
KA1984	MB (Monkey) ELISA Kit	1.95 to 125 ng/mL	N/A
KA1947	Mb (Mouse) ELISA Kit	4.69 to 300 ng/mL	N/A
KA1986	MB (Pig) ELISA Kit	0.49 to 31.25 ng/mL	N/A
KA2053	MB (Rabbit) ELISA Kit	1.95 to 125 ng/mL	N/A
KA1954	Mb (Rat) ELISA Kit	3.125 to 200 ng/mL	N/A
KA1183	MCT (Human) ELISA Kit	1.56 to 100 ng/mL	0.39 ng/mL
KA0028	MDK (Human) ELISA Kit	2 to 10 ng/mL	0.33 ng/mL
KA1166	Melatonin (Human) ELISA Kit	3.12 to 200 pg/mL	1.56 pg/mL
KA2005	MERTK (Human) ELISA Kit	0.312 to 20 ng/mL	< 10 pg/mL
KA1889	Metanephrine ELISA Kit	41 to 2100 pg/mL	10 pg/mL
KA1890	Metanephrine ELISA Kit	40 to 1600 ng/mL	13 ng/mL
KA1878	Metanephrine/Normetanephrine ELISA Kit	40 to 2100 pg/mL ¹¹	10 pg/mL ¹²
KA1879	Metanephrine/Normetanephrine ELISA Kit	40 to 1600 pg/mL ¹³	5 ng/mL ¹⁴
KA1802	Mfge8 (Mouse) ELISA Kit	N/A	1 pg/mL
KA0557	MICA (Human) ELISA Kit	62.5 to 4000 pg/mL	10 pg/mL
KA1746	MICB (Human) ELISA Kit	N/A	75 pg/mL
KA1978	Microalbumin (Dog) ELISA Kit	12.5 to 400 ng/mL	N/A
KA0558	MIF (Human) ELISA Kit	156 to 10000 pg/mL	20 pg/mL
KA0390	MMP1 (Human) ELISA Kit	156 to 10000 pg/mL	20 pg/mL
KA1750	MMP10 (Human) ELISA Kit	N/A	3 pg/mL
KA0182	MMP13 (Human) ELISA Kit	0.31 to 20 ng/mL	0.18 ng/mL
KA1600	MMP13 (Human) ELISA Kit	32 to 2000 pg/mL	7 pg/mL
KA0391	MMP2 (Human) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA0392	Mmp2 (Mouse) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA0393	Mmp2 (Rat) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA0394	MMP3 (Human) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA1807	Mmp3 (Mouse) ELISA Kit	N/A	10 pg/mL
KA0395	MMP7 (Human) ELISA Kit	156 to 10000 pg/mL	6 pg/mL
KA0396	MMP8 (Human) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA1833	Mmp8 (Rat) ELISA Kit	N/A	50 pg/mL
KA0179	MMP9 (Human) ELISA Kit	0.23 to 15 ng/mL	0.05 ng/mL
KA0398	Mmp9 (Mouse) ELISA Kit	156 to 10000 pg/mL	20 pg/mL
KA0537	MPO (Human) ELISA Kit	0.78 to 50 ng/mL	0.78 ng/mL
KA1752	MST1 (Human) ELISA Kit	N/A	8 pg/mL
KA1180	MUC4 (Human) ELISA Kit	0.312 to 20 ng/mL	0.08 ng/mL
KA1988	Murinoglobulin (Rat) ELISA Kit	3.95 to 125 ng/mL	N/A
KA2063	Myl1 (Mouse, Rat) ELISA Kit	1.56 to 50 ng/mL	N/A
KA2062	Myl3 (Mouse, Rat) ELISA Kit	0.078 to 2.5 ng/mL	N/A
KA1865	NAMPT (Human/Mouse/Rat) ELISA Kit	0.1 to 1000 ng/mL	379 pg/mL
KA0095	Nampt (Mouse) ELISA Kit	0.75 to 60 ng/mL	0.24 ng/mL
KA1704	NBL1 (Human) ELISA Kit	N/A	35 pg/mL
KA2003	NCAM1 (Human) ELISA Kit	78 to 5000 pg/mL	< 10 pg/mL
KA1478	NEFL (Human) ELISA Kit	N/A	50 pg/mL
KA1684	Nesfatin (Human, Mouse, Rat) ELISA Kit	0.1 to 1000 ng/mL	147 pg/mL
KA1208	NGB (Human) ELISA Kit	3.12 to 200 ng/mL	0.78 ng/mL
KA0399	NGF (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0400	Ngf (Mouse) ELISA Kit	31.2 to 2000 pg/mL	1 pg/mL
KA0401	Ngf (Rat) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0445	Nitrotyrosine ELISA Kit	2.1 nM to 1500 nM	2 nM
KA1165	Norepinephrine (Human) ELISA Kit	31.2 to 2,000 pg/mL	7.8 pg/mL
KA1891	Norepinephrine ELISA Kit	20 to 339 ng/mL ¹⁵	1.5 ng/mL ¹⁶
KA1892	Normetanephrine ELISA Kit	58 to 5800 pg/mL	23 pg/mL
KA1893	Normetanephrine ELISA Kit	40 to 5200 ng/mL	23 ng/mL
KA1173	NPHS1 (Human) ELISA Kit	0.156 to 10 ng/mL	0.078 ng/mL
KA1680	NPPA (Human/Mouse/Rat) ELISA Kit	0 to 100 ng/mL	230 pg/mL
KA0503	Nppa (Rat) ELISA Kit	0.25 to 4 ng/mL	0.2 ng/mL
KA1861	NPPB (Human, Mouse, Rat) ELISA Kit	0.1 to 1000 pg/mL	1.45 pg/mL
KA0504	Nppb (Rat) ELISA Kit	0.031 to 2 ng/mL	0.03 ng/mL
KA0979	Nppb (Rat) ELISA Kit	0.031 to 1 ng/mL	0.03 ng/mL
KA1864	NPY (Human/Mouse/Rat) ELISA Kit	0.1 to 1000 ng/mL	180 pg/mL
KA1754	NRG1-beta 1 (Human) ELISA Kit	N/A	50 pg/mL
KA0402	NTF3 (Human) ELISA Kit	15.6 to 1000 pg/mL	2 pg/mL




KA2005




KA0391




KA0395



KA0396



KA0399




KA0400


Note:

11. Metanephrine: 40 to 2100 pg/mL, Normetanephrine: 58 to 5800 pg/mL
 13. Metanephrine: 40 to 1600 pg/mL, Normetanephrine: 40 to 5200 pg/mL
 15. 20 to 339 ng/mL urine, 318 to 2436 pg/mL plasma


12. Metanephrine: 10 pg/mL, Normetanephrine: 20 pg/mL
 14. Metanephrine: 5 ng/mL, Normetanephrine: 20 ng/mL
 16. 1.5 ng/mL urine, 50 pg/mL plasma




KA0404




KA1507




KA1676



KA0409



KA1509




KA0425

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA1860	NTF3 (Human) ELISA Kit	N/A	4 pg/mL
KA0403	Ntf3 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	2 pg/mL
KA0404	Ntf3 (Rat) ELISA Kit	15.6 to 1000 pg/mL	2 pg/mL
KA1755	NTF4 (Human) ELISA Kit	N/A	2 pg/mL
KA1507	Ntrk1 (Rat) ELISA Kit	156 to 10000 pg/mL	N/A
KA0996	Ntrk2 (Mouse) ELISA Kit	93.8 to 6000 pg/mL	5 pg/mL
KA1218	NTX (Human) ELISA Kit	0.47 to 30 nM BCE	0.24 nM BCE
KA1393	OLR1 (Human) ELISA Kit	N/A	12.4 pg/mL
KA1676	Olr1 (Mouse) ELISA Kit	31.2 to 2000 pg/mL	1 pg/mL
KA1023	ORM (Human) ELISA Kit	3.13 to 200 ng/mL	5 ng/mL
KA1961	ORM1 (Chicken) ELISA Kit	2.3 to 150 ng/mL	N/A
KA2069	ORM1 (Chicken) ELISA Kit	6.25 to 200 ng/mL	N/A
KA2066	ORM1 (Monkey) ELISA Kit	3.13 to 100 ng/mL	N/A
KA2075	ORM1 (Turkey) ELISA Kit	4.69 to 300 ng/mL	N/A
KA0480	ORM2 (Human) ELISA Kit	0.063 to 4 ug/mL	0.08 ug/mL
KA1756	OSM (Human) ELISA Kit	N/A	1 pg/mL
KA1147	OTC (Human) ELISA Kit	3.12 to 200 U/L	1.5 U/L
KA1856	Pancreatic Amylase (Human) ELISA Kit	N/A	0.3 mU/mL
KA1148	PAPPA (Human) ELISA Kit	6.25 to 400 pg/mL	3.12 pg/mL
KA0100	PARK7 (Human) ELISA Kit	1.56 to 100 ng/mL	0.55 ng/mL
KA1072	PCSK9 (Human) ELISA Kit	0.16 to 10 ng/mL	0.154 ng/mL
KA1057	Pcsk9 (Mouse/Rat) ELISA Kit	94 to 6000 pg/mL	38.1 pg/mL
KA1227	PCT (Human) ELISA Kit	N/A	3.9 pg/mL
KA1240	PCX (Human) ELISA Kit	156 to 10000 pg/mL	78 pg/mL
KA1758	PDGFA (Human) ELISA Kit	N/A	< 40 pg/mL
KA0408	PDGF-AB (Human) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA0409	Pdgf-ab (Rat) ELISA Kit	31.2 to 2000 pg/mL	3 pg/mL
KA1759	PDGFB (Human) ELISA Kit	N/A	1 pg/mL
KA1760	PDGFB (Human) ELISA Kit	N/A	1 pg/mL
KA0136	PECAM1 (Human) ELISA Kit	0.48 to 30 ng/mL	0.06 ng/mL
KA1508	PECAM1 (Human) ELISA Kit	156 to 10000 pg/mL	N/A
KA0267	PECAM1 (Monkey) ELISA Kit	0.48 to 30 U/mL	0.06 U/mL
KA1509	Pecam1 (Mouse) ELISA Kit	156 to 10000 pg/mL	N/A
KA1761	PF4 (Human) ELISA Kit	N/A	20 pg/mL
KA1809	Pf4 (Mouse) ELISA Kit	N/A	6 pg/mL
KA0410	PGF (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1762	PGF (Human) ELISA Kit	N/A	2 pg/mL
KA1479	Phosphorylated NEFH (Human) ELISA Kit	N/A	31 pg/mL
KA1771	PI3 (Human) ELISA Kit	N/A	5 pg/mL
KA0157	PLAT (Human) ELISA Kit	15.6 to 1000 pg/mL	6.4 pg/mL
KA0425	PLAU (Human) ELISA Kit	62.5 to 4000 pg/mL	5 pg/mL
KA0426	PLAUR (Human) ELISA Kit	62.5 to 4000 pg/mL	4 pg/mL
KA1851	PLG (Bovine) ELISA Kit	N/A	3 ng/mL
KA1858	PLG (Bovine) ELISA Kit	N/A	300 ng/mL
KA1939	PLG (Dog) ELISA Kit	9.38 to 600 ng/mL	N/A
KA0494	PLG (Human) ELISA Kit	1.25 to 80 ng/mL	1 ng/mL
KA1845	Plg (Mouse) ELISA Kit	N/A	250 ng/mL
KA1846	Plg (Mouse) ELISA Kit	N/A	0.18 ng/mL
KA1959	PLG (Pig) ELISA Kit	3.125 to 200 ng/mL	N/A
KA1948	PLG (Rabbit) ELISA Kit	6.25 to 200 ng/mL	N/A
KA1960	Plg (Rat) ELISA Kit	3.13 to 200 ng/mL	N/A
KA0163	PMN-Elastase (Human) ELISA Kit	0.16 to 10 ng/mL	1.98 ng/mL
KA1134	PNLIP (Human) ELISA Kit	4.69 to 300 U/L	1.17 U/L
KA1678	PNPLA3 (Human/Mouse/Rat) ELISA Kit	0.1 to 1000 ng/mL	21 pg/mL
KA0112	Poly-Ubiquitinated Protein (Human) ELISA Kit	1.56 to 100 U/mL	1.182 U/mL
KA1215	PP (Human) ELISA Kit	3.12 to 200 ng/mL	0.78 ng/mL
KA1753	PPBP (Human) ELISA Kit	N/A	8.5 pg/mL
KA1176	PPIA (Human) ELISA Kit	0.156 to 10 ng/mL	0.039 ng/mL
KA0536	PRDX1 (Human) ELISA Kit	0.5 to 32 ng/mL	0.5 ng/mL
KA0531	PRDX3 (Human) ELISA Kit	2 to 32 ng/mL	2 ng/mL
KA2070	Prealbumin (Mouse) ELISA Kit	15.6 to 500 ng/mL	N/A
KA1912	Pregnenolone ELISA Kit	N/A	0.05 ng/mL
KA0217	PRL (Human) ELISA Kit	0 to 200 ng/mL	2 ng/mL
KA0923	Prl (Mouse/Rat) ELISA Kit	0.2 to 200 ng/mL	0.2 ng/mL
KA1834	Prlr (Rat) ELISA Kit	N/A	0.4 ng/mL

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA1839	PROC (Human) ELISA Kit	N/A	900 ng/mL
KA1840	PROC (Human) ELISA Kit	N/A	3 ng/mL
KA1763	Procalcitonin (Human) ELISA Kit	N/A	30 pg/mL
KA0235	Progesterone (Human) ELISA Kit	0 to 50 ng/mL	0.3 ng/mL
KA0299	Progesterone ELISA Kit	15.62 to 500 pg/mL	8.57 pg/mL
KA0031	Proguanylin (Human) ELISA Kit	0.3 to 10 ng/mL	0.09 ng/mL
KA1675	PROK1 (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0304	Prostacyclin ELISA Kit	7.81 to 2000 pg/mL	6.58 pg/mL
KA0296	Prostaglandin E1 ELISA Kit	4.88 to 5000 pg/mL	5.58 pg/mL
KA0324	Prostaglandin E2 ELISA Kit	7.81 to 1000 pg/mL	6.03 pg/mL
KA0325	Prostaglandin E2 ELISA Kit	1562 to 100000 pg/mL	684 pg/mL
KA0326	Prostaglandin E2 ELISA Kit	7.8 to 1000 pg/mL	8.26 pg/mL
KA0310	Prostaglandin F2 alpha ELISA Kit	3.05 to 50000 pg/mL	6.71 pg/mL
KA0327	Prostaglandin F2 alpha ELISA Kit	1.95 to 2000 pg/mL	0.98 pg/mL
KA0496	Prothrombin (Human) ELISA Kit	0.13 to 32 ug/mL	0.13 ug/mL
KA0208	PSA (Human) ELISA Kit	0 to 120 ng/mL	1 ng/mL
KA0268	P-selectin (Monkey) ELISA Kit	0.8 to 50 U/mL	1.3 U/mL
KA1228	PTGS2 (Human) ELISA Kit	1.56 to 100 ng/mL	0.39 ng/mL
KA0924	PTH (Human) ELISA Kit	0 to 210 pg/mL	0.44 ng/mL
KA1154	PTH LH (Human) ELISA Kit	15.6 to 1,000 pg/mL	3.9 pg/mL
KA1007	PTX3 (Human) ELISA Kit	0.5 to 20 ng/mL	0.1 ng/mL
KA1181	PYD (Human) ELISA Kit	6.25 to 400 nmol/L	2.5 nmol/L
KA1686	PYY (Human/Mouse/Rat) ELISA Kit	0.1 to 1000 pg/mL	2.84 pg/mL
KA2006	RAGE (Human) ELISA Kit	78 to 5000 pg/mL	< 10 pg/mL
KA1810	Rage (Mouse) ELISA Kit	N/A	15 pg/mL
KA2027	Rbp (Rat) ELISA Kit	6.25 to 200 ng/mL	N/A
KA0499	RBP4 (Human) ELISA Kit	6.25 to 400 ng/mL	6 ng/mL
KA0989	RBP4 (Human) ELISA Kit	312 to 20000 pg/mL	10 pg/mL
KA1053	RBP4 (Human) ELISA Kit	1.56 to 100 ng/mL	1.08 ng/mL
KA2008	Rbp4 (Mouse) ELISA Kit	0.132 to 20 ng/mL	< 10 pg/mL
KA0500	RCVRN (Human) ELISA Kit	0.078 to 20 ng/mL	0.05 ng/mL
KA1163	REN (Human) ELISA Kit	15.6 to 1,000 pg/mL	3.9 pg/mL
KA0498	RENBP (Human) ELISA Kit	0.031 to 2 ug/mL	0.008 ug/mL
KA1838	RENBP (Human) ELISA Kit	N/A	100 ng/mL
KA1140	RNASE3 (Human) ELISA Kit	0.625 to 40 ng/mL	0.156 ng/mL
KA0091	S100A12 (Human) ELISA Kit	78 to 5000 pg/mL	56 pg/mL
KA1392	S100A13 (Human) ELISA Kit	N/A	66.2 pg/mL
KA0092	S100A4 (Human) ELISA Kit	0.39 to 25 ng/mL	0.28 ng/mL
KA1054	S100A7 (Human) ELISA Kit	0.12 to 90 ng/mL	0.12 ng/mL
KA1049	S100A8 (Human) ELISA Kit	78.1 to 5000 pg/mL	43.4 pg/mL
KA1050	S100A9 (Human) ELISA Kit	50 to 3200 pg/mL	6.55 pg/mL
KA0037	S100B (Human) ELISA Kit	0.05 to 2 ng/mL	0.005 ng/mL
KA0093	S100P (Human) ELISA Kit	39.0 to 2500 pg/mL	29 pg/mL
KA0518	SAA (Human) ELISA Kit	1.1 to 80 ng/mL	1.1 ng/mL
KA1162	SAA (Human) ELISA Kit	3.12 to 200 ng/mL	0.78 ng/mL
KA2014	SAA (Monkey) ELISA Kit	0.78 to 25 ng/mL	N/A
KA1928	Saa (Mouse) ELISA Kit	31.25 to 2000 ng/mL	N/A
KA0020	SCGB1A1 (Human) ELISA Kit	2 to 100 ng/mL	0.02 ng/mL
KA1875	SDF-1alpha (Human) ELISA Kit	N/A	10 pg/mL
KA1876	SDF-1beta (Human) ELISA Kit	N/A	35 pg/mL
KA1058	SELL (Human) ELISA Kit	62.5 to 4000 pg/mL	5 pg/mL
KA1059	Sell (Mouse) ELISA Kit	62.5 to 4000 pg/mL	5 pg/mL
KA1830	Sell (Rat) ELISA Kit	N/A	6 pg/mL
KA0548	SELP (Human) ELISA Kit	62.5 to 4000 pg/mL	5 pg/mL
KA0413	Selp (Mouse) ELISA Kit	62.5 to 4000 pg/mL	5 pg/mL
KA0154	SELPLG (Human) ELISA Kit	1.6 to 50 U/mL	0.99 U/mL
KA1151	Serotonin (Human) ELISA Kit	93.75 to 6,000 pg/mL	23.4 pg/mL
KA1894	Serotonin ELISA Kit	40 to 860 ng/mL	5 ng/mL
KA1936	SERPINA1 (Dog) ELISA Kit	12.5 to 800 ng/mL	N/A
KA0458	SERPINA1 (Human) ELISA Kit	0.039 to 40 ug/mL	0.020 ug/mL
KA0459	SERPINA1 (Human) ELISA Kit	0.39 to 100 ng/mL	0.4 ng/mL
KA1979	Serpina1 (Mouse) ELISA Kit	6.25 to 200 ng/mL	N/A
KA1688	SERPINA12 (Human/Mouse/Rat) ELISA Kit	1 to 100 ng/mL	2.62 ng/mL
KA0456	SERPINC1 (Human) ELISA Kit	3.125 to 200 ng/mL	2 ng/mL
KA1032	SERPINC1 (Human) ELISA Kit	4 to 4000 ng/mL	40 ng/mL




KA1053




KA1054




KA0518



KA1058




KA0548




KA0413


ELISA Kits




KA0547




KA1068




KA0416



KA0667



KA0437




KA0549

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA1836	SERPINC1 (Human) ELISA Kit	N/A	62.5 ng/mL
KA0490	Serpinc1 (Mouse) ELISA Kit	1.56 to 400 ng/mL	1 ng/mL
KA1980	Serpinc1 (Rat) ELISA Kit	12.5 to 400 ng/mL	N/A
KA0190	SERPINE1 (Human) ELISA Kit	78 to 5000 pg/mL	29 pg/mL
KA1138	SFTPB (Human) ELISA Kit	7.8 to 500 ng/mL	1.95 ng/mL
KA1139	SFTPC (Human) ELISA Kit	0.47 to 30 ng/mL	0.12 ng/mL
KA1251	SFTPD (Human) ELISA Kit	1.56 to 100 ng/mL	0.01 ng/mL
KA1913	SHBG (Human) ELISA Kit	N/A	0.1 nmol/L
KA1765	SIGLEC5 (Human) ELISA Kit	N/A	8 pg/mL
KA1144	SLC2A4 (Human) ELISA Kit	0.156 to 10 ng/mL	0.039 ng/mL
KA1191	SNCA (Human) ELISA Kit	0.78 to 50 ng/mL	0.195 ng/mL
KA0528	SOD2 (Human) ELISA Kit	78.12 to 5000 pg/mL	25 pg/mL
KA0529	SOD3 (Human) ELISA Kit	2 to 128 ng/mL	2 ng/mL
KA1246	SPARC (Human) ELISA Kit	3.12 to 200 ng/mL	N/A
KA1256	Spermatozoa Ab ELISA Kit	0 to 60 U/mL	N/A
KA1257	Spermatozoa Ab ELISA Kit	0 to 60 U/mL	N/A
KA1258	Spermatozoa IgG ELISA Kit	0 to 60 Units/mL	N/A
KA0407	SPP1 (Human) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA0547	Spp1 (Mouse) ELISA Kit	156 to 10000 pg/mL	10 pg/mL
KA0312	Spp1 (Rat) ELISA Kit	3.13 to 100 ng/mL	3.03 ng/mL
KA0302	Substance P ELISA Kit	9.76 to 10000 pg/mL	8.04 pg/mL
KA0199	T3 (Free) (Human) ELISA Kit	0 to 19 pg/mL	0.05 pg/mL
KA0198	T3 (Total) (Human) ELISA Kit	0 to 10 ng/mL	0.2 ng/mL
KA0925	T3 (Total) (Mouse/Rat) ELISA Kit	0.25 to 7.5 ng/mL	0.125 ng/mL
KA0200	T4 (Total) (Human) ELISA Kit	0 to 25 ug/dL	0.4 ug/dL
KA0508	TAT Complexes (Human) ELISA Kit	1.48 to 120 ng/mL	1 ng/mL
KA1857	TAT Complexes (Mouse) ELISA Kit	N/A	0.125 ng/mL
KA0193	TEK (Human) ELISA Kit	0.16 to 10 ng/mL	0.03 ng/mL
KA0236	Testosterone (Human) ELISA Kit	0 to 18 ng/mL	0.05 ng/mL
KA0309	Testosterone ELISA Kit	7.81 to 2000 pg/mL	5.67 pg/mL
KA2047	TF (Bovine) ELISA Kit	9.375 to 600 ng/mL	N/A
KA1924	TF (Dog) ELISA Kit	6.25 to 400 ng/mL	N/A
KA0509	TF (Human) ELISA Kit	0.098 to 25 ug/mL	100 ng/mL
KA0510	TF (Human) ELISA Kit	1.56 to 100 ng/mL	1.5 ng/mL
KA1847	Tf (Rat) ELISA Kit	N/A	150 ng/mL
KA1848	Tf (Rat) ELISA Kit	N/A	0.3 ng/mL
KA1237	TFF3 (Human) ELISA Kit	6.25 to 400 pg/mL	1.0 pg/mL
KA0507	TFPI (Human) ELISA Kit	0.313 to 20 ng/mL	0.3 ng/mL
KA0039	TFRC (Human) ELISA Kit	0.05 to 2 ng/mL	10 ng/mL
KA1068	TGFA (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0151	TGFB1 (Human) ELISA Kit	31.3 to 2000 pg/mL	9 pg/mL
KA0256	Tgfb1 (Mouse) ELISA Kit	31.3 to 2000 pg/mL	9 pg/mL
KA0279	Tgfb1 (Rat) ELISA Kit	0.48 to 30 ng/mL	0.048 ng/mL
KA0416	Tgfb1 (Rat) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0153	TGFB2 (Human) ELISA Kit	31.3 to 1000 pg/mL	6.6 pg/mL
KA1210	THBS1 (Human) ELISA Kit	1.56 to 100 ng/mL	0.78 ng/mL
KA0667	THBS2 (Human) ELISA Kit	156 to 10000 pg/mL	5 pg/mL
KA0511	Thrombin (Human) ELISA Kit	0.313 to 20 ng/mL	0.3 ng/mL
KA1131	Thromboxane A2 (Human) ELISA Kit	3.12 to 200 pg/mL	0.78 pg/mL
KA0294	Thromboxane B2 ELISA Kit	13.7 to 10000 pg/mL	10.54 pg/mL
KA1118	Thyroid peroxidase Ab ELISA Kit	33 to 3000 IU/mL	5 U/mL
KA1767	TIE1 (Human) ELISA Kit	N/A	0.6 ng/mL
KA0181	TIMP1 (Human) ELISA Kit	39 to 2500 pg/mL	10 pg/mL
KA1199	TIMP1 (Human) ELISA Kit	0.156 to 10 ng/mL	0.078 ng/mL
KA1814	Timp1 (Mouse) ELISA Kit	N/A	3 pg/mL
KA1815	Timp1 (Mouse) ELISA Kit	N/A	3 pg/mL
KA0437	Timp1 (Rat) ELISA Kit	31.2 to 2000 pg/mL	3 pg/mL
KA0549	TIMP2 (Human) ELISA Kit	156 to 10000 pg/mL	2 pg/mL
KA1816	Timp2 (Mouse) ELISA Kit	N/A	80 pg/mL
KA0420	TIMP3 (Human) ELISA Kit	156 to 10000 pg/mL	2 pg/mL
KA1768	TIMP4 (Human) ELISA Kit	N/A	25 pg/mL
KA0506	Tissue Factor (Human) ELISA Kit	25 to 800 pg/mL	20 pg/mL
KA1222	TLR2 (Human) ELISA Kit	0.47 to 30 ng/mL	0.12 ng/mL
KA1238	TLR4 (Human) ELISA Kit	0.312 to 20 ng/mL	N/A
KA1142	TLR5 (Human) ELISA Kit	0.312 to 20 ng/mL	0.156 ng/mL

ELISA Kits


Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA0191	TNF (Human) ELISA Kit	23 to 1500 pg/mL	5 pg/mL
KA0648	TNF (Human) ELISA Kit	0 to 500 pg/mL	7 pg/mL
KA0257	Tnf (Mouse) ELISA Kit	15.6 to 1000 pg/mL	3 pg/mL
KA1823	TNF (Pig) ELISA Kit	N/A	40 pg/mL
KA0280	Tnf (Rat) ELISA Kit	39 to 2500 pg/mL	11.2 pg/mL
KA0270	TNF-R I (Monkey) ELISA Kit	0.09 to 6 U/mL	0.05 U/mL
KA1770	TNFRSF10D (Human) ELISA Kit	N/A	10 pg/mL
KA1255	TNFRSF11A (Human) ELISA Kit	62.5 to 4000 pg/mL	2 pg/mL
KA2007	Tnfrsf11a (Mouse) ELISA Kit	62.5 to 4000 pg/mL	< 2 pg/mL
KA0030	TNFRSF11B (Human) ELISA Kit	1.5 to 60 pmol/L	0.13 pmol/L
KA0406	Tnfrsf11b (Mouse) ELISA Kit	93.8 to 6000 pg/mL	5 pg/mL
KA1719	TNFRSF18 (Human) ELISA Kit	N/A	15 pg/mL
KA1720	TNFRSF18 (Human) ELISA Kit	N/A	15 pg/mL
KA0114	TNFRSF1A (Human) ELISA Kit	0.08 to 5 ng/mL	0.053 ng/mL
KA0121	TNFRSF1A (Human) ELISA Kit	0.16 to 10 ng/mL	0.09 ng/mL
KA1254	TNFRSF1A (Human) ELISA Kit	7.8 to 500 pg/mL	1 pg/mL
KA0424	TNFRSF1B (Human) ELISA Kit	7.8 to 500 pg/mL	2 pg/mL
KA0145	TNFRSF8 (Human) ELISA Kit	1.6 to 100 ng/mL	0.33 ng/mL
KA1779	Tnfrsf8 (Mouse) ELISA Kit	N/A	10 pg/mL
KA0171	TNFRSF9 (Human) ELISA Kit	0.8 to 50 ng/mL	0.2 ng/mL
KA0176	TNFSF10 (Human) ELISA Kit	15.6 to 1000 pg/mL	5.4 pg/mL
KA0990	TNFSF10 (Human) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1808	Tnfsf11 (Mouse) ELISA Kit	N/A	1 pg/mL
KA1817	Tnfsf11 (Mouse) ELISA Kit	N/A	4 pg/mL
KA1818	Tnfsf12 (Mouse) ELISA Kit	N/A	0.3 pg/mL
KA0177	TNFSF13 (Human) ELISA Kit	0.8 to 50 ng/mL	0.4 ng/mL
KA1499	TNFSF13B (Human) ELISA Kit	62.5 to 4000 pg/mL	N/A
KA1500	Tnfsf13b (Mouse) ELISA Kit	46.9 to 3000 pg/mL	N/A
KA1741	TNFSF14 (Human) ELISA Kit	N/A	15 pg/mL
KA1721	TNFSF18 (Human) ELISA Kit	N/A	25 pg/mL
KA0174	TNFSF8 (Human) ELISA Kit	1.6 to 100 ng/mL	0.5 ng/mL
KA0550	Tnfsf8 (Mouse) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA1992	Tnni1/Tnni2 (Mouse) ELISA Kit	3.125 to 50 ng/mL	N/A
KA2010	Tnni1/Tnni2 (Rat) ELISA Kit	3.125 to 50 ng/mL	N/A
KA2055	TNNI2 (Monkey) ELISA Kit	3.125 to 200 ng/mL	N/A
KA2074	TNNI3 (Bovine) ELISA Kit	0.156 to 10 ng/mL	N/A
KA1991	TNNI3 (Dog) ELISA Kit	0.156 to 10 ng/mL	N/A
KA2067	Tnni3 (Guinea Pig) ELISA Kit	0.20 to 12.5 ng/mL	N/A
KA2054	TNNI3 (Monkey) ELISA Kit	0.156 to 10 ng/mL	N/A
KA2056	Tnni3 (Mouse) ELISA Kit	0.156 to 10 ng/mL	N/A
KA2057	Tnni3 (Mouse) ELISA Kit	0.078 to 5 ng/mL	N/A
KA2058	TNNI3 (Pig) ELISA Kit	0.039 to 2.5 ng/mL	N/A
KA2059	TNNI3 (Pig) ELISA Kit	0.039 to 2.5 ng/mL	N/A
KA2061	TNNI3 (Rabbit) ELISA Kit	0.156 to 10 ng/mL	N/A
KA1994	Tnni3 (Rat) ELISA Kit	0.156 to 10 ng/mL	N/A
KA1996	Tnni3 (Rat) ELISA Kit	0.078 to 5 ng/mL	N/A
KA0108	Total p53 (Human) ELISA Kit	1.56 to 100 U/mL	1.335 U/mL
KA0155	TP53 (Human) ELISA Kit	0.8 to 50 U/mL	0.33 U/mL
KA0110	TP53 (phospho S392) (Human) ELISA Kit	1.56 to 100 U/mL	1.413 U/mL
KA0109	TP53 (phospho S46) (Human) ELISA Kit	1.56 to 100 U/mL	1.450 U/mL
KA1769	TPO (Human) ELISA Kit	N/A	25 pg/mL
KA0418	Tpo (Mouse) ELISA Kit	31.2 to 2000 pg/mL	5 pg/mL
KA1772	TREM1 (Human) ELISA Kit	N/A	8 pg/mL
KA1929	Trf (Mouse) ELISA Kit	1.56 to 100 ng/mL	N/A
KA1977	TRF (Rabbit) ELISA Kit	6.25 to 100 ng/mL	N/A
KA0233	Troponin I (Human) ELISA Kit	0 to 75 ng/mL	1 ng/mL
KA1916	Tryptophan ELISA Kit	N/A	1.2 ug/mL
KA0196	TSH (Human) ELISA Kit	0 to 25 uIU/mL	0.2 uIU/mL
KA0197	TSH (Human) ELISA Kit	0 to 10 uIU/mL	0.07 uIU/mL
KA0515	TSH (Human) ELISA Kit	0.2 to 40 uIU/mL	0.2 uIU/mL
KA0516	TSH (Human) ELISA Kit	0.77 to 150 uIU/mL	0.77 uIU/mL
KA0517	TSH (Human) ELISA Kit	0.05 to 12 uIU/mL	0.05 uIU/mL
KA0495	TTR (Human) ELISA Kit	0.03 to 31.25 ng/mL	0.02 ng/mL
KA0535	TXN (Human) ELISA Kit	0.39 to 25 ng/mL	0.39 ng/mL
KA0534	TXNRD1 (Human) ELISA Kit	0.39 to 25 ng/mL	0.39 ng/mL




KA1255




KA2007




KA0406




KA0424




KA0196




KA0517



KA0428



KA0430



KA0666

ELISA Kits

Cat. Num.	Product Name	Calibration Range	Limit of Detection
KA0138	VCAM1 (Human) ELISA Kit	3.2 to 100 ng/mL	0.59 ng/mL
KA0428	Vcam1 (Mouse) ELISA Kit	156 to 10000 pg/mL	5 pg/mL
KA0650	VEGF (Human) ELISA Kit	9 to 700 pg/mL	5 pg/mL
KA0166	VEGFA (Human) ELISA Kit	15.6 to 1000 pg/mL	7.9 pg/mL
KA0258	Vegfa (Mouse) ELISA Kit	62.5 to 4000 pg/mL	20 pg/mL
KA0430	Vegfa (Rat) ELISA Kit	15.6 to 1000 pg/mL	1 pg/mL
KA0173	VEGFC (Human) ELISA Kit	0.23 to 15 ng/mL	0.057 ng/mL
KA0666	VEGFC (Human) ELISA Kit	31.2 to 2000 pg/mL	2 pg/mL
KA1170	Vitamin B12 (Human) ELISA Kit	15.6 to 1000 pmol/L	7.8 pmol/L
KA1168	Vitamin D (Human) ELISA Kit	0.78 to 50 ng/mL	N/A
KA1169	Vitamin D3 (Human) ELISA Kit	0.312 to 20 nmol/L	0.078 nmol/L
KA0512	VWF (Human) ELISA Kit	2.5 to 80 mU/mL	2.5 mU/mL
KA1742	XCL1 (Human) ELISA Kit	N/A	150 pg/mL
KA1800	Xcl1 (Mouse) ELISA Kit	N/A	40 pg/mL
KA0314	11-dehydro-TXB2 ELISA Kit	9.8 to 10000 pg/mL	4.31 pg/mL
KA0306	12(S)-HETE ELISA Kit	195 to 50000 pg/mL	146 pg/mL
KA0317	13(S)-HODE ELISA Kit	3.9 to 1000 ng/mL	1.6 ng/mL
KA0307	15(S)-HETE ELISA Kit	78.1 to 20000 pg/mL	69.21 pg/mL
KA0303	15-deoxy-Delta12,14-PGJ2 ELISA Kit	195 to 200000 pg/mL	36.8 pg/mL
KA1249	17-hydroxycorticosteroid (Human) ELISA Kit	0.625 to 40 nmol/mL	0.312 nmol/mL
KA1196	17-Ketosteroid (Human) ELISA Kit	78 to 5000 ng/mL	19.5 ng/mL
KA1896	17-OH-Progesterone ELISA Kit	N/A	3.6 pg/mL
KA1177	2,3-DPG (Human) ELISA Kit	0.625 to 40 nmol/mL	0.156 nmol/mL
KA1897	5a-Androstane-3a, 17b-Diol Glucuronide ELISA Kit	N/A	0.1 ng/mL
KA1881	5-hydroxyindoleacetic acid ELISA Kit	2.4 to 24.3 mg/L	0.17 mg/L
KA1229	6-keto-PGF1-alpha (Human) ELISA Kit	12.5 to 800 pg/mL	6.25 pg/mL
KA0295	6-keto-PGF1-alpha ELISA Kit	3.2 to 50000 pg/mL	1.40 pg/mL
KA0443	8-Epi-PGF2-alpha ELISA Kit	0.1 to 100 ng/mL	0.1 ng/mL
KA0444	8-hydroxy-2'-deoxyguanosine ELISA Kit	0.5 to 200 ng/mL	N/A
KA0298	8-iso-PGF2-alpha ELISA Kit	6.1 to 100000 pg/mL	16.3 pg/mL
KA0313	8-iso-PGF2-alpha ELISA Kit	160 to 100000 pg/mL	103.2 pg/mL
KA1221	8-OHdG (Human) ELISA Kit	156 to 10000 pg/mL	39 pg/mL

Disease ELISA Kits

Cat. Num.	Product Name	Suitable Samples	Label
KA1268	ACPA ELISA Kit	Serum, Plasma	HRP-conjugate
KA1048	ADIPOQ (Dog) ELISA Kit	Plasma, Serum, Tissue Culture Medium, Other Biological Media	HRP-conjugate
KA0016	ADIPOQ (Human) ELISA Kit	Plasma (Citrate, EDTA, Heparin), Serum	HRP-conjugate
KA0017	ADIPOQ (Human) ELISA Kit	Cell Culture Supernatant, Cerebrospinal Fluid, Milk, Plasma (Citrate, EDTA, Heparin), Serum, Urine	HRP-conjugate
KA0089	ADIPOQ (Human) ELISA Kit	Plasma, Serum	HRP-conjugate
KA0978	ADIPOQ (Human) ELISA Kit	Cell Culture Supernatant, Serum, Plasma, Urine, Milk, Saliva	Biotin-conjugate
KA0018	Adipoq (Mouse) ELISA Kit	Cell Culture Supernatant, Serum	HRP-conjugate
KA0090	Adipoq (Mouse) ELISA Kit	Serum, Plasma, Tissue Culture Medium	HRP-conjugate
KA1025	Adipoq (Mouse) ELISA Kit	Cell Culture Supernatant, Plasma, Serum, Urine	HRP-conjugate
KA1026	Adipoq (Rat) ELISA Kit	Cell Culture Supernatant, Plasma, Serum, Urine	HRP-conjugate
KA1047	Adipoq (Rat) ELISA Kit	Plasma, Serum, Tissue Culture Medium, Other Biological Media	HRP-conjugate
KA1087	Alpha Fodrin Ab IgG/IgA ELISA Kit	Serum, Plasma	HRP-conjugate
KA1078	AMA-M2 Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0926	Amphetamine ELISA Kit	Whole Blood, Oral Fluid, Serum, Plasma, Urine	HRP-conjugate
KA1080	ANA Hep Screen ELISA Kit	Serum, Plasma	HRP-conjugate
KA1079	ANA Profile ELISA Kit	Serum, Plasma	HRP-conjugate
KA0939	ANA Screen ELISA Kit	Plasma, Serum	HRP-conjugate
KA1081	ANA Screen ELISA Kit	Serum, Plasma	HRP-conjugate
KA1082	ANCA Profile ELISA Kit	Serum, Plasma	HRP-conjugate
KA1083	ANCA Screen ELISA Kit	Serum, Plasma	HRP-conjugate
KA1084	ANCA-C (PR3) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1086	ANCA-P (MPO) ELISA Kit	Serum, Plasma	HRP-conjugate
KA0953	Anthrax protective antigen IgG ELISA Kit	Serum	HRP-conjugate
KA1075	ASCA Combi ELISA Kit	Serum	HRP-conjugate
KA1076	ASCA IgA ELISA Kit	Serum	HRP-conjugate

Disease ELISA Kits

Cat. Num.	Product Name	Suitable Samples	Label
KA1077	ASCA IgG ELISA Kit	Serum	HRP-conjugate
KA1270	ASCA IgG/IgA ELISA Kit	Serum, Plasma	HRP-conjugate
KA1413	Atenolol ELISA Kit	Urine	AP-conjugate
KA0927	Barbiturate ELISA Kit	Whole Blood, Serum, Plasma, Urine	HRP-conjugate
KA1273	beta2-Glycoprotein I Ab Screen ELISA Kit	Serum, Plasma	HRP-conjugate
KA1271	beta2-Glycoprotein I IgA ELISA Kit	Serum, Plasma	HRP-conjugate
KA1272	beta2-Glycoprotein I IgG/IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA2092	Blastocystis spp. ELISA Kit	Stool	HRP-conjugate
KA2091	Bordetella pertussis IgA/IgM ELISA Kit	Serum	HRP-conjugate
KA2090	Bordetella pertussis IgG ELISA Kit	Serum	HRP-conjugate
KA1282	Borrelia burgdorferi IgG ELISA Kit	Serum, Plasma	HRP-conjugate
KA1283	Borrelia burgdorferi IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA1465	Borrelia IgG ELISA Kit	Serum	HRP-conjugate
KA1466	Borrelia IgG ELISA Kit	Serum, Cerebrospinal Fluid, Synovial Fluid	HRP-conjugate
KA1467	Borrelia IgM ELISA Kit	Serum	HRP-conjugate
KA1468	Borrelia IgM ELISA Kit	Serum, Cerebrospinal Fluid, Synovial Fluid	HRP-conjugate
KA1088	BPI Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0954	Brucella IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0955	Brucella IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA1274	C1q Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0928	Cannabinoids ELISA Kit	Whole Blood, Serum, Plasma, Urine	HRP-conjugate
KA0940	Cardiolipin IgA ELISA Kit	Plasma, Serum	HRP-conjugate
KA1096	Cardiolipin IgA ELISA Kit	Serum, Plasma	HRP-conjugate
KA0941	Cardiolipin IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0942	Cardiolipin IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA0943	Cardiolipin Total Ab ELISA Kit	Plasma, Serum	HRP-conjugate
KA1265	Cathepsin G Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1275	Centromere B Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1477	CFH IgG ELISA Kit	Serum	HRP-conjugate
KA2086	Chlamydia IgA ELISA Kit	Serum	HRP-conjugate
KA2084	Chlamydia IgG ELISA Kit	Serum	HRP-conjugate
KA2085	Chlamydia IgM ELISA Kit	Serum	HRP-conjugate
KA1464	Chlamydia pneumoniae IgA ELISA Kit	Serum	HRP-conjugate
KA2081	Chlamydia pneumoniae IgA ELISA Kit	Serum	HRP-conjugate
KA2083	Chlamydia pneumoniae IgA ELISA Kit	Serum	HRP-conjugate
KA0956	Chlamydia pneumoniae IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA2079	Chlamydia pneumoniae IgG ELISA Kit	Serum	HRP-conjugate
KA2082	Chlamydia pneumoniae IgG ELISA Kit	Serum	HRP-conjugate
KA0957	Chlamydia pneumoniae IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA2080	Chlamydia pneumoniae IgM ELISA Kit	Serum	HRP-conjugate
KA0958	Chlamydia trachomatis IgA ELISA Kit	Plasma, Serum	HRP-conjugate
KA2078	Chlamydia trachomatis IgA ELISA Kit	Serum	HRP-conjugate
KA0959	Chlamydia trachomatis IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA2077	Chlamydia trachomatis IgG ELISA Kit	Serum	HRP-conjugate
KA0960	Chlamydia trachomatis IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA0098	CML Ab (Human) ELISA Kit	Plasma, Serum	HRP-conjugate
KA0097	CML Ab (Mouse) ELISA Kit	Plasma, Serum	HRP-conjugate
KA0099	CML Ab (Rat) ELISA Kit	Plasma, Serum	HRP-conjugate
KA0929	Cocaine Metabolite ELISA Kit	Whole Blood, Serum, Plasma, Urine	HRP-conjugate
KA0930	Cotinine ELISA Kit	Serum, Urine	HRP-conjugate
KA1416	Cotinine ELISA Kit	Urine	AP-conjugate
KA1259	C-Peptide ELISA Kit	Serum, Plasma, Urine	HRP-conjugate
KA2094	Cryptosporidium spp. ELISA Kit	Stool	HRP-conjugate
KA1453	Cytomegalovirus IgA ELISA Kit	Serum, cerebrospinal fluid	HRP-conjugate
KA0227	Cytomegalovirus IgG ELISA Kit	Serum	HRP-conjugate
KA1452	Cytomegalovirus IgG ELISA Kit	Serum, cerebrospinal fluid	HRP-conjugate
KA0228	Cytomegalovirus IgM ELISA Kit	Serum	HRP-conjugate
KA0961	Dengue virus IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0962	Dengue virus IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA1097	DNase ELISA Kit	Serum, Plasma	HRP-conjugate
KA1098	dsDNA IgA ELISA Kit	Serum, Plasma	HRP-conjugate
KA0944	dsDNA IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA1100	dsDNA IgG ELISA Kit	Serum, Plasma	HRP-conjugate
KA1099	dsDNA IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA1101	dsDNA Screen ELISA Kit	Serum, Plasma	HRP-conjugate

Disease ELISA Kits

Cat. Num.	Product Name	Suitable Samples	Label
KA1102	Elastase Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1103	ENA Combi ELISA Kit	Serum, Plasma	HRP-conjugate
KA1106	ENA Screen ELISA Kit	Serum, Plasma	HRP-conjugate
KA1104	ENA-4 Profile ELISA Kit	Serum, Plasma	HRP-conjugate
KA1105	ENA-6 Profile ELISA Kit	Serum, Plasma	HRP-conjugate
KA1450	Epstein Barr Virus EA(D) IgG ELISA Kit	Serum, Plasma	HRP-conjugate
KA1451	Epstein Barr Virus EA(D) IgM ELISA Kit	Serum	HRP-conjugate
KA1446	Epstein Barr Virus EBNA-1 IgG ELISA Kit	Serum, Plasma	HRP-conjugate
KA1448	Epstein Barr Virus EBNA-1 IgG ELISA Kit	Serum	HRP-conjugate
KA1449	Epstein Barr Virus EBNA-1 IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA0653	Epstein Barr Virus IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0654	Epstein Barr Virus IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA0963	Epstein Barr Virus VCA IgA ELISA Kit	Plasma, Serum	HRP-conjugate
KA1443	Epstein Barr Virus VCA IgG ELISA Kit	Serum	HRP-conjugate
KA1444	Epstein Barr Virus VCA IgG ELISA Kit	Serum, Cerebrospinal Fluid	HRP-conjugate
KA1445	Epstein Barr Virus VCA IgG ELISA Kit	Serum	HRP-conjugate
KA0932	Fentanyl (Canine) ELISA Kit	Whole Blood, Oral Fluid, Serum, Plasma, Urine	HRP-conjugate
KA0931	Fentanyl (Human) ELISA Kit	Whole Blood, Serum, Plasma, Urine	HRP-conjugate
KA0933	Flunitrazepam ELISA Kit	Whole Blood, Serum, Plasma, Urine	HRP-conjugate
KA1266	GBM Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1863	Ghrelin (Human/Mouse/Rat) ELISA Kit	Serum, Plasma, Cell Culture Supernatant, Cell Lysate	Biotin-conjugate
KA2093	Giardia lamblia ELISA Kit	Stool	HRP-conjugate
KA2095	Giardia lamblia/Cryptosporidium spp. ELISA Kit	Stool	HRP-conjugate
KA1286	Gliadin IgA ELISA Kit	Serum, Plasma	HRP-conjugate
KA1285	Gliadin IgG ELISA Kit	Serum, Plasma	HRP-conjugate
KA1287	Gliadin IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA0290	HBeAg and HBeAg Ab ELISA Kit	Serum, Plasma (EDTA, Citrate, Heparin)	HRP-conjugate
KA0964	Helicobacter pylori IgA ELISA Kit	Plasma, Serum	HRP-conjugate
KA0219	Helicobacter pylori IgG ELISA Kit	Serum	HRP-conjugate
KA0220	Helicobacter pylori IgG ELISA Kit	Serum	HRP-conjugate
KA0221	Helicobacter pylori IgM ELISA Kit	Serum	HRP-conjugate
KA0284	Hepatitis A virus Ab ELISA Kit	Plasma, Serum	HRP-conjugate
KA0285	Hepatitis A virus IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA0287	Hepatitis B surface antigen Ab ELISA Kit	Plasma, Serum	HRP-conjugate
KA0286	Hepatitis B surface antigen ELISA Kit	Plasma, Serum	HRP-conjugate
KA0288	Hepatitis B virus core antigen Ab ELISA Kit	Plasma, Serum	HRP-conjugate
KA0289	Hepatitis B virus IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA0291	Hepatitis C virus Ab ELISA Kit	Plasma, Serum	HRP-conjugate
KA0450	Hepatitis E virus IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA1009	Hepatitis E virus IgM (Human) ELISA Kit	Serum, Plasma	HRP-conjugate
KA0451	Hepatitis E virus IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA1046	Herpes Simplex Virus I IgG (Mouse/Rat) ELISA Kit	Serum	HRP-conjugate
KA0229	Herpes Simplex Virus I IgG ELISA Kit	Serum	HRP-conjugate
KA0230	Herpes Simplex Virus I IgM ELISA Kit	Serum	HRP-conjugate
KA0651	Herpes Simplex Virus I/II IgG ELISA Kit	Serum, Plasma (EDTA, Heparin)	HRP-conjugate
KA1454	Herpes Simplex Virus I/II IgG ELISA Kit	Serum, Cerebrospinal Fluid	HRP-conjugate
KA0652	Herpes Simplex Virus I/II IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA0231	Herpes Simplex Virus II IgG ELISA Kit	Serum	HRP-conjugate
KA0232	Herpes Simplex Virus II IgM ELISA Kit	Serum	HRP-conjugate
KA1457	HHV-6 IgG ELISA Kit	Serum	HRP-conjugate
KA1120	Histone Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0452	HIV-1 p24 ELISA Kit	Tissue, Cell Culture Supernatant	HRP-conjugate
KA0292	HIV-1/2 Ab ELISA Kit	Plasma, Serum	HRP-conjugate
KA0453	HIV-1/2 Ab ELISA Kit	Plasma, Serum	HRP-conjugate
KA1244	HPV (Human) ELISA Kit	Serum, Tissue Homogenate, Other Biological Fluid	HRP-conjugate
KA1245	HPV IgG (Human) ELISA Kit	Cell Culture Supernate, Serum, Plasma	HRP-conjugate
KA0921	Insulin (Human) ELISA Kit	Plasma, Serum	HRP-conjugate
KA1107	Insulin Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1089	Intrinsic Factor Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1276	Jo-1 Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1267	Lactoferrin Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0025	LEP (Human) ELISA Kit	Cell Culture Supernatant, Plasma (Citrate, EDTA, Heparin), Serum	HRP-conjugate
KA0649	LEP (Human) ELISA Kit	Cell Culture Supernatant, Plasma, Serum, Breast Milk	HRP-conjugate
KA0994	Lep (Mouse) ELISA Kit	Sera, Plasma, Body Fluid, Tissue Lysate, Cell Culture Supernatant	Biotin-conjugate
KA0026	Lep (Mouse/Rat) ELISA Kit	Cell Culture Supernatant, Plasma (Citrate, EDTA, Heparin), Serum	HRP-conjugate

Disease ELISA Kits

Cat. Num.	Product Name	Suitable Samples	Label
KA0027	LEPR (Human) ELISA Kit	Cell Culture Supernatant, Plasma (Citrate, EDTA, Heparin), Serum	HRP-conjugate
KA1280	LKM-1 Ab ELISA Kit	Serum	HRP-conjugate
KA1090	Lysozyme Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1206	MDA (Human) ELISA Kit	Serum, Plasma, Other Biological Fluid	HRP-conjugate
KA0655	Measles IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0656	Measles IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA0934	Methamphetamine ELISA Kit	Whole Blood, Oral Fluid, Serum, Plasma, Urine	HRP-conjugate
KA0945	Mitochondrial IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0936	Morphine (Canine) ELISA Kit	Whole Blood, Serum, Plasma, Urine	HRP-conjugate
KA0935	Morphine ELISA Kit	Whole Blood, Serum, Plasma, Urine	HRP-conjugate
KA0657	Mumps IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0658	Mumps IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA1463	Mycoplasma pneumoniae IgA ELISA Kit	Serum	HRP-conjugate
KA2089	Mycoplasma pneumoniae IgA ELISA Kit	Serum	HRP-conjugate
KA0659	Mycoplasma pneumoniae IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA2087	Mycoplasma pneumoniae IgG ELISA Kit	Serum	HRP-conjugate
KA0660	Mycoplasma pneumoniae IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA2088	Mycoplasma pneumoniae IgM ELISA Kit	Serum	HRP-conjugate
KA1091	Nucleosome Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0937	Opiates ELISA Kit	Whole Blood, Oral Fluid, Serum, Plasma, Urine	HRP-conjugate
KA1260	Ovarian Ab ELISA Kit	Serum	HRP-conjugate
KA1281	PAPPA ELISA Kit	Serum, Plasma	HRP-conjugate
KA1108	Parietal Cell Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0938	PCP ELISA Kit	Whole Blood, Oral Fluid, Serum, Plasma, Urine	HRP-conjugate
KA1269	Phosphatidic Acid IgG/IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA1092	Phosphatidyl Inositol IgG/IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA1277	Phosphatidyl Serine Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1109	Phospholipid Screen IgG/IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA1085	PR3 IgG ELISA Kit	Serum, Plasma	HRP-conjugate
KA1685	Proglucagon (Human, Mouse, Rat) ELISA Kit	Serum, Plasma, Culture Supernatant, Cell Lysate	Biotin-conjugate
KA0922	Pro-insulin (Human) ELISA Kit	Plasma, Serum	HRP-conjugate
KA1094	Prothrombin Ab IgA ELISA Kit	Serum, Plasma	HRP-conjugate
KA1093	Prothrombin Ab IgG/IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA1095	Prothrombin Screen ELISA Kit	Serum, Plasma	HRP-conjugate
KA0034	RETN (Human) ELISA Kit	Cell Culture Supernatant, Plasma (Citrate, EDTA, Heparin), Serum	HRP-conjugate
KA0192	RETN (Human) ELISA Kit	Cell Culture Supernatant, Serum, Plasma (EDTA, Citrate, Heparin)	Biotin-conjugate
KA0035	Retn (Mouse) ELISA Kit	Cell Culture Supernatant, Plasma (Citrate, EDTA), Serum	HRP-conjugate
KA0995	Retn (Mouse) ELISA Kit	Sera, Plasma, Body Fluid, Tissue Lysate, Cell Culture Supernatant	Biotin-conjugate
KA0036	Retn (Rat) ELISA Kit	Cell Culture Supernatant, Serum	HRP-conjugate
KA1687	RETNLB (Human/Mouse/Rat) ELISA Kit	Serum, Plasma, Culture Supernatant, Cell Lysate	Biotin-conjugate
KA1442	Rheumatoid Factor Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1289	Rheumatoid factor IgA ELISA Kit	Serum, Plasma	HRP-conjugate
KA1288	Rheumatoid Factor IgG ELISA Kit	Serum, Plasma	HRP-conjugate
KA1290	Rheumatoid factor IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA1110	Rib-P Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1278	RNP 70 Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1111	RNP/Sm Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0946	RNP/Sm IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0223	Rubella virus IgG ELISA Kit	Serum	HRP-conjugate
KA1458	Rubella virus IgG ELISA Kit	Serum	HRP-conjugate
KA0224	Rubella virus IgM ELISA Kit	Serum	HRP-conjugate
KA1459	Rubella virus IgM ELISA Kit	Serum	HRP-conjugate
KA0965	Salmonella IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0966	Salmonella IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA0947	Sci-70 IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA1112	Sm Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0948	Sm IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0949	SS-A (Ro) IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA1113	SS-A 52 Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1114	SS-A 60 Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1115	SS-A Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA0950	SS-B (La) (Human) IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA1116	SS-B Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1117	ss-DNA Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1425	Sulfamethazine ELISA Kit	Swine Sample	AP-conjugate

Disease ELISA Kits

Cat. Num.	Product Name	Suitable Samples	Label
KA1008	Syphilis (Human) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1279	Thrombosis Profile IgG/IgM ELISA Kit	Serum, Plasma	HRP-conjugate
KA1121	Thyroglobulin (Human) ELISA Kit	Serum, Plasma	HRP-conjugate
KA0951	Thyroglobulin Ab ELISA Kit	Plasma, Serum	HRP-conjugate
KA1119	Thyroglobulin Ab ELISA Kit	Serum, Plasma	HRP-conjugate
KA1899	Thyroglobulin IgG ELISA Kit	Serum	HRP-conjugate
KA1469	Tick borne encephalitis virus IgG ELISA Kit	Serum	HRP-conjugate
KA1470	Tick borne encephalitis virus IgG ELISA Kit	Serum	HRP-conjugate
KA1471	Tick borne encephalitis virus IgM ELISA Kit	Serum	HRP-conjugate
KA1263	Tissue Transglutaminase IgG ELISA Kit	Serum, Plasma	HRP-conjugate
KA1264	Tissue Transglutaminase Screen ELISA Kit	Serum, Plasma	HRP-conjugate
KA0967	Toxoplasma gondii IgA ELISA Kit	Plasma, Serum	HRP-conjugate
KA1462	Toxoplasma gondii IgA ELISA Kit	Serum	HRP-conjugate
KA0225	Toxoplasma gondii IgG ELISA Kit	Serum	HRP-conjugate
KA1460	Toxoplasma gondii IgG ELISA Kit	Serum	HRP-conjugate
KA0226	Toxoplasma gondii IgM ELISA Kit	Serum	HRP-conjugate
KA1461	Toxoplasma gondii IgM ELISA Kit	Serum	HRP-conjugate
KA1900	TPO IgG ELISA Kit	Serum	HRP-conjugate
KA0968	Treponema pallidum IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA0969	Treponema pallidum IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA0970	Treponema pallidum Total ELISA Kit	Plasma, Serum	HRP-conjugate
KA1474	Tubb3 Ab (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1472	TUBB3 IgG (Human) ELISA Kit	Serum, Cerebrospinal Fluid	HRP-conjugate
KA1473	TUBB3 IgM (Human) ELISA Kit	Serum, Cerebrospinal Fluid	HRP-conjugate
KA0661	Varicella zoster virus IgG ELISA Kit	Plasma, Serum	HRP-conjugate
KA1455	Varicella zoster virus IgG ELISA Kit	Serum	HRP-conjugate
KA1456	Varicella zoster virus IgG ELISA Kit	Serum, Cerebrospinal Fluid	HRP-conjugate
KA0662	Varicella zoster virus IgM ELISA Kit	Plasma, Serum	HRP-conjugate
KA1475	VIM IgG (Human) ELISA Kit	Serum	HRP-conjugate
KA1476	VIM IgM (Human) ELISA Kit	Serum	HRP-conjugate

Immunoglobulin ELISA Kits

Cat. Num.	Product Name	Suitable Samples	Label
KA1010	HAMA (Human) ELISA Kit	Serum	HRP-conjugate
KA2028	IgA (Bovine) ELISA Kit	Biological Sample	HRP-conjugate
KA2031	IgA (Chicken) ELISA Kit	Biological Sample	HRP-conjugate
KA1966	IgA (Dog) ELISA Kit	Biological Sample	HRP-conjugate
KA2032	IgA (Guinea Pig) ELISA Kit	Biological Sample	HRP-conjugate
KA2035	IgA (Horse) ELISA Kit	Biological Sample	HRP-conjugate
KA1971	IgA (Monkey) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1943	IgA (Mouse) ELISA Kit	Serum, Plasma	HRP-conjugate
KA2038	IgA (Pig) ELISA Kit	Biological Sample	HRP-conjugate
KA2050	IgA (Rabbit) ELISA Kit	Biological Sample	HRP-conjugate
KA1950	IgA (Rat) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1967	IgD (Human) ELISA Kit	Biological Sample	HRP-conjugate
KA1940	IgE (Cat) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1937	IgE (Dog) ELISA Kit	Serum, Plasma	HRP-conjugate
KA2013	IgE (Horse) ELISA Kit	Biological Sample	HRP-conjugate
KA0216	IgE (Human) ELISA Kit	Serum	HRP-conjugate
KA2024	IgE (Monkey) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1944	IgE (Mouse) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1951	IgE (Rat) ELISA Kit	Serum, Plasma	HRP-conjugate
KA2029	IgG (Bovine) ELISA Kit	Biological Sample	HRP-conjugate
KA2046	IgG (Bovine) ELISA Kit	Biological Sample	HRP-conjugate
KA2048	IgG (Cat) ELISA Kit	Biological Sample	HRP-conjugate
KA2043	IgG (Chicken) ELISA Kit	Serum, Plasma, Egg-yolk Extract	HRP-conjugate
KA1938	IgG (Dog) ELISA Kit	Biological Sample	HRP-conjugate
KA2033	IgG (Guinea Pig) ELISA Kit	Biological Sample	HRP-conjugate
KA2036	IgG (Horse) ELISA Kit	Biological Sample	HRP-conjugate
KA1729	IgG (Human) ELISA Kit	Serum, Plasma	Biotin-conjugate

Immunoglobulin ELISA Kits

Cat. Num.	Product Name	Suitable Samples	Label
KA1963	IgG (Monkey) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1945	IgG (Mouse) ELISA Kit	Serum, Plasma	HRP-conjugate
KA2016	IgG (Pig) ELISA Kit	Biological Sample	HRP-conjugate
KA2017	IgG (Rabbit) ELISA Kit	Biological Sample	HRP-conjugate
KA1952	IgG (Rat) ELISA Kit	Serum, Plasma	HRP-conjugate
KA2042	IgG (Turkey) ELISA Kit	Serum, Plasma, Egg-yolk Extract	HRP-conjugate
KA1730	IgG1 (Human) ELISA Kit	Serum, Plasma	Biotin-conjugate
KA2025	IgG1 (Monkey) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1011	IgG1 (Mouse) ELISA Kit	Culture Supernatant, Serum	HRP-conjugate
KA1013	IgG2a (Mouse) ELISA Kit	Culture Supernatant, Serum	HRP-conjugate
KA1012	IgG2b (Mouse) ELISA Kit	Culture Supernatant, Serum	HRP-conjugate
KA1970	IgG2c (Mouse) ELISA Kit	Biological Sample	HRP-conjugate
KA1854	IgG3 (Human) ELISA Kit	Cell Culture Supernatants, Milk, Plasma, Saliva, Serum, Urine	Biotin-conjugate
KA1946	IgG3 (Mouse) ELISA Kit	Serum, Plasma	HRP-conjugate
KA2030	IgM (Bovine) ELISA Kit	Biological Sample	HRP-conjugate
KA1956	IgM (Chicken) ELISA Kit	Serum, Plasma	HRP-conjugate
KA2045	IgM (Chicken) ELISA Kit	Serum, Plasma, Egg-yolk Extract	HRP-conjugate
KA1958	IgM (Dog) ELISA Kit	Biological Sample	HRP-conjugate
KA2034	IgM (Guinea Pig) ELISA Kit	Biological Sample	HRP-conjugate
KA2037	IgM (Horse) ELISA Kit	Biological Sample	HRP-conjugate
KA1855	IgM (Human) ELISA Kit	Cell Culture Supernatants, Milk, Plasma, Saliva, Serum, Urine	Biotin-conjugate
KA1962	IgM (Monkey) ELISA Kit	Serum, Plasma	HRP-conjugate
KA1014	IgM (Mouse) ELISA Kit	Culture Supernatant, Serum	HRP-conjugate
KA2039	IgM (Pig) ELISA Kit	Biological Sample	HRP-conjugate
KA2040	IgM (Rabbit) ELISA Kit	Biological Sample	HRP-conjugate
KA1953	IgM (Rat) ELISA Kit	Serum, Plasma	HRP-conjugate
KA2044	IgM (Turkey) ELISA Kit	Serum, Plasma, Egg-yolk Extract	HRP-conjugate
KA1957	IgY (Chicken) ELISA Kit	Serum, Plasma	HRP-conjugate


Environmental Risk ELISA Kits

Cat. Num.	Product Name	Suitable Samples	Label
KA1480	Acetnaftalen Ab (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1481	Antracen IgG (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1414	Atrazine ELISA Kit	Water sample leach from soil	AP-conjugate
KA1490	Benzo(a)pyrene Ab (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1493	Benzo(a)pyrene ELISA Kit	Soil, Water, Sediments	HRP-conjugate
KA1494	Benzo(a)pyrene ELISA Kit	Water sample	HRP-conjugate
KA1489	Benzo(a)pyrene IgA (Human) ELISA Kit	Serum	HRP-conjugate
KA1488	Benzo(a)pyrene IgG (Human) ELISA Kit	Serum	HRP-conjugate
KA1495	Bisphenol A ELISA Kit	Soil, Water	HRP-conjugate
KA1491	BPDE-DNA Ab (Human) ELISA Kit	Serum	HRP-conjugate
KA1492	BPDE-DNA Ab (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1485	Chrysene IgG (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1483	Fluoranthene IgG (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1484	Fluorene IgG (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1496	Microcystin-LR ELISA Kit	Water	HRP-conjugate
KA1486	Naphthalene IgG (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1424	Paraquat ELISA Kit	Plasma	HRP-conjugate
KA1482	Phenanthrene IgG (Mouse) ELISA Kit	Serum	HRP-conjugate
KA1487	Pyrene IgG (Mouse) ELISA Kit	Serum	HRP-conjugate


Food Safety ELISA Kits

Cat. Num.	Product Name	Suitable Samples	Label
KA1418	Fumonisin ELISA Kit	Rice, corn, cornmeal	HRP-conjugate
KA1422	Monensin ELISA Kit	Animal feeds	AP-conjugate
KA1423	Orchratoxin-A ELISA Kit	Grains, cereals, coffee, animal feeds	HRP-conjugate
KA1427	Tylosin ELISA Kit	Feed extract	AP-conjugate
KA1428	Zearalenone ELISA Kit	Grains, milk	AP-conjugate


Activity Assay Kits




KA0073




KA0066




KA0064



KA0058



KA0061




KA0057


Cat. Num.	Product Name
KA1607	Acetylcholinesterase Assay Kit
KA1497	Aggrecanase Activity Assay Kit
KA0070	AKT/PKB Kinase (Human) Assay/Inhibitor Screening Kit
KA1294	Alanine Transaminase Assay Kit
KA1295	Aldehyde Site Detection Kit
KA1608	alpha-Glucosidase Assay Kit
KA1622	Antioxidant Assay Kit
KA1296	Apoptotic Blebs Assay Kit
KA1609	Arginase Assay Kit
KA1297	Ascorbate Assay
KA1429	Asparaginase Activity Assay Kit
KA1625	Aspartate Transaminase Assay Kit
KA1298	ATF2 (phospho T69/T71) Transcription Factor Assay
KA0067	Aurora A Kinase (Human) Assay/Inhibitor Screening Kit
KA0074	Aurora Family Kinase (Human) Assay/Inhibitor Screening Kit
KA1299	Autophagy/Cytotoxicity Dual Staining Kit
KA1300	BACE Inhibitor Screening Assay Kit
KA1612	BCG Albumin Assay Kit
KA1613	BCP Albumin Assay Kit
KA1614	Bilirubin Assay Kit
KA0723	Calpain Activity Assay Kit
KA0073	CaM-kinase II (Human) Assay Kit
KA1626	Catalase Assay Kit
KA0066	Cdc2-Cyclin B Kinase (Human) Assay Kit
KA1301	Cell Cycle Phase Determination Kit
KA0101	Cellular Histone Acetylation (Human) Assay Kit
KA1302	CFSE Cell Division Assay Kit
KA0064	Checkpoint Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA1544	Chk 2 Gene Knockdown Quantification Kit
KA1303	Cholesterol Assay Kit
KA1304	Cholesterol Cell-Based Detection Assay Kit
KA1305	ChREBP Cell-Based Translocation Assay Kit
KA1306	ChREBP Transcription Factor Assay Kit
KA0071	CK2 Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA1615	Copper Assay Kit
KA0323	COX Activity Kit
KA1308	CREB1 (phospho S133) Transcription Factor Assay Kit
KA0058	c-Src Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA0063	Cyclic GMP dependent protein kinase (cGK) (Human) Assay Kit
KA1546	Cytochrome p450 Gene Knockdown Quantification Kit
KA1639	DNA Assay Kit
KA1440	DNA Methyltransferase 1 Activity/Inhibitor Screening Assay Core Kit
KA1604	DNA Methyltransferase 1 Activity/Inhibitor Screening Assay Core Kit
KA1547	DNA Methyltransferase 3B Activity/Inhibitor Screening Assay Core Kit
KA1515	DNMT Activity/Inhibition Assay Ultra Kit (Colorimetric)
KA1516	DNMT Activity/Inhibition Assay Ultra Kit (Fluorometric)
KA1311	DPP (IV) Inhibitor Screening Assay Kit
KA1545	E2F4 Gene Knockdown Quantification Kit
KA0061	EphA2 Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA0972	F10 (Human) Chromogenic Activity Assay Kit
KA0971	F7 (Human) Chromogenic Activity Assay Kit
KA0057	FGFR2 Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA1315	Formaldehyde Assay Kit
KA1542	General Gene Knockdown Quantification Kit
KA1536	Global Acetyl Histone H4K12 Quantification Kit (Colorimetric)
KA1537	Global Acetyl Histone H4K12 Quantification Kit (Fluorometric)
KA1534	Global Acetyl Histone H4K16 Quantification Kit (Colorimetric)
KA1535	Global Acetyl Histone H4K16 Quantification Kit (Fluorometric)
KA1530	Global Acetyl Histone H4K5 Quantification Kit (Colorimetric)
KA1531	Global Acetyl Histone H4K5 Quantification Kit (Fluorometric)
KA1532	Global Acetyl Histone H4K8 Quantification Kit (Colorimetric)
KA1533	Global Acetyl Histone H4K8 Quantification Kit (Fluorometric)
KA1519	Global Di-Methyl Histone H4K20 Quantification Kit (Colorimetric)
KA1520	Global Di-Methyl Histone H4K20 Quantification Kit (Fluorometric)
KA1517	Global Mono-Methyl Histone H4K20 Quantification Kit (Colorimetric)

Activity Assay Kits


Cat. Num.	Product Name
KA1518	Global Mono-Methyl Histone H4K20 Quantification Kit (Fluorometric)
KA1523	Global Pan-Methyl Histone H4K20 Quantification Kit (Colorimetric)
KA1524	Global Pan-Methyl Histone H4K20 Quantification Kit (Fluorometric)
KA1521	Global Tri-Methyl Histone H4K20 Quantification Kit (Colorimetric)
KA1522	Global Tri-Methyl Histone H4K20 Quantification Kit (Fluorometric)
KA1316	Glutathione S-Transferase Assay Kit
KA1317	Glycerol Cell-Based Assay Kit
KA1318	HAT Inhibitor Screening Assay Kit
KA1319	HDAC Cell-Based Assay Kit
KA1320	HDAC1 Inhibitor Screening Assay Kit
KA0107	HDAC8 Deacetylase Fluorometric (Human) Assay Kit
KA0104	HDACs Deacetylase Fluorometric (Human) Assay Kit
KA1617	Heme Assay Kit
KA1616	Hemoglobin Assay Kit
KA1321	HIF1A Transcription Factor Assay Kit
KA1441	Histone Demethylase LSD1 Inhibitor Screening Assay Core Kit
KA1548	Histone Demethylase LSD1 Inhibitor Screening Assay Core Kit
KA1322	Hydrogen Peroxide Cell-Based Assay Kit
KA1512	Hydroxymethylated DNA Quantification Kit (Colorimetric)
KA1513	Hydroxymethylated DNA Quantification Kit (Fluorometric)
KA0078	IKK-alpha/beta Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA1323	In Vitro Angiogenesis Assay Kit
KA1629	Invertase Assay Kit
KA1527	JARID Demethylase Activity/Inhibition Assay Kit (Colorimetric)
KA1529	JARID Demethylase Activity/Inhibition Assay Kit (Fluorometric)
KA1324	JC-1 Mitochondrial Membrane Potential Assay Kit
KA1526	JMJD2 Demethylase Activity/Inhibition Assay Kit (Colorimetric)
KA1630	Ketone Body Assay Kit
KA1631	Kinase Assay Kit
KA1327	LDL Uptake Cell-Based Assay Kit
KA1328	Lipid Hydroperoxide (LPO) Assay Kit
KA1329	Lipoxygenase Inhibitor Screening Assay Kit
KA1525	LSD1 Demethylase Activity/Inhibition Assay Kit (Colorimetric)
KA1603	LSD1 Demethylase Activity/Inhibition Assay Kit (Fluorometric)
KA1330	Lysosome/Cytotoxicity Dual Staining Kit
KA1331	Malachite Green Phosphate Assay Kit
KA1312	MAPK1 (phospho T202/Y204) Cell-Based Phosphorylation/Translocation Assay Kit
KA1347	MAPK14 (phospho T180/Y182) Cell-Based Phosphorylation/Translocation Assay Kit
KA0068	MAPKAP-kinase 2 (Human) Assay/Inhibitor Screening Assay Kit
KA0055	Met Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA1510	Methylated DNA Quantification Kit (Colorimetric)
KA1511	Methylated DNA Quantification Kit (Fluorometric)
KA1632	Monoamine Oxidase Assay Kit
KA1336	MPO Chlorination Assay Kit
KA1337	MPO Inhibitor Screening Assay Kit
KA1338	MPO Peroxidation Assay Kit
KA0079	Mps1/TTK Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA1334	MTT Cell Proliferation Assay Kit
KA1543	NCOR2 Gene Knockdown Quantification Kit
KA1633	Neuraminidase Assay Kit
KA1339	NF-kB (human p50) Transcription Factor Assay Kit
KA1340	NF-kB (human p50/p65) Combo Transcription Factor Assay Kit
KA1341	NF-kB (p65) Transcription Factor Assay Kit
KA1342	Nitrate/Nitrite Colorimetric Assay Kit
KA1343	Nitrate/Nitrite Colorimetric Assay Kit (LDH method)
KA1344	Nitrate/Nitrite Fluorometric Assay Kit
KA1345	NOS Activity Assay Kit
KA1346	Nuclear Extraction Kit
KA0077	p38 Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA0059	p56Lck Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA1354	PAF Acetylhydrolase Assay Kit
KA1355	PAF Acetylhydrolase Inhibitor Screening Assay Kit
KA0080	PKD1 Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA1620	Peroxidase Assay Kit
KA1618	Peroxide Assay Kit




KA0078




KA0055




KA0079



KA0077




KA0059




KA0080


Activity Assay Kits




KA0069




KA0075




KA0076



KA0056



KA1349



KA0060

Cat. Num.	Product Name
KA1307	Phospholipases A2 (calcium-dependent cytosolic) Assay Kit
KA0069	Pim-1 Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA0075	PKC Super Family Kinase (Human) Assay Kit
KA0974	PLAT (Human) Chromogenic Activity Assay Kit
KA0976	PLAU (Human) Chromogenic Activity Assay Kit (Direct Assay)
KA0977	PLAU (Human) Chromogenic Activity Assay Kit (Indirect Assay)
KA1637	pNPP Phosphatase Assay Kit
KA0065	Polo-like kinase 1 (Human) Assay/Inhibitor Screening Assay Kit
KA0076	Polo-like kinase 3 (Human) Assay/Inhibitor Screening Assay Kit
KA1356	PPARA Transcription Factor Assay Kit
KA1357	PPARA/PPARD/PPARG Complete Transcription Factor Assay Kit
KA1359	PPARD Transcription Factor Assay Kit
KA1358	PPARG Transcription Factor Assay Kit
KA1431	Proteasome Activity Assay Kit
KA1638	Protein Assay Kit
KA0086	Protein Phosphatase Cdc25 Combo Fluorometric (Human) Assay Kit
KA0083	Protein Phosphatase Cdc25A Fluorometric (Human) Assay Kit
KA0084	Protein Phosphatase Cdc25B Fluorometric (Human) Assay Kit
KA0085	Protein Phosphatase Cdc25C Fluorometric (Human) Assay Kit
KA0087	Protein Phosphatase Cdi1/KAP Fluorometric (Human) Assay Kit
KA0088	Protein Phosphatase LMW-PTP/ACP1 Fluorometric (Human) Assay Kit
KA0081	Protein Tyrosine Phosphatase 1B (PTP1B) Fluorometric (Human) Assay Kit
KA0056	Pyk2 Kinase (Human) Assay/Inhibitor Screening Assay Kit
KA1361	Renin Inhibitor Screening Assay Kit
KA0062	Rho-kinase (Human) Assay Kit
KA1362	Secreted Alkaline Phosphatase Reporter Gene Assay Kit (Luminescence)
KA1498	Sensitive Aggrecanase Activity Assay Kit
KA0973	SERPINE1 (Human) Chromogenic Activity Assay Kit
KA1366	SIRT1 Direct Fluorescent Screening Assay Kit
KA1367	SIRT1 FRET-Based Screening Assay Kit
KA0105	SIRT1/Sir2 Deacetylase Fluorometric (Human) Assay Kit
KA1368	SIRT2 Direct Fluorescent Screening Assay Kit
KA0106	SIRT3 Deacetylase Fluorometric (Human) Assay Kit
KA1369	SIRT3 Direct Fluorescent Screening Assay Kit
KA1370	SIRT6 Direct Fluorescent Screening Assay Kit
KA1373	Sphingomyelin Assay Kit
KA1374	Sphingomyelinase Assay Kit
KA1377	SREBF1 Transcription Factor Assay Kit
KA1378	SREBF2 Cell-Based Translocation Assay Kit
KA1379	SREBF2 Transcription Factor Assay Kit
KA1380	Steatosis Colorimetric Assay Kit
KA1621	Sulfate Assay Kit
KA0082	T Cell Protein Tyrosine Phosphatase (TC-PTP) Fluorometric (Human) Assay Kit
KA1381	TBARS Assay Kit
KA0975	TF (Human) Chromogenic Activity Assay Kit
KA1382	Thiol Detection Assay Kit
KA1538	Total Histone H3 Acetylation Detection Fast Kit (Colorimetric)
KA1539	Total Histone H3 Acetylation Detection Fast Kit (Fluorometric)
KA1540	Total Histone H4 Acetylation Detection Fast Kit (Colorimetric)
KA1541	Total Histone H4 Acetylation Detection Fast Kit (Fluorometric)
KA1383	Total Phosphatidic Acid Assay Kit
KA1348	TP53 Cell-Based Activation/Translocation Assay Kit
KA1349	TP53 Designer Transcription Factor Assay Kit
KA1350	TP53 Total and TP53 (phospho S392) Dual Staining Assay Kit
KA1351	TP53 Transcription Factor Assay Kit
KA0060	TrkA Kinase (Human) Assay/Inhibitor Screening Kit
KA0779	Trypsin Activity Assay Kit
KA1623	Urease Assay Kit
KA0218	Visual hCG Pregnancy Kit
KA0072	Wee1 Kinase (Human) Assay/Inhibitor Screening Kit
KA1384	WST-1 Cell Proliferation Assay Kit
KA1385	WST-8 Cell Proliferation Assay Kit
KA1386	Xanthine Oxidase Assay Kit
KA1387	XTT Cell Proliferation Assay Kit
KA1619	Zinc Assay Kit


Activity Assay Kits

Cat. Num.	Product Name
KA1352	20S Proteasome Assay Kit
KA1353	7-AAD Cell Viability Assay Kit
KA1293	7-AAD/CFSE Cell-Mediated Cytotoxicity Assay Kit


Apoptosis Assay Kits

Cat. Num.	Product Name	Cat. Num.	Product Name
KA0720	Apoptotic DNA Ladder Detection Kit	KA0750	Caspase-5 Inhibitor Screening Kit
KA0721	Apoptotic DNA Ladder Detection Kit (Premium)	KA0751	Caspase-6 Assay Kit (Colorimetric)
KA0722	Apoptotic DNA Ladder Isolation Kit	KA0752	Caspase-6 Assay Kit (Fluorometric)
KA0724	Caspase Screening Kit	KA0753	Caspase-6 Inhibitor Screening Kit
KA0725	Caspase Staining Kit (Green)	KA0754	Caspase-7 Immunoassay Kit
KA0726	Caspase Staining Kit (Red)	KA0755	Caspase-7 Inhibitor Screening Kit
KA0727	Caspase-1 Assay Kit (Colorimetric)	KA0756	Caspase-8 Assay Kit (Colorimetric)
KA0728	Caspase-1 Assay Kit (Fluorometric)	KA0757	Caspase-8 Assay Kit (Fluorometric)
KA0729	Caspase-1 Inhibitor Screening Kit	KA0758	Caspase-8 Inhibitor Screening Kit
KA0730	Caspase-10 Assay Kit (Colorimetric)	KA0759	Caspase-8 Staining Kit (Green)
KA0731	Caspase-10 Assay Kit (Fluorometric)	KA0760	Caspase-8 Staining Kit (Red)
KA0732	Caspase-10 Inhibitor Screening Kit	KA0761	Caspase-9 Assay Kit (Colorimetric)
KA0733	Caspase-12 Assay Kit (Fluorometric)	KA0762	Caspase-9 Assay Kit (Fluorometric)
KA0734	Caspase-12 Staining Kit (Green)	KA0763	Caspase-9 Inhibitor Screening Kit
KA0735	Caspase-2 Assay Kit (Colorimetric)	KA0764	Caspase-9 Staining Kit (Green)
KA0736	Caspase-2 Assay Kit (Fluorometric)	KA0765	Caspase-9 Staining Kit (Red)
KA0737	Caspase-2 Inhibitor Screening Kit	KA0766	Cathepsin B Activity Assay Kit
KA0738	Caspase-2 Staining Kit (Green)	KA0767	Cathepsin D Activity Assay Kit
KA0739	Caspase-3 Immunoassay Kit	KA0768	Cathepsin H Activity Assay Kit
KA0740	Caspase-3 Assay Kit (Colorimetric)	KA0769	Cathepsin K Activity Assay Kit
KA0741	Caspase-3 Assay Kit (Fluorometric)	KA0770	Cathepsin L Activity Assay Kit
KA0742	Caspase-3 Inhibitor Screening Kit	KA0771	Cathepsin S Activity Assay Kit
KA0743	Caspase-3 Staining Kit (Green)	KA0772	Cytochrome c Apoptosis Detection Kit
KA0744	Caspase-3 Staining Kit (Red)	KA0773	In Situ BrdU DNA Fragmentation Assay Kit
KA0745	Caspase-4 Assay Kit (Colorimetric)	KA0774	In Situ DNA Fragmentation Assay Kit
KA0746	Caspase-4 Assay Kit (Fluorometric)	KA0775	In Situ BrdU DNA Fragmentation Assay Kit (IHC)
KA0747	Caspase-4 Inhibitor Screening Kit	KA0776	In Situ BrdU DNA Fragmentation Assay Kit (Red)
KA0748	Caspase-5 Assay Kit (Colorimetric)	KA0778	Mitochondrial Apoptosis Detection Kit
KA0749	Caspase-5 Assay Kit (Fluorometric)	KA1335	Multi-Parameter Apoptosis Assay Kit


Apoptosis




Normal cell



Condensation



Fragmentation




Apoptotic bodies

Assay Kits

Cell Damage and Oxidative Stress Assay Kits


Cat. Num.	Product Name	Cat. Num.	Product Name
KA0784	DNA Damage Quantification Kit	KA0789	HDAC Inhibitor Drug Screening Kit
KA0782	HAT Activity Assay Kit	KA0785	LDH-Cytotoxicity Assay Kit
KA0788	HDAC Activity Assay Kit (Colorimetric)	KA0786	LDH-Cytotoxicity Assay Kit II
KA0787	HDAC Activity Assay Kit (Fluorometric)	KA0783	Superoxide Dismutase Activity Assay Kit



KA0783

Cell Proliferation and Cytotoxicity Assay Kits

Cat. Num.	Product Name	Cat. Num.	Product Name
KA0898	Bioluminescence Cytotoxicity Assay Kit	KA0897	Fluorescence Cell Viability Assay
KA1605	Blue Cell Viability Assay Kit	KA0901	Live-Dead Cell Staining Kit
KA0781	Cell Proliferation Assay Kit II	KA1606	MTT Cell Viability Assay Kit



KA0897


Cell Proliferation and Cytotoxicity Assay Kits

Cat. Num.	Product Name
KA0899	Senescence Detection Kit

Epigenetic Assay Kits


Cat. Num.	Product Name
KA0293	Cellular UV DNA-Damage Detection Kit
KA0677	DNA Demethylase Activity/Inhibition Assay Kit
KA0680	DNA Methyltransferase Activity/Inhibition Assay Kit (Colorimetric)
KA0681	DNA Methyltransferase Activity/Inhibition Assay Kit (Fluorometric)
KA0682	DNMT1 Assay Kit
KA0683	DNMT3A Assay Kit
KA0684	DNMT3B Assay Kit
KA0637	Global Acetyl Histone H3-K14 Quantification Kit (Colorimetric)
KA0638	Global Acetyl Histone H3-K14 Quantification Kit (Fluorometric)
KA0639	Global Acetyl Histone H3-K18 Quantification Kit (Colorimetric)
KA0640	Global Acetyl Histone H3-K18 Quantification Kit (Fluorometric)
KA0641	Global Acetyl Histone H3-K23 Quantification Kit (Colorimetric)
KA0642	Global Acetyl Histone H3-K23 Quantification Kit (Fluorometric)
KA0643	Global Acetyl Histone H3-K36 Quantification Kit (Colorimetric)
KA0644	Global Acetyl Histone H3-K36 Quantification Kit (Fluorometric)
KA0645	Global Acetyl Histone H3-K56 Quantification Kit (Colorimetric)
KA0646	Global Acetyl Histone H3-K56 Quantification Kit (Fluorometric)
KA0635	Global Acetyl Histone H3-K9 Quantification Kit (Colorimetric)
KA0636	Global Acetyl Histone H3-K9 Quantification Kit (Fluorometric)
KA0602	Global Di-Methyl Histone H3-K27 Quantification Kit (Colorimetric)
KA0603	Global Di-Methyl Histone H3-K27 Quantification Kit (Fluorometric)
KA0610	Global Di-Methyl Histone H3-K36 Quantification Kit (Colorimetric)
KA0611	Global Di-Methyl Histone H3-K36 Quantification Kit (Fluorometric)
KA0584	Global Di-Methyl Histone H3-K4 Quantification Kit (Colorimetric)
KA0585	Global Di-Methyl Histone H3-K4 Quantification Kit (Fluorometric)
KA0618	Global Di-Methyl Histone H3-K79 Quantification Kit (Colorimetric)
KA0619	Global Di-Methyl Histone H3-K79 Quantification Kit (Fluorometric)
KA0594	Global Di-Methyl Histone H3-K9 Quantification Kit (Colorimetric)
KA0595	Global Di-Methyl Histone H3-K9 Quantification Kit (Fluorometric)
KA0633	Global Histone H3 Acetylation Assay Kit
KA0703	Global Histone H3 Phosphorylation (Ser10) Assay Kit (Colorimetric)
KA0704	Global Histone H3 Phosphorylation (Ser10) Assay Kit (Fluorometric)
KA0705	Global Histone H3 Phosphorylation (Ser28) Assay Kit (Colorimetric)
KA0706	Global Histone H3 Phosphorylation (Ser28) Assay Kit (Fluorometric)
KA0582	Global Histone H3-K27 Methylation Assay Kit
KA0583	Global Histone H3-K27 Tri-Methylation Assay Kit
KA0580	Global Histone H3-K4 Methylation Assay Kit
KA0581	Global Histone H3-K9 Methylation Assay Kit
KA0634	Global Histone H4 Acetylation Assay Kit
KA0600	Global Mono-Methyl Histone H3-K27 Quantification Kit (Colorimetric)
KA0601	Global Mono-Methyl Histone H3-K27 Quantification Kit (Fluorometric)
KA0608	Global Mono-Methyl Histone H3-K36 Quantification Kit (Colorimetric)
KA0609	Global Mono-Methyl Histone H3-K36 Quantification Kit (Fluorometric)
KA0586	Global Mono-Methyl Histone H3-K4 Quantification Kit (Colorimetric)
KA0587	Global Mono-Methyl Histone H3-K4 Quantification Kit (Fluorometric)
KA0616	Global Mono-Methyl Histone H3-K79 Quantification Kit (Colorimetric)
KA0617	Global Mono-Methyl Histone H3-K79 Quantification Kit (Fluorometric)
KA0592	Global Mono-Methyl Histone H3-K9 Quantification Kit (Colorimetric)
KA0593	Global Mono-Methyl Histone H3-K9 Quantification Kit (Fluorometric)
KA0606	Global Pan-Methyl Histone H3-K27 Quantification Kit (Colorimetric)
KA0607	Global Pan-Methyl Histone H3-K27 Quantification Kit (Fluorometric)
KA0614	Global Pan-Methyl Histone H3-K36 Quantification Kit (Colorimetric)
KA0615	Global Pan-Methyl Histone H3-K36 Quantification Kit (Fluorometric)
KA0590	Global Pan-Methyl Histone H3-K4 Quantification Kit (Colorimetric)
KA0591	Global Pan-Methyl Histone H3-K4 Quantification Kit (Fluorometric)

Epigenetic




Epigenetic Assay Kits


Cat. Num.	Product Name
KA0622	Global Pan-Methyl Histone H3-K79 Quantification Kit (Colorimetric)
KA0623	Global Pan-Methyl Histone H3-K79 Quantification Kit (Fluorometric)
KA0598	Global Pan-Methyl Histone H3-K9 Quantification Kit (Colorimetric)
KA0599	Global Pan-Methyl Histone H3-K9 Quantification Kit (Fluorometric)
KA0604	Global Tri-Methyl Histone H3-K27 Quantification Kit (Colorimetric)
KA0605	Global Tri-Methyl Histone H3-K27 Quantification Kit (Fluorometric)
KA0612	Global Tri-Methyl Histone H3-K36 Quantification Kit (Colorimetric)
KA0613	Global Tri-Methyl Histone H3-K36 Quantification Kit (Fluorometric)
KA0588	Global Tri-Methyl Histone H3-K4 Quantification Kit (Colorimetric)
KA0589	Global Tri-Methyl Histone H3-K4 Quantification Kit (Fluorometric)
KA0620	Global Tri-Methyl Histone H3-K79 Quantification Kit (Colorimetric)
KA0621	Global Tri-Methyl Histone H3-K79 Quantification Kit (Fluorometric)
KA0596	Global Tri-Methyl Histone H3-K9 Quantification Kit (Colorimetric)
KA0597	Global Tri-Methyl Histone H3-K9 Quantification Kit (Fluorometric)
KA0626	HAT Activity/Inhibition Assay Kit (Colorimetric)
KA0628	HDAC Activity/Inhibition Assay Kit (Colorimetric)
KA0627	HDAC Activity/Inhibition Assay Kit (Fluorometric)
KA0629	HDAC1 Assay Kit
KA0688	HDAC2 Assay Kit
KA0631	HDAC8 Assay Kit
KA0573	Histone Demethylase (H3-K4 Specific) Activity/Inhibition Assay Kit
KA0576	Histone Demethylase (H3-K9 Specific) Activity/Inhibition Fast Assay Kit
KA0575	Histone Demethylase LSD1 Activity/Inhibition Assay Kit
KA0624	Histone H3 (Total) Quantification Kit (Colorimetric)
KA0625	Histone H3 (Total) Quantification Kit (Fluorometric)
KA0572	Histone Methyltransferase Activity/Inhibition Assay Kit (H3-K27)
KA0570	Histone Methyltransferase Activity/Inhibition Assay Kit (H3-K4)
KA0571	Histone Methyltransferase Activity/Inhibition Assay Kit (H3-K9)
KA0687	In Situ H3-K9 Acetylation Assay Kit
KA0707	In Situ Histone H3 Phosphorylation (Ser10) Assay Kit
KA0577	In Situ Histone H3-K27 Tri-Methylation Assay Kit
KA0689	In Situ Histone H3-K4 Methylation Assay Kit
KA0690	In Situ Histone H3-K9 Methylation Assay Kit
KA0675	MBD2 Binding Activity/Inhibition Assay Kit
KA0691	Methylated DNA Capture Kit
KA0676	Methylated DNA Quantification Kit




KA0911




KA0817




KA1643



KA0872



KA0810



KA0875


Fluorescent Protein Assay Kits

Cat. Num.	Product Name	Cat. Num.	Product Name
KA0911	GFP Quantitation Kit	KA1640	Luciferase Reporter Gene Assay Kit


Metabolism Assay Kits

Cat. Num.	Product Name	Cat. Num.	Product Name
KA1624	Acetylcholine Assay Kit	KA0817	Alkaline Phosphatase Assay Kit
KA0803	Acetyl-CoA Assay Kit	KA1642	Alkaline Phosphatase Assay Kit
KA0816	Acid Phosphatase Assay Kit (Colorimetric)	KA0820	Alkaline Phosphatase Assay Kit (Fluorometric)
KA0819	Acid Phosphatase Assay Kit (Fluorometric)	KA1643	alpha-Amylase Assay Kit
KA0876	Aconitase Assay Kit	KA0872	alpha-Ketoglutarate Assay Kit
KA0902	Adipogenesis Assay Kit	KA0810	Ammonia Assay Kit
KA1659	Adipolysis Assay Kit	KA0875	Amylase Assay Kit
KA1658	ADP Assay Kit	KA1660	Ascorbic Acid Assay Kit
KA0807	ADP Assay Kit (Colorimetric/Fluorometric)	KA1661	ATP Assay Kit
KA1673	ADP/ATP Ratio Assay Kit	KA0806	ATP Assay Kit (Colorimetric/Fluorometric)
KA1656	AF HDL and LDL/VLDL Assay Kit	KA0794	ATP Cell Viability Assay Kit


Metabolism Assay Kits




KA0870




KA0855




KA0832



KA0834



KA0836



KA0887


Cat. Num.	Product Name	Cat. Num.	Product Name
KA1610	ATPase/GTPase Assay Kit	KA0799	GST Activity Assay Kit (Colorimetric)
KA1611	beta-Glucosidase Assay Kit	KA0796	GST Activity Assay Kit (Fluorometric)
KA0854	beta-Hydroxybutyrate Assay Kit	KA1671	HDL and LDL/VLDL Assay Kit
KA0824	Branched Chain Amino Acid (Leu, Ile, Val) Assay Kit	KA0871	Hemin Assay Kit
KA1644	Calcium Assay Kit	KA0801	Hydrogen Peroxide Assay Kit
KA0812	Calcium Assay Kit (Colorimetric)	KA0814	Iron Assay Kit
KA0884	Catalase Assay Kit	KA1647	Iron Assay Kit
KA0790	CETP Activity Assay Kit	KA0865	Isocitrate Assay Kit
KA0828	CETP Inhibitor Drug Screening Kit	KA0833	Lactate Assay Kit
KA0823	Chloride Assay Kit	KA0834	Lactate Assay Kit II
KA1645	Chloride Assay Kit	KA0878	Lactate Dehydrogenase Assay Kit
KA0829	Cholesterol/Cholesteryl Ester Quantitation Kit	KA1653	Lactate Dehydrogenase Assay Kit
KA1662	Choline Assay Kit	KA1672	Lactose Assay Kit
KA0840	Choline/Acetylcholine Quantification Kit	KA0863	L-Alanine Assay Kit
KA0864	Citrate Assay Kit	KA0859	L-Amino Acid Quantitation Kit
KA0821	Cobalt Assay Kit	KA0860	L-Carnitine Assay Kit
KA0809	Coenzyme A Assay Kit	KA0877	Lipase Assay Kit
KA1664	Coenzyme A Assay Kit	KA1654	Lipase Assay Kit
KA0856	Creatine Assay Kit	KA0813	Magnesium Assay Kit
KA1666	Creatine Assay Kit	KA1650	Magnesium Assay Kit
KA1665	Creatine Kinase Assay Kit	KA0858	Malate Assay Kit
KA0849	Creatinine Assay Kit	KA0843	Maltose and Glucose Assay Kit
KA0869	D-Lactate Assay Kit (Colorimetric)	KA0851	Maltose Assay Kit
KA1646	Ethanol Assay Kit	KA1657	NAD+/NADH Assay Kit
KA0808	FAD Assay Kit	KA1663	NADP+/NADPH Assay Kit
KA0870	Ferric Reducing Ascorbate Assay Kit	KA0822	Nickel Assay Kit
KA0818	Fluorescent Phosphate Assay Kit	KA1641	Nitric Oxide Assay Kit
KA1667	Free Fatty Acid Assay Kit	KA0793	Nitric Oxide Assay Kit (Fluorometric)
KA0853	Free Glycerol Assay Kit	KA1634	Nitric Oxide Synthase Assay Kit
KA1668	Fructose Assay Kit	KA0867	Oxaloacetate Assay Kit
KA0855	Fumarate Assay Kit	KA0826	Phenylalanine Assay Kit (Fluorometric)
KA0842	Galactose and Lactose Assay Kit	KA0815	Phosphate Assay Kit (Colorimetric)
KA1669	Galactose Assay Kit	KA0827	Phosphatidylcholine Assay Kit
KA0841	Glucose and Sucrose Assay Kit	KA1636	Phospholipase D Assay Kit
KA0831	Glucose Assay Kit	KA1635	Phospholipid Assay Kit
KA1648	Glucose Assay Kit	KA0791	PLTP Activity Assay Kit
KA0832	Glucose Assay Kit II	KA0830	PLTP Inhibitor Drug Screening Kit
KA0866	Glucose-6-Phosphate Assay Kit	KA0836	Pyruvate Assay Kit
KA1430	Glucose-6-Phosphate Assay Kit	KA1674	Pyruvate Assay Kit
KA0880	Glucose-6-Phosphate Dehydrogenase Assay Kit	KA0873	Pyruvate Kinase Assay Kit
KA1670	Glutamate Assay Kit	KA0857	Sarcosine Assay Kit
KA0879	Glutamate Dehydrogenase Assay Kit	KA1655	Sialic Acid Assay Kit
KA1627	Glutamine Assay Kit	KA0862	Starch Assay Kit
KA0800	Glutathione Assay Kit	KA0850	Sucrose Assay Kit
KA1649	Glutathione Assay Kit	KA0883	Thioredoxin Reductase Assay Kit
KA0797	Glutathione Assay Kit (Colorimetric)	KA0802	Total Antioxidant Capacity Assay Kit
KA0882	Glutathione Peroxidase Assay Kit	KA0847	Triglyceride Quantification Kit
KA1628	Glutathione Peroxidase Assay Kit	KA1652	Urea Assay Kit
KA0881	Glutathione Reductase Assay Kit	KA1651	Uric Acid Assay Kit
KA0861	Glycogen Assay Kit	KA0874	Xanthine Oxidase Assay Kit

Signal Transduction Assay Kits

Cat. Num.	Product Name	Cat. Num.	Product Name
KA0885	Akt Activity Assay Kit	KA0888	JNK Activity Screening Kit
KA0900	beta-Secretase Activity Assay Kit	KA0903	PDE5A Assay Kit
KA0886	cAMP Direct Immunoassay Kit	KA0904	PI3 Kinase Assay Kit
KA0887	cGMP Direct Immunoassay Kit	KA0905	PTPN1 Assay Kit
KA0889	JNK Activity Assay Kit	KA0906	SPHK1 Assay Kit

Others Assay Kits


Cat. Num.	Product Name	Cat. Num.	Product Name
KA2097	Europium Calibration Kit	KA0895	Mitochondrial DNA Isolation Kit
KA1016	Hydrogen Peroxide Chemiluminescent Detection Kit	KA1074	PCSK9-LDLR in vitro Binding Assay Kit
KA1017	Hydrogen Peroxide Colorimetric Detection Kit	KA1015	sPLA2 activity Kit
KA0777	Mitochondria/Cytosol Fractionation Kit		



KA1017

Ab Concentration Kits


Cat. Num.	Product Name
KA1598	Antibody Concentration/Clean Up Kit




KA1563

Ab Labeling Kits


Cat. Num.	Product Name	Cat. Num.	Product Name
KA0001	Alkaline Phosphatase Labeling Kit-NH ₂	KA1575	FastLink Biotin Labeling Kit
KA0002	Alkaline Phosphatase Labeling Kit-SH	KA1576	FastLink Biotin Labeling Kit
KA0005	Allophycocyanin Labeling Kit-NH ₂	KA1574	FastLink B-Phycoerythrin Labeling Kit
KA0006	Allophycocyanin Labeling Kit-SH	KA1592	FastLink FluoProbes® 547H Labeling Kit
KA0003	Biotin Labeling Kit-NH ₂	KA1591	FastLink FluoProbes® 594 Labeling Kit
KA0004	Biotin Labeling Kit-SH	KA1593	FastLink FluoProbes® 647H Labeling Kit
KA0007	B-Phycoerythrin Labeling Kit-NH ₂	KA1594	FastLink FluoProbes® 682H Labeling Kit
KA0008	B-Phycoerythrin Labeling Kit-SH	KA1590	FastLink FluoProbes752 Labeling Kit
KA1549	FastLink Alkaline Phosphatase Labeling Kit	KA1571	FastLink Fluorescein Labeling Kit
KA1554	FastLink Allophycocyanin - XL Labeling Kit	KA1550	FastLink Glucose Oxidase Labeling Kit
KA1559	FastLink AMCA Labeling Kit	KA1551	FastLink HRP Labeling Kit
KA1560	FastLink Atto390 Labeling Kit	KA1581	FastLink PE-Atto594 Tandem Labeling Kit
KA1561	FastLink Atto425 Labeling Kit	KA1552	FastLink PerCP Labeling Kit
KA1562	FastLink Atto465 Labeling Kit	KA1582	FastLink PE-Texas Red Tandem Labeling Kit
KA1563	FastLink Atto488 Labeling Kit	KA1572	FastLink Rhodamine Labeling Kit
KA1564	FastLink Atto532 Labeling Kit	KA1553	FastLink R-PE Labeling Kit
KA1565	FastLink Atto565 Labeling Kit	KA1556	FastLink Streptavidin Labeling Kit
KA1566	FastLink Atto590 Labeling Kit	KA1573	FastLink Texas Red Labeling Kit
KA1567	FastLink Atto594 Labeling Kit	KA0009	Fluorescein Labeling Kit - NH ₂
KA1568	FastLink Atto633 Labeling Kit	KA0010	HiLyte Fluor™ 555 Labeling Kit-NH ₂
KA1583	FastLink Atto637 Labeling Kit	KA0011	HiLyte Fluor™ 647 Labeling Kit-NH ₂
KA1584	FastLink Atto655 Labeling Kit	KA0014	Peroxidase Labeling Kit-NH ₂
KA1569	FastLink Atto680 Labeling Kit	KA0015	Peroxidase Labeling Kit-SH
KA1570	FastLink Atto700 Labeling Kit	KA0012	R-Phycoerythrin Labeling Kit-NH ₂
KA1555	FastLink Avidin Labeling Kit	KA0013	R-Phycoerythrin Labeling Kit-SH



KA1568




KA1584




KA1591

Ab Purification Kits

Cat. Num.	Product Name	Cat. Num.	Product Name
KA1595	Antibody Purification Kit	KA1597	TCS Antibody Purification Kit
KA1596	Serum Antibody Purification Kit		



KA1593



KA0564



KA0567



KA0563



KA0565



KA0566



KA0562

IHC Kits

Cat. Num.	Product Name	Cat. Num.	Product Name
KA0564	BLACK Peroxidase Substrate Kit	KA0565	RealBlue Peroxidase Substrate Kit
KA0567	BLUE Phosphatase Substrate Kit	KA0566	RED Phosphatase Substrate Kit
KA0563	ORANGE Peroxidase Substrate Kit	KA0562	ReliableDAB Peroxidase Substrate Kit

IP & ChIP Kits

Cat. Num.	Product Name
KA0669	Acetyl-Histone H3 ChIP Kit
KA0670	Acetyl-Histone H4 ChIP Kit
KA0685	Chromatin Immunoprecipitation (ChIP) Kit
KA0692	Methylated DNA Immunoprecipitation Kit
KA0694	Methyl-CpG Binding domain Protein 2 ChIP Kit
KA0696	Methyl-Histone H3-K27 ChIP Kit
KA0697	Methyl-Histone H3-K4 ChIP Kit
KA0698	Methyl-Histone H3-K9 ChIP Kit
KA1514	One-Step ChIP Kit
KA0568	Protein A Agarose Immunoprecipitation Kit
KA0569	Protein G Agarose Immunoprecipitation Kit
KA0671	Tissue Acetyl-Histone H3 ChIP Kit
KA0672	Tissue Acetyl-Histone H4 ChIP Kit
KA0686	Tissue Chromatin Immunoprecipitation (ChIP) Kit
KA0693	Tissue Methylated DNA Immunoprecipitation Kit
KA0695	Tissue Methyl-CpG Binding domain Protein 2 ChIP Kit
KA0699	Tissue Methyl-Histone H3-K27 ChIP Kit
KA0700	Tissue Methyl-Histone H3-K4 ChIP Kit
KA0701	Tissue Methyl-Histone H3-K9 ChIP Kit
KA0702	Tri-Methyl-Histone H3-K9 ChIP Kit

Protein Expression Kits

Cat. Num.	Product Name
KA1434	Protein Expression Kit I
KA1435	Protein Expression Kit I (GST-Fusion Protein)
KA1436	Protein Expression Kit I (His-Fusion Protein)
KA1437	Protein Expression Kit II
KA1438	Protein Expression Kit II (GST-Fusion Protein)
KA1439	Protein Expression Kit II (His-Fusion Protein)

