



**PROGEN**  
*passion for research*

**SPEED UP!**

**His-tag  
detection**



**proflo**

His-tag quick test

IMMUNOCHROMATOGRAPHIC LATERAL FLOW TEST

Order No.: PRLFHIS



proflo  
HIS-TAG QUICK TEST

### Speed up detection of His-tagged proteins!


The **proflo** His-tag quick test detects histidine-tagged recombinant proteins in cell lysates or purified protein preparations. The semi-quantitative assay determines quickly & reliably the concentration of His-tagged proteins in your sample. Thus, it saves time for downstream steps and avoids time-consuming methods like SDS-PAGE or Western Blot.

### Use for universal applications!

- ✓ for small-scale screening
- ✓ during scale-up
- ✓ in purification fractions

### 2-step protocol – results in only 10 minutes!

- ✓ robust, immunochromatographic test principle
- ✓ competitive, semiquantitative assay
- ✓ concentration range  
2 – 50  $\mu\text{g/ml}$



### CONTACT

+49 6221 8278-0 • [info@progen.com](mailto:info@progen.com)  
[www.progen.com](http://www.progen.com)