



# A LEAP FORWARD

in standardized AAV titer determination

- Dependable & reproducible AAV titer quantification
- Specific recognition of intact rAAV particles
- Available for serotypes 1, 2, 3, 5, 6, 8, 9
- Validated assays for standardized approaches
- Damage-free rAAV purification with OptiPrep™



## SELECT FROM OUR PORTFOLIO OF AAV TITRATION ELISAs!

Serotype	Cat. No.
AAV1	PRAAV1
AAV2	PRATV
<b>NEW!</b> AAV2	PRAAV2R
<b>NEW!</b> AAV3 <i>new serotype!</i>	PRAAV3R
AAV5	PRAAV5
AAV6	PRAAV6
AAV8	PRAAV8
AAV9	PRAAV9

PROGEN offers a line of AAV quantification ELISAs for different AAV serotypes (1, 2, 3, 5, 6, 8, 9), utilizing its portfolio of capsid-specific AAV antibodies.

The assays for the determination of rAAV2 and rAAV8 titers have been validated in international studies\* and have been classified as superior method for a reliable, standardized rAAV particle titer quantification.

\*Ayuso E et al., 2014, Human Gene Therapy 25:977-987  
Lock M et al., 2010, Human Gene Therapy 21:1273-1285

PROGEN Biotechnik GmbH • Maaßstr. 30 • D-69123 Heidelberg • T +49 6221 8278-0

Contact us at [AAV@progen.com](mailto:AAV@progen.com)  
[www.progen.com](http://www.progen.com)

