



proteintech®

FlexAble Antibody Labeling Kits

Human kidney: anti-Calbindin (14479-1-AP) + FlexAble CoraLite Plus 550 Kit (KFA002, yellow), anti-ACE2 (66699-1-Ig) + FlexAble CoraLite Plus 650 Kit (KFA023, magenta), CoraLite@488-conjugated Podocalyxin antibody (CL488-18150, green), and DAPI (blue).

ptglab.com

FlexAble

Any antibody. Any color. Any time.



FlexAble is a novel antibody labeling kit for immunofluorescence, western blot and flow cytometry. It uses an affinity linker to conjugate fluorochromes, enzymes and molecules in any buffer condition.

Compatible with:

- Antibodies from any supplier
- Any antibody concentration
- Antibodies that contain BSA and glycerol

CONJUGATES

CoraLite® 488
CoraLite Plus 550
CoraLite Plus 650
CoraLite Plus 750

SPECIFICITIES

Rabbit IgG
Mouse IgG1

ptglab.com

