



New!

ImmPACT™ Vector® Red

Alkaline Phosphatase Substrate (SK-5105)

ULTRA SENSITIVE and **CONVENIENT**



Vector® Red Alkaline Phosphatase Substrate in our ImmPACT™ Formulation:

- Produces a bright red, insoluble, heat stable reaction product that is also highly fluorescent and non-fading
- Two to four times more sensitive than the original Vector® Red chromogen
- Supplied with a ready-to-use (R.T.U.), optimized diluent and the chromogen in a dropper bottle
- Can be used on tissue sections or cells, or on nitrocellulose, PVDF, or nylon membranes
- Ideal for single and multiple antigen labeling
- Can be non-aqueously or aqueously mounted

FIGURE LEGEND

Prostate: Anti-PSA (mouse monoclonal, VP-P981), ImmPRESS™ AP Anti-Mouse Ig Reagent (MP-5402), ImmPACT™ Vector® Red Alkaline Phosphatase Substrate (SK-5105), viewed under brightfield and fluorescence microscopy.