



# MERCODIA TOTAL GIP NL-ELISA NORTHERN LIGHTS

ORDER NOW

- **Range:** 2.7-1000 pmol/L
- **Sample volume:** 25  $\mu$ L
- **Species:** human
- **Detection:** chemiluminescence

## Mercodia's 3<sup>rd</sup> Northern Lights Assay

- High sensitivity
- LLOQ 2.7 pmol/L
- Broad measuring range
- No dilution of samples
- No cross-reactivity to any proglucagon derivate

## Why measure total GIP?

- Majority of circulating GIP is the DPP4-degraded inactive GIP(3-42).
- To get an overview of secretion, measure both GIP(1-42) and GIP(3-42).

