

FuGENE® SI Transfection Reagent

Multicomponent Transfection reagent for the delivery of siRNA into Mammalian Cells.

For Research Use Only. Not for use in therapeutic or diagnostic procedures

SKU	Size
SI-1000	1mL
SI-5000	5mL
SI-10000	10mL
SI-25000	25mL
SI-50000	50mL
SI-250000	250mL

Certificate of Analysis

Attribute	Requirement	Result
Ethanol Concentration	60-70%	PASS
Water Quality	16-18.2 Ω	PASS
Final Filtration	0.1 μ m filter	PASS
Appearance	Clear, colorless solution	PASS
Composition	Free from any ingredients of animal or human origin	PASS
Transfection Efficiency: Gene Knockdown as measured by Flow Cytometry	\geq 70% of Reference lot	PASS
Cell viability as measured by Flow Cytometry	\geq 90% of Reference lot	PASS

Storage Conditions

Product to be stored for long-term storage at 4 degrees Celsius for up to 2 years from date of delivery.

Product can be shipped at ambient temperature.

Limited Use Label License

You shall use this Product in accordance with all applicable laws, rules and regulations. You may not attempt to reverse engineer this Product. If you are not willing to accept these conditions, the distributor from whom the Product was purchased will accept the return of the unopened Product and provide you with a full refund. However, in the event this Product is opened, you agree to be bound by the conditions of this limited use statement.

FUGENT warrants and represents to that the Product shall, under the proper care and handling, be free from defects in materials and workmanship, and shall conform to the specifications on Datasheet for a period of 2 years from date of delivery. If any Products fail to meet the foregoing warranties, then, at the discretion of FUGENT, FUGENT shall either replace the non-conforming Products or refund the purchase price for such Products. FUGENT shall have the sole discretion to determine whether the non-conforming Products shall be returned to FUGENT.

The limited warranty set forth above is End User's sole and exclusive remedy against FUGENT for any claims, damages or other liabilities to the Product. EXCEPT FOR SUCH LIMITED WARRANTY, FUGENT MAKES NO OTHER WARRANTY, EXPRESS OR IMPLIED, AND SPECIFICALLY DISCLAIMS ALL **WARRANTIES OF NON-INFRINGEMENT, MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. FUGENT SHALL NOT BE LIABLE TO END USER FOR ANY CONSEQUENTIAL, SPECIAL, INDIRECT, PUNITIVE, EXEMPLARY OR INCIDENTAL DAMAGES ARISING OUT OF OR RELATED TO THE PRODUCT OR ITS PERFORMANCE, INCLUDING WITHOUT LIMITATION, ANY PERSONAL INJURY, PROPERTY DAMAGE, PRODUCTS LIABILITY, OR LOST PROFITS, EVEN IF FUGENT HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. IN NO EVENT SHALL FUGENT'S LIABILITY, IF ANY, FOR LOSS OR DAMAGES RELATING TO OR ARISING OUT OF THE PRODUCT EXCEED THE PURCHASE PRICE (OR PORTION THEREOF) ACTUALLY PAID BY END USER FOR THE DEFECTIVE PRODUCT(S) GIVING RISE TO SUCH LOSS OR DAMAGE.**

FuGENE® is licensed trademark of FuGENT LLC
1600 Aspen Commons #101
Middleton, WI 53562
www.fugene.com
For inquiries, contact us at contact@fugene.com

